




RUSCH® LARYNGOSCOPE REFERENCE GUIDE

We provide Rusch quality



MINIMIZE RISK. MAXIMIZE OUTCOMES.

Infection Control ... It's in your hands.

Teleflex is your partner in finding solutions. We believe that the tools you use should keep you and your patients safe, while providing uncompromised quality and reliability.

Disposable Blades

EQUIPLITE™, GREEN RUSCH® LITE™, AND RUSCH® LITE SLIM™

Single-use disposable blades that offer the performance of reusables, while presenting the opportunity to substantially reduce the potential for infection associated with reusable devices.

Quick Clean Handles

TRULED™ HANDLES

Easy-to-remove light and battery cartridge allows the handle to be quickly and easily sterilized between cases.

DOLPHIN™ HANDLES

Water resistant handles that you can leave intact during disinfection. Easier cleaning — better compliance.



Protecting patients, protecting you.

TELEFLEX®

Our mission is to enhance patient outcomes by providing products that help clinicians protect against infection and improve patient and provider safety.

RUSCH®

Teleflex is proud to continue a more than 125 year partnership between our Rusch brand and clinicians to develop laryngoscopes that make intubation as simple, convenient, and safe as possible. The result is one of the broadest ranges of laryngoscopes available under one brand, including options to manage difficult airways, and improve procedural efficiency. Each world-class product is precisely engineered to meet the needs of Anesthesia professionals and their patients.

FIBER OPTICS

FIBER OPTIC HANDLES

- Strong, durable stainless steel construction
- Bright Xenon lamp
- Textured to enhance grip



TRULED™ HANDLES

- Easy-to-remove battery and light cartridge for quick sterilization between cases
- Ergonomic cushioned grip
- LED lamp provides optimal blue-white light
- Rechargeable battery provides 10 hours of continuous use between charges
- 5-year warranty



SNAPLIGHT® BLADES

- Magnetized light pipe makes cleaning a “snap”

REUSABLE BLADES

SNAPLIGHT® FIBER OPTIC BLADES

- Unique snap-on light bundle eliminates screws
- Magnetized light bundle can be removed quickly for easy cleaning
- Integrated oxygen port available on infant sizes



EMERALD® FIBER OPTIC BLADES

- 5.0 mm diameter fiber optic bundle
- Integrated bundle design for easy cleaning
- Matte finish to reduce glare



DISPOSABLE BLADES

GREEN RUSCH LITE™ DISPOSABLE FIBER OPTIC BLADES

- Disposable fiber optic blade
- Solid metal construction from tip to heel
- Reliable “click” locking mechanism
- Anti-glare sheath on light pipe



STANDARD/CONVENTIONAL

METAL HANDLES




CROSS CUT FINISH

- Strong, durable construction
- Textured to enhance grip



MATTE FINISH

- Strong, durable construction
- Unique ribbed texture



DOLPHIN™ HANDLES

- Water resistant
- Handle can stay intact during disinfection
- Textured to enhance grip

REUSABLE BLADES

STANDARD BLADES

- Full range of styles and sizes
- Manufactured from surgical grade stainless steel
- Halogen light



MACINTOSH



MILLER



PHILLIPS



WISCONSIN




IMPROVED VIEW
MACINTOSH

DISPOSABLE BLADES

EQUIPLITE™ DISPOSABLE BLADES

- Disposable standard blades
- Solid metal construction from tip to heel
- LED bulb provides optimal blue-white light
- Color-coded by size
- “Click” locking mechanism



RUSCH LITE SLIM™ DISPOSABLE BLADES

- Disposable standard blades
- Rigid plastic construction, tensile tested for strength
- Bright Xenon light



SPECIALTY PRODUCTS

MRI KIT

- Independently certified to support MR conditions with a static magnetic field of 3 Tesla or less and highest gradient field of 720 Gauss/cm or less
- Unique design — black aluminum handle and gold-plated blades allow the set to be easily distinguishable
- Lightweight handle and blades
- Specially designed carbon zinc dry batteries



HYBRID BLADES



FOCS™ BLADES

- Fiber optic blades that work on standard handles
- Allows conversion of a standard handle to a fiber optic system without buying new handles
- Bright, cool focused light
- Easily removed light pipe allows for quick cleaning
- Protected water resistant Halogen in the base of the blade
- Integrated oxygen port available on infant sizes



FLIPPER™ BLADE

- FOCS style blade — fiber optic blade that works on standard handles
- Articulating tip aids in difficult intubations
- Reduces the force required to visualize the larynx
- Sold as a kit with a Dolphin™ handle, spare bulb, and carrying case

ORDERING INFORMATION

LARYNGOSCOPE BLADES										
DESCRIPTION	QTY. (BOX)	SIZE								
		00	0	1	1 1/2	2	3	3 1/2	4	5
		SMALL PREMATURE	PREMATURE	INFANT	INFANT	CHILD	MEDIUM ADULT	MEDIUM ADULT	LARGE ADULT	LARGE ADULT
FIBER OPTIC REUSABLE										
Emerald®										
Miller	1	004449000	004450050	004451100	004451150	004452200	004453300		004454400	
Macintosh	1			004431100		004432200	004433300	004433500	004434400	
SnapLight®										
Miller	1	002200000	002215000	002216100		002217200	002218300		002219400	
Macintosh	1			002201100		002202200	002203300		002204400	
Miller OxiPort	1	004457700	004458800							
FIBER OPTIC DISPOSABLE										
Green Rusch® Lite™										
Miller	20	004550010	004550000	004550001		004550002	004550003		004550004	
Macintosh	20			004551001		004551002	004551003	004551035	004551004	
Standard Reusable										
Miller	1	008600000	008615050	008616100		008617200	008618300		008619400	
Macintosh	1			008601100		008602200	008603300	008603350	008604400	
Improved View Mac	1		001967000	001967100		001965200	001963300	001966350	001964400	
Phillips	1			008607100		008608200				
Wisconsin	1			008611100		008612200	008613300		008614400	
STANDARD DISPOSABLE										
EquipLite™										
Miller	20	004650010	004650000	004650001		004650002	004650003		004650004	
Macintosh	20		004651000	004651001		004651002	004651003	004651035	004651004	004651005
Rusch® Lite Slim™										
Miller	20		004850050	004851100		004852200	004853300			
Macintosh	20			004801100		004802200	004803300		004804400	
HYBRID (FIBER OPTIC LIGHT FOR A STANDARD HANDLE)										
FOCS™										
Miller	1	005648000	005650050	005651100		005652200	005653300		005654400	
Macintosh	1			005640100		005641200	005642300	005645350	005643400	
Miller OxiPort	1	005648700	005650710	005651720						
Flipper™	1 Kit (blade, medium handle, spare bulb, carrying case)							005640000		

LARYNGOSCOPE HANDLES					
DESCRIPTION	QTY. (BOX)	SIZE			
		PEDIATRIC/ SHORT	SMALL/ PENLIGHT	STUBBY/ SHORT	MEDIUM/ STANDARD
FIBER OPTIC					
Fiber optic	1		004412200	004413300	004411100
TruLED™	1				0055500
TruLED™ Replacement Battery					0055510
TruLED™ Replacement Charger					0055530
STANDARD					
Cross Cut Finish	1	008710000	008620100	008720050	008621000
Matte Finish	1		142600		142500
Dolphin™	1		145100		145000

RUSCH® MRI		
DESCRIPTION	QTY.	PRODUCT #
Rusch® MRI Kit	1	005852300
Macintosh Blades 2, 3, 4		
Miller Blades 1, 2, 3		
1 Medium Fiber Optic Handle		
2 Battery Packs		
Rusch® MRI Extra Batteries	1	005853300

STANDARD/CONVENTIONAL BULBS		
PRODUCT #	FITS PRODUCTS	SIZE
008623100	Miller	
008600000		00
008615050		0
008616100		1
008617200		2
008618300		3
Macintosh		
008601100		1
008602200		2
008603300		3
Wisconsin		
008611100		1
Improved View		
001967000		0
001967100		1
001965200		2
Phillips		
008607100		1
008622300	Miller	
008619400		4
Macintosh		
008603350		3.5
008604400		4
Wisconsin		
008612200		2
008613300		3
008614400		4
Improved View		
001963300		3
001966350		3.5
001964400		4
Phillips		
008608200		2

FIBER OPTIC LAMPS		
PRODUCT #	FITS PRODUCTS	SIZE
005630000	004412200	Small/Penlight
	004411100	Medium/Standard
	004413300	Stubby/Short

FOCS™ ILS FIBER BUNDLES		
PRODUCT #	FITS PRODUCTS	SIZE
	Miller	
005601000	005648000	00
005601050	005650050	0
005601100	005651100	1
005601200	005652200	2
005601300	005653300	3
005601400	005644000	4
Macintosh		
005600100	005640100	1
005600200	005641200	2
005600300	005642300	3
005600350	005645350	3.5
005600400	005643400	4
Miller-OxiPort		
005602000	005648700	00
005602050	005650710	0
005602100	005651720	1
Flipper™		
005600350	005640000	3.5

SNAPLIGHT® FIBER BUNDLES		
PRODUCT #	FITS PRODUCTS	SIZE
	Macintosh	
002501100	002201100	1
002502200	002202200	2
002503300	002203300	3
002504400	002204400	4
Miller		
002500000	002200000	00
002515000	002215000	0
002516100	002216100	1
002517200	002217200	2
002518300	002218300	3
002519400	002219400	4
Miller-OxiPort		
002500000	004457700	00
002515000	004458800	0

Rx Only
Not manufactured with natural rubber latex.



VISIT MYRUSCH.COM FOR OUR FULL LINE OF RUSCH® AIRWAY MANAGEMENT PRODUCTS
TO PLACE AN ORDER CONTACT CUSTOMER SERVICE AT 1.866.246.6990

TruLED and EquipLite are trademarks of Truphatek International Ltd.
Teleflex, Dolphin, Emerald, Flipper, FOCS, Green Rusch Lite, Rusch, Rusch Lite Slim and SnapLight are trademarks or registered trademarks of Teleflex Incorporated or its affiliates.

Teleflex is a global provider of medical products designed to enable healthcare providers to protect against infections and improve patient and provider safety. The company specializes in products and services for vascular access, respiratory, general and regional anesthesia, cardiac care, urology and surgery. Teleflex also provides specialty products for device manufacturers.

© 2013 Teleflex Incorporated. All rights reserved. 2012-1580



TELEFLEX MEDICAL PO Box 12600 Research Triangle Park, NC 27709
Toll Free: 866.246.6990 Phone: 919.544.8000 Intl: 919.433.8088
TELEFLEX.COM