

**ETHICON**

**NUROLON\***

**Braided Polyamide 6/6 Suture**



# NUROLON\*

## Braided Polyamide 6/6 Suture

### Braided

For easy handling and secure knot tying.

### Coated

With wax to minimise capillarity.

### Synthetic

For minimal tissue reaction.

### Non-Absorbable

Provides prolonged tensile strength retention in tissue.

### Colour

Black.

### Range










Gauge sizes 5-0 to 1 (U.S.P.). Supplied as needled sutures and ligatures.

### Sterilisation

By Irradiation.

These sutures comply with the requirements of U.S. Pharmacopoeia (U.S.P.) and The European Pharmacopoeia (PhEUR).

# NUROLON\* Braided Polyamide 6/6

Needle length and description		Length	Colour	Gauge						
<b>⊙ 1/2 circle Taper Point</b>				3-0 2	2-0 3	0 3.5	1 4	2 5	4 6	5 7
	26mm	8x45cm CR	black	<b>C513G</b>						
	26mm	8x75cm CR	black	<b>C5436G</b>						
<b>⊖ 1/2 circle Taper Point plus</b>				3-0 2	2-0 3	0 3.5	1 4	2 5	4 6	5 7
	26mm	75cm	black		<b>NL833H</b>					
<b>⊙ 1/2 circle Taper Point (heavier wire diameters)</b>				3-0 2	2-0 3	0 3.5	1 4	2 5	4 6	5 7
	31mm	75cm	black				<b>W5723</b>			
<b>⊙ J Needle Taper Point</b>				3-0 2	2-0 3	0 3.5	1 4	2 5	4 6	5 7
	31mm	75cm	black				<b>W5985</b>			
<b>⊖ 1/2 circle TAPERCUT*</b>				3-0 2	2-0 3	0 3.5	1 4	2 5	4 6	5 7
	40mm	100cm	black				<b>W5415</b>			
<b>▼ 1/2 circle Reverse Cutting</b>				3-0 2	2-0 3	0 3.5	1 4	2 5	4 6	5 7
	22mm	75cm	black				<b>W5957</b>			
Needle length and description		Length	Colour	Gauge						
<b>▲ 3/8 circle Conventional Cutting</b>				3-0 2	2-0 3	0 3.5	1 4	2 5		
	26mm	45cm	black		<b>W5534</b>					
<b>▼ Straight Cutting</b>				3-0 2	2-0 3	0 3.5	1 4	2 5		
	55mm	100cm	black		<b>W5775</b>					

# NUROLON\* Braided Polyamide 6/6

Ligatures and Non-needled Sutures			
SUTUPAK* Pre-cut Sutures			
Gauge		black 17 x 45cm lengths 12 sachets per box	black 10 x 75cm lengths 12 sachets per box
<b>3-0</b>	2		<b>W5222</b>
<b>2-0</b>	3	<b>W5203</b>	<b>W5223</b>
<b>0</b>	3.5		<b>W5224</b>
<b>1</b>	4	<b>W5205</b>	<b>W5225</b>