

Introducing the **OxiMax**[®] Pulse Oximetry Platform

New Max Sensors

OXIMAX
DIGITAL
MEMORY CHIP

SENSOR
MESSAGES
AND
SENSOR
EVENT
REPORT

IMPROVED
PERFORMANCE



tyco

Healthcare

NELLCOR[®]

Max Sensor Product Summary

	OxiMax Model #	Equivalent RCal Model #	OxiMax Features
	MAX-A-I	D-25	<ul style="list-style-type: none"> OxiMax digital memory chip in MAX sensors Sensor ID: Lot Code/Sensor Model, Signature Verification Sensor Messages and Sensor Event Report Tear resistant bandage and ambient light performance (opaque layer) in MAX-A-I, MAX-P-I, MAX-AL-I Forwards and backwards compatible Dark lavender plug
	MAX-N-I	N-25	
	MAX-P-I	D-20	
	MAX-I-I	I-20	
	MAX-AL-I	D-25L	
	MAX-R-I	R-15	
Oxisensor II			<ul style="list-style-type: none"> Will remain R-Cal compatible only.
	OC-3	OC-3	<ul style="list-style-type: none"> OxiMax digital memory chip in OC-3 Sensor ID: Lot Code/Sensor Model, Signature Verification Sensor Messages (no Sensor Event Report) Forwards and backwards compatible Dark lavender plug on OC-3
	DS100A	DS100A	<ul style="list-style-type: none"> OxiMax digital memory chip in DS100A, Dura-Y and Oxibands Sensor ID: Lot Code/Sensor Model, Signature Verification Sensor Messages (no Sensor Event Report) Forwards and backwards compatible Dark lavender plug on OC-3
	D-YS	D-YS	
	OXI A/N	OXI A/N	
	OXI P/I	OXI P/I	
RS-10	N/A	RS-10	<ul style="list-style-type: none"> Will remain R-Cal compatible only.
	MAX-FAST	N/A	<ul style="list-style-type: none"> OxiMax digital memory chip Sensor ID: Lot Code/Sensor Model, Signature Verification Sensor Messages and Sensor Event Report Exclusively OxiMax compatible White plug
Patient Cables	DEC4	EC4	<ul style="list-style-type: none"> Forwards and backwards compatible Dark lavender connectors on both ends Exclusively OxiMax compatible Dark lavender connectors on both ends
	DEC8	EC8	
	DOC10	N/A	
	OEM10		

All products are trademarks of Nellcor.

© 2001 Tyco Healthcare. All rights reserved.

M-MAXSENS/GB

tyco
Healthcare

NELCOR

Tyco Healthcare UK Limited
154 Fareham Road
GOSPORT
Hampshire
PO13 0AS
UK
Tel: +44 (0) 1329 224226
Fax: +44 (0) 1329 224334

Tyco Healthcare
Birchwood Court
Montrose Road
P.O. Box 8108
Halfway House
1685 MIDRAND
SOUTH AFRICA
Tel: +27 (0) 11 265 2000
Fax: +27 (0) 11 265 2050

Tyco Healthcare Regional Export
154 Fareham Road
GOSPORT
Hampshire
PO13 0AS
UK
Tel: +44 (0) 1329 224258
Fax: +44 (0) 1329 224086

Tyco Healthcare Norden AB
Gårdsvägen 18
SE-169 67 SOLNA
SWEDEN
Tel: +46 (0) 8 585 605 05
Fax: +46 (0) 8 585 605 06

Tyco Healthcare Nederland BV
Hogeweg 105
5301 LL ZALTBOMMEL
The NETHERLANDS
Tel: +31 418 57 66 00
Fax: +31 418 57 67 93
TycoHealthNed@TycoHealth.com
www.tycohealthcare.nl

Tyco Healthcare Belgium N.V.-S.A.
Generaal de Wittelaan 9/5
B-2800 MECHELEN
BELGIUM
Tel: +32 (0) 15 29 44 50
Fax: +32 (0) 15 29 44 55
tycohealthcare.belgium@tycohealth.be