


RÜSCH SURGERY AND GASTROENTEROLOGY

Quality with a fine tradition



QUALITY WITH A FINE TRADITION

Trusted brands make Teleflex a reliable and strong partner. Built on a solid tradition of innovation, Teleflex is a global leader in superior medical supplies designed to help providers minimise risk and maximise outcomes for their patients. Our understanding of the importance to our customers of a full range of products has led to the development of a unique line of medical devices, all of which complement one another.

Teleflex's RÜSCH instruments for surgery and gastroenterology supply the medical customer with user-friendly products which meet general market demands and provide safety and efficiency for patient and medical staff. Whether made of soft rubber, latex, silicone or PVC, choose from our large range of drains, T-tubes, feeding and stomach tubes, rectal tubes, tubings and various patient care products to find the right instrument for your requirements. Our comprehensive quality management system ensures maximum reliability and functionality of all products.

With 125 years of know-how, product development and experience, the RÜSCH brand has gained a worldwide reputation for quality and performance.


You can find further details and the technical specifications of our products within this brochure.

TELEFLEX – HIGH QUALITY MEDICAL SUPPLIES FROM A SINGLE SOURCE



| | |
|-----------------------------|----|
| Drains | 4 |
| T-tubes | 6 |
| Feeding tubes | 7 |
| Stomach tubes | 8 |
| Rectal tubes | 10 |
| Tubing | 12 |
| Patient care products | 14 |
| Pessaries | 15 |

DRAINS



PENROSE DRAIN
made of Silkotex Rüsich Gold®

- length approx. 41 cm
- wall thickness approx. 0.13 mm
- sterile

| PENROSE DRAIN | | | | | Ⓢ |
|---------------|------------|-----------|------------------|--------------|-----|
| REF. | ORDER SIZE | SIZE INCH | SIZE | DIAMETER | QTY |
| 473600 | size 1 | 1/4" | approx. 6 1/3 mm | Ø 10 mm flat | 10 |
| | size 2 | 5/16" | approx. 8 mm | Ø 13 mm flat | |
| | size 3 | 3/8" | approx. 10 mm | Ø 16 mm flat | |
| | size 4 | 1/2" | approx. 13 mm | Ø 20 mm flat | |
| | size 5 | 5/8" | approx. 16 mm | Ø 25 mm flat | |
| | size 6 | 3/4" | approx. 19 mm | Ø 30 mm flat | |
| | size 7 | 7/8" | approx. 22 mm | Ø 34 mm flat | |
| | size 8 | 1" | approx. 25 mm | Ø 39 mm flat | |



CORRUGATED DRAIN
made of soft rubber

- red
- autoclavable
- non-sterile

| CORRUGATED DRAIN | | | Ⓢ |
|------------------|--------------------|--|-----|
| REF. | SIZE | | QTY |
| 473800* | approx. 30 x 15 cm | | 10 |



CORRUGATED DRAIN
made of silicone

- white
- radiopaque
- latex-free
- sterile, double-packed

| CORRUGATED DRAIN | | | Ⓢ |
|------------------|----------------------|--|-----|
| REF. | SIZE | | QTY |
| 473901 | approx. 30 x 12.5 cm | | 10 |

* This product contains natural rubber latex which may cause allergic reactions.

ADHESIVE DRAINAGE BAG



| REF. | VOLUME OF BAG | QTY |
|--------|---------------|-----|
| 173604 | 200 ml | 5 |



ADHESIVE DRAINAGE BAG

- Latex-free adhesive drainage bag with skin protection and micro-skin adhesive surface:
- anti-reflux valve
 - outlet
 - sterile, double-packed

PASSIVE WOUND DRAINAGE SYSTEM



| REF. | ORDER SIZE (CH) | | | | DESCRIPTION | QTY |
|--------|-----------------|----|----|----|--|-----|
| 179410 | 12 | 15 | 18 | 21 | length of drain approx. 100 cm volume of bag 350 ml | 5 |
| | 24 | 27 | 30 | | | |



PASSIVE WOUND DRAINAGE SYSTEM

- Latex-free circular drain made of silicone:
- atraumatic drainage tip
 - five lateral eyes, staggered
 - X-ray marker throughout entire length

- Drainage bag:
- anti-reflux valve
 - outlet valve
 - universal suspension
 - sterile, double-packed



LIGHTFLOW CAPILLARY DRAIN



| REF. | ORDER SIZE (O.D.) MM | | | | LENGTH | QTY |
|--------|----------------------|---|----|----|---------------|-----|
| 173600 | 6 | 8 | 10 | 12 | approx. 30 cm | 10 |



LIGHTFLOW CAPILLARY DRAIN
made of silicone

- transparent
- asymmetrical ribs on inner wall
- X-ray marker throughout entire length
- latex-free
- sterile, double-packed

PASSIVE CIRCULAR DRAIN



| REF. | ORDER SIZE (CH) | | | | | | | | LENGTH | QTY | | |
|--------|-----------------|----|----|----|----|----|----|----|--------|-----|---------------|----|
| 179401 | 10 | 12 | 15 | 18 | 21 | 24 | 27 | 30 | 36 | 39 | approx. 50 cm | 10 |

PASSIVE CIRCULAR DRAIN
made of silicone


- atraumatic drainage tip
- five lateral eyes, staggered
- X-ray marker throughout entire length
- latex-free
- sterile, double-packed

T-TUBES

RÜSCH GOLD® T-TUBE, KEHR TYPE

made of Silkolatex


- without funnel
- sterile, double-packed



T-TUBE, KEHR TYPE

made of soft rubber

- without funnel
- sterile, double-packed



DETACHABLE FUNNEL

made of PVC

- for T-Tubes
- sterile

RÜSCH GOLD® T-TUBE

| REF. | ORDER SIZE (O.D.) MM | | | | | LENGTH | QTY |
|---------|----------------------|-----|-----|-----|-----|------------|-----|
| 178700* | 3.0 | 3.5 | 4.0 | 4.5 | 5.0 | 50 x 16 cm | 10 |
| | 5.5 | 6.0 | 7.0 | 8.0 | | | |
| 178900* | 3.5 | 4.0 | 4.5 | 5.0 | 5.5 | 70 x 16 cm | 10 |
| | 6.0 | 7.0 | 8.0 | | | | |

T-TUBE

| REF. | ORDER SIZE (O.D.) MM | | | | | | LENGTH | QTY |
|---------|-----------------------|-----|-----|-----|-----|------|------------|-----|
| 423600* | 2.5 | 3.0 | 3.5 | 4.0 | 4.5 | 5.0 | 70 x 16 cm | 10 |
| | 5.5 | 6.0 | 7.0 | 8.0 | 9.0 | 10.0 | | |
| 423601* | 3.0 | 4.0 | 5.0 | 6.0 | | | 70 x 16 cm | 5 |
| | open gutter cross bar | | | | | | | |

DETACHABLE FUNNEL

| REF. | ORDER SIZE FUNNEL SIZES (O.D.) | IN COMBINATION WITH T-TUBE MADE OF SILKOLATEX REF. 178700* / 178900* | IN COMBINATION WITH T-TUBE MADE OF SOFT RUBBER REF. 423600* / 423601* | QTY |
|--------|--------------------------------|--|---|-----|
| 333169 | 1.8 mm | O.D. 3.0 mm | O.D. 2.5 - 3.0 mm | 10 |
| | 2.5 mm | O.D. 3.5 - 4.0 mm | O.D. 3.5 - 4.0 mm | |
| | 3.0 mm | - | O.D. 4.5 mm | |
| | 3.5 mm | O.D. 4.5 - 5.0 mm | O.D. 5.0 mm | |
| | 4.0 mm | - | O.D. 5.5 - 6.0 mm | |
| | 4.5 mm | O.D. 5.5 - 6.0 mm | O.D. 6.5 mm | |
| | 5.0 mm | - | O.D. 7.0 mm | |
| | 5.5 mm | O.D. 7.0 - 8.0 mm | O.D. 8.0 / 9.0 / 10.0 mm | |

* These products contain natural rubber latex which may cause allergic reactions.

FEEDING TUBES



FEEDING TUBE



| REF. | ORDER SIZE | 1 | 2 | 3 | 4 | 5 | LENGTH | QTY |
|--------|--------------|--|---|---|--|--|---------------|-----|
| 224000 | I.D./O.D. MM | 1.0 / 1.5 | 1.5 / 2.1 | 1.7 / 2.5 | 2.0 / 3.0 | 2.0 / 3.5 | approx. 50 cm | 50 |
| | COLOUR |  yellow |  green |  red |  blue |  orange | | |

FEEDING TUBE

made of PVC, latex-free

- transparent
- plastic fitting with stopper, suitable for Luer adaptor
- sizes colour-coded
- cylindrical, 2 eyes
- siliconised
- sterile



FEEDING TUBE



| REF. | ORDER SIZE | 1 | 2 | 3 | 4 | 5 | LENGTH | QTY |
|--------|------------------------------|---|---|---|---|---------------|--------|-----|
| 334000 | colour-coded see Ref. 224000 | | | | | approx. 50 cm | 10 | |

FEEDING TUBE

made of Wiruthan, latex-free

- for prolonged periods of use
- funnel with stopper
- cylindrical, 2 eyes
- sterile



FEEDING TUBE



| REF. | ORDER SIZE | 1 | 2 | 3 | 4 | 5 | LENGTH | QTY |
|--------|------------------------------|---|---|---|---|---------------|--------|-----|
| 235500 | colour-coded see Ref. 224000 | | | | | approx. 50 cm | 100 | |

FEEDING TUBE

made of PVC, latex-free

- transparent
- plastic fitting with stopper, suitable for Luer adaptor
- cylindrical, 2 eyes
- siliconised
- sterile



FEEDING TUBE



| REF. | ORDER SIZE | CH. 12 | CH. 16 | CH. 18 | LENGTH | QTY |
|--------|------------|---|--|---|----------------|-----|
| 334102 | ORDER SIZE | CH. 12 | CH. 16 | CH. 18 | approx. 120 cm | 10 |
| | COLOUR |  white |  orange |  red | | |

FEEDING TUBE

made of Wiruthan, latex-free

- for prolonged periods of use
- with X-ray marker throughout entire length
- double funnel for Luer and universal adaptor
- cylindrical, 4 eyes
- sterile



FEEDING TUBE



| REF. | ORDER SIZE | CH. 12 | CH. 14 | CH. 16 | CH. 18 | LENGTH | QTY |
|--------|------------|---|---|--|---|----------------|-----|
| 334101 | ORDER SIZE | CH. 12 | CH. 14 | CH. 16 | CH. 18 | approx. 120 cm | 10 |
| | COLOUR |  white |  green |  orange |  red | | |

FEEDING TUBE

made of Wiruthan, latex-free

- for prolonged periods of use
- with X-ray marker throughout entire length
- funnel without stopper
- cylindrical, 2 eyes
- sterile

STOMACH TUBES



SENGSTAKEN TUBE
made of soft rubber

- balloons made of Silkolatex Rüschi-Gold®, with pilot balloons
- Luer syringe tips
- radiopaque ring between the balloons
- shaft completely radiopaque
- graduated
- cylindrical tip
- 4 suction eyes at the distal end of the tube
- sterile

SENGSTAKEN TUBE

| REF. | ORDER SIZE (CH.) | | | DESCRIPTION | LENGTH | QTY |
|---------|------------------|----|----|-------------|---------------------------------|-----|
| 204800* | 16 | 18 | 21 | three-way | for adults in length of 100 cm | 2 |
| 204801* | 14 | | | three-way | for children in length of 65 cm | 2 |



SENGSTAKEN-BLAKEMORE TUBE
made of soft rubber

- balloons made of Silkolatex Rüschi-Gold®, with pilot balloons
- Luer syringe tips
- shaft completely radiopaque
- graduated
- cylindrical tip
- 4 suction eyes at the distal end of the tube
- sterile

SENGSTAKEN-BLAKEMORE TUBE

| REF. | ORDER SIZE (CH.) | | | DESCRIPTION | OVERALL LENGTH | QTY |
|--|------------------|----|----|-------------|----------------|-----|
| 204802* | 14 | 16 | 18 | three-way | approx. 100 cm | 2 |
| 204803* | | | 18 | four-way | approx. 100 cm | 2 |
| with additional suction option through the oesophageal balloon | | | | | | |



LINTON-NACHLAS TUBE
made of soft rubber

- balloons made of Silkolatex Rüschi-Gold®, with pilot balloons
- Luer syringe tips
- shaft completely radiopaque
- sterile

LINTON NACHLAS TUBE

| REF. | ORDER SIZE (CH.) | DESCRIPTION | LENGTH | QTY |
|---------|------------------|-------------|----------------|-----|
| 204700* | 21 | three-way | approx. 100 cm | 1 |

* These products contain natural rubber latex which may cause allergic reactions.



OESOPHAGEAL TUBE

| REF. | ORDER SIZE (O.D.) MM | | | | | | | | LENGTH | QTY |
|---------|----------------------|------|------|------|------|------|------|------|----------------|-----|
| 445300* | 3.0 | 4.0 | 5.0 | 6.0 | 7.0 | 8.0 | 9.0 | 10.0 | approx. 80 cm | 10 |
| | 11.0 | 12.0 | 13.0 | 14.0 | 16.0 | 18.0 | | | | |
| 447100* | 5.0 | 6.0 | 7.0 | 8.0 | 9.0 | 10.0 | 11.0 | 12.0 | approx. 150 cm | 1 |

OESOPHAGEAL TUBE
made of soft rubber

- with funnel
- cylindrical solid tip, 2 eyes
- siliconised



OESOPHAGEAL TUBE

| REF. | ORDER SIZE (O.D.) MM | | | | | | | | LENGTH | QTY |
|--------|----------------------|-----|-----|-----|-----|-----|------|--|---------------|-----|
| 225300 | 4.0 | 5.0 | 6.0 | 7.0 | 8.0 | 9.0 | 10.0 | | approx. 75 cm | 10 |

OESOPHAGEAL TUBE
made of PVC, latex-free

- transparent
- with funnel
- cylindrical hollow tip, 2 eyes
- siliconised
- sterile



OESOPHAGEAL BOUGIE

| REF. | ORDER SIZE (CH.) | | | | | | | | LENGTH | QTY |
|--------|------------------|----|----|----|----|----|----|----|---------------|-----|
| 282300 | 7 | 8 | 9 | 10 | 11 | 12 | 13 | 14 | approx. 72 cm | 1 |
| | 15 | 16 | 17 | 18 | 19 | 20 | 21 | 22 | | |
| | 23 | 24 | 25 | 26 | 27 | 28 | 29 | 30 | | |
| | 32 | 34 | 36 | 38 | 40 | | | | | |

OESOPHAGEAL BOUGIE
made of PVC, latex-free

- brown
- with olive tip



MILLER-ABBOTT TUBE

| REF. | ORDER SIZE (CH.) | | | LENGTH | QTY |
|---------|------------------|----|----|----------------|-----|
| 205000* | 12 | 14 | 16 | approx. 310 cm | 2 |

MILLER-ABBOTT TUBE
made of soft rubber

- with latex Rüschi Gold® balloon
- radiopaque
- graduated
- with metal fitting



LEVIN DUODENAL TUBE

| REF. | ORDER SIZE (CH.) | | | | | | | LENGTH | QTY |
|--------|------------------|----|----|----|----|----|--|----------------|-----|
| 235200 | 8 | 10 | 12 | 14 | 16 | 18 | | approx. 125 cm | 50 |

LEVIN DUODENAL TUBE
made of PVC, latex-free

- transparent
- cylindrical, 4 eyes
- graduated
- sterile

RECTAL TUBES




NORDMANN BALLOON RECTAL TUBE

made of soft rubber

- Silkolatex coating
- Silkolatex Rüsich Gold® balloon and pilot balloon
- with funnel
- two-way
- universal syringe adaptor
- cylindrical hollow tip, 2 eyes

NORDMANN BALLOON RECTAL TUBE

| REF. | ORDER SIZE (O.D.) MM | | | | LENGTH | QTY |
|---------|--|-----|------|------|---------------|-----|
| 209000* | 6.0 | 8.0 | 10.0 | 12.0 | approx. 30 cm | 5 |
| | balloon content 75 ml (max. capacity 300 ml) | | | | | |



GAUL DOUBLE BALLOON RECTAL TUBE

made of soft rubber

- Silkolatex coating
- Silkolatex Rüsich Gold® balloon and pilot balloon
- with funnel
- two-way
- universal syringe adaptor
- cylindrical hollow tip, 2 eyes

GAUL DOUBLE BALLOON RECTAL TUBE

| REF. | ORDER SIZE (O.D.) MM | LENGTH | QTY |
|---------|----------------------|---------------|-----|
| 209200* | 10.0 | approx. 30 cm | 1 |



BALLOON RECTAL TUBE

made of PVC

- transparent
- two-way
- high-volume low-pressure cuff
- X-ray marker
- latex-free
- sterile

BALLOON RECTAL TUBE

| REF. | ORDER SIZE (O.D.) MM | LENGTH | QTY |
|--------|----------------------|---------------|-----|
| 253475 | 10.0 | approx. 30 cm | 5 |

* These products contain natural rubber latex which may cause allergic reactions.



RECTAL TUBE

| REF. | ORDER SIZE (O.D.) MM | | | | | | | | | LENGTH | QTY |
|---------|----------------------|------|-----|-----|------|------|------|------|------|---------------|-----|
| 437500* | 6.0 | 7.0 | 8.0 | 9.0 | 10.0 | 11.0 | 12.0 | 13.0 | 14.0 | approx. 30 cm | 10 |
| | 15.0 | 16.0 | | | | | | | | approx. 30 cm | 5 |

RECTAL TUBE

made of soft rubber

- with funnel
- cylindrical solid tip, 2 eyes



RECTAL TUBE

| REF. | ORDER SIZE (O.D.) MM | | | | | | | | | LENGTH | QTY |
|---------|----------------------|-----|-----|-----|------|------|------|--|--|---------------|-----|
| 431100* | 4.0 | 5.0 | 6.0 | 7.0 | 8.0 | | | | | approx. 10 cm | 10 |
| 431600* | 6.0 | 7.0 | 8.0 | 9.0 | 10.0 | 11.0 | 12.0 | | | approx. 20 cm | 10 |

RECTAL TUBE

made of soft rubber

- with funnel
- central opening



RECTAL TUBE

| REF. | ORDER SIZE (CH.) | | | | | | | | | | LENGTH | QTY |
|---------|----------------------|-----|-----|------|------|------|----|----|----|---------------|---------------|-----|
| 436700* | 12 | 14 | 16 | 18 | 20 | 22 | 24 | 26 | 28 | | approx. 50 cm | 10 |
| | 30 | 32 | 34 | 36 | | | | | | | | |
| REF. | ORDER SIZE (O.D.) MM | | | | | | | | | LENGTH | QTY | |
| 436400* | 7.0 | 8.0 | 9.0 | 10.0 | 11.0 | 12.0 | | | | approx. 30 cm | 10 | |

RECTAL TUBE

made of soft rubber

- with funnel
- central opening, 1 eye



RECTAL TUBE


| REF. | ORDER SIZE (O.D.) MM | | | | | LENGTH | QTY |
|--------|---|--|--|---|--|---------------|-----|
| 235800 | 6.0 | 7.0 | 8.0 | 9.0 | 10.0 | approx. 30 cm | 50 |
| 235900 |  red |  violet |  blue |  green |  grey | approx. 40 cm | |

RECTAL TUBE

made of PVC

- transparent
- funnel according to colour code
- cylindrical, 2 eyes
- latex-free
- sterile

TUBING



SILKOMED® TUBING

made of silicone

- pyrogen-free
- autoclavable
- latex-free
- non-sterile



SILKOMED® TUBING

| REF. | ORDER SIZE LENGTH | I.D. MM | O.D. MM | SIZE LENGTH | I.D. MM | O.D. MM |
|--------|-------------------|---------|---------|--------------|---------|---------|
| 471600 | approx. 15 m | 2.0 | x 4.0 | approx. 15 m | 5.0 | x 9.0 |
| | | 3.0 | x 5.0 | | 6.0 | x 9.0 |
| | | 3.0 | x 7.0 | | 6.0 | x 12.0 |
| | | 4.0 | x 6.0 | | 7.0 | x 10.0 |
| | | 4.0 | x 7.0 | | 8.0 | x 12.0 |
| | | 5.0 | x 7.0 | | 10.0 | x 14.0 |
| QTY | 1 coil | 5.0 | x 8.0 | | 11.0 | x 18.0 |

DRAINAGE TUBING

made of soft rubber

- autoclavable
- non-sterile



DRAINAGE TUBING

| REF. | ORDER SIZE LENGTH | I.D. MM | O.D. MM | SIZE LENGTH | I.D. MM | O.D. MM | |
|---------|-------------------|--------------|---------|--------------|--------------|--------------|--------|
| 474100* | approx. 55 m | 1.5 | x 3.0 | approx. 30 m | 6.0 | x 10.0 | |
| | | 2.0 | x 4.0 | | 7.0 | x 10.0 | |
| | | 3.0 | x 5.0 | | approx. 25 m | 7.0 | x 11.0 |
| | | 3.5 | x 6.0 | | 8.0 | x 11.0 | |
| | | 4.0 | x 6.0 | | 8.0 | x 12.0 | |
| | | approx. 45 m | 4.0 | | x 7.0 | approx. 20 m | 9.0 |
| QTY | 1 coil | 5.0 | x 8.0 | | 9.0 | x 14.0 | |
| | approx. 30 m | 6.0 | x 9.0 | | 10.0 | x 14.0 | |
| | | 6.0 | x 9.5 | | | | |

DRAINAGE TUBING


made of soft rubber


- extra smooth
- in specific lengths
- autoclavable
- non-sterile



DRAINAGE TUBING

| REF. | ORDER SIZE LENGTH | I.D. MM | O.D. MM | SIZE LENGTH | I.D. MM | O.D. MM | | |
|---------|-------------------|---------|---------|---------------|---------|---------|--------|--------|
| 474200* | approx. 1.25 m | 1.5 | x 3.0 | approx. 1.6 m | 6.0 | x 10.0 | | |
| | | 2.0 | x 4.0 | | 7.0 | x 9.0 | | |
| | approx. 1.6 m | 3.0 | x 5.0 | | 7.0 | x 10.0 | | |
| | | 3.5 | x 5.5 | | 7.0 | x 10.5 | | |
| | | 3.5 | x 6.0 | | 7.0 | x 11.0 | | |
| | | 4.0 | x 6.0 | | 8.0 | x 11.0 | | |
| | QTY | 1 coil | 4.0 | | x 7.0 | | 8.0 | x 12.0 |
| | | 5.0 | x 7.0 | | | 9.0 | x 13.0 | |
| | | 5.0 | x 8.0 | | | 10.0 | x 14.0 | |
| | | 5.0 | x 8.5 | | | 10.0 | x 15.0 | |
| | | 6.0 | x 8.0 | | | 11.0 | x 15.0 | |
| | | 6.0 | x 9.0 | | | 12.0 | x 17.0 | |





LATEX TUBING

| REF. | ORDER SIZE LENGTH | I.D. MM | O.D. MM | SIZE LENGTH | I.D. MM | O.D. MM |
|---------------|-------------------|---------|---------|--------------|---------|---------|
| 471500* | approx. 15 m | 2.0 | x 4.0 | approx. 15 m | 5.0 | x 8.0 |
| | | 3.0 | x 5.0 | | 5.0 | x 9.0 |
| | | 3.0 | x 7.0 | | 6.0 | x 9.0 |
| | | 4.0 | x 6.0 | | 7.0 | x 10.0 |
| | | 4.0 | x 7.0 | | 8.0 | x 12.0 |
| | | 4.5 | x 10.0 | | 10.0 | x 14.0 |
| QTY 1 coil | | 5.0 | x 7.0 | | | |



LATEX TUBING

- autoclavable
- non-sterile

LATEX TUBING

| REF. | ORDER SIZE LENGTH | I.D. MM | O.D. MM | SIZE LENGTH | I.D. MM | O.D. MM |
|---------------|-------------------|---------|---------|--------------|---------|---------|
| 471501* | approx. 15 m | 2.0 | x 4.0 | approx. 15 m | 5.0 | x 8.0 |
| | | 3.0 | x 5.0 | | 5.0 | x 9.0 |
| | | 3.0 | x 7.0 | | 6.0 | x 9.0 |
| | | 4.0 | x 6.0 | | 7.0 | x 10.0 |
| | | 4.0 | x 7.0 | | 8.0 | x 12.0 |
| | | 5.0 | x 7.0 | | 10.0 | x 14.0 |
| QTY 1 coil | | | | | | |



LATEX TUBING

- autoclavable
- non-sterile

PVC TUBING

| REF. | ORDER SIZE LENGTH | I.D. MM | O.D. MM | SIZE LENGTH | I.D. MM | O.D. MM |
|---------------|-------------------|---------|---------|--------------|---------|---------|
| 475000 | approx. 250 m | 1.5 | x 2.1 | approx. 50 m | 6.0 | x 8.0 |
| | approx. 200 m | 3.5 | x 4.5 | | 7.0 | x 9.0 |
| | approx. 100 m | 4.0 | x 5.4 | | 8.0 | x 11.3 |
| | approx. 50 m | 3.0 | x 6.0 | | 11.0 | x 15.3 |
| | | | 5.0 | x 7.0 | | |
| QTY 1 coil | | | | | | |



PVC TUBING

made of PVC

- transparent
- sterilisable with ethylene oxide
- latex-free
- non-sterile

SILKOMED® VACUUM TUBING

| REF. | ORDER SIZE LENGTH | I.D. MM | O.D. MM |
|---------------|-------------------|---------|---------|
| 471800 | approx. 10 m | 7.0 | x 13.0 |
| | | 8.0 | x 14.0 |
| | | 10.0 | x 16.0 |
| QTY 1 coil | | | |



SILKOMED® VACUUM TUBING

made of silicone

- pyrogen-free
- autoclavable
- latex-free
- non-sterile

* These products contain natural rubber latex which may cause allergic reactions.

PATIENT CARE PRODUCTS

ESMARCH BANDAGE
made of soft rubber

- plain
- 0.8 mm in thickness



ESMARCH BANDAGE

| REF. | ORDER SIZE (WIDTH) CM | | | LENGTH | QTY |
|---------|-----------------------|-----|------|-------------|-----|
| 607000* | 6.0 | 8.0 | 10.0 | approx. 5 m | 2 |

MARTIN BANDAGE
made of soft rubber

- 0.6 mm in thickness



MARTIN BANDAGE

| REF. | ORDER SIZE (WIDTH) CM | | | | LENGTH | QTY |
|---------|-----------------------|-----|------|------|----------------|-----|
| 601000* | 6.0 | 8.0 | 10.0 | 12.0 | approx. 3.25 m | 2 |

KIDNEY ROLL
made of soft rubber

- with valve (neck roll)



KIDNEY ROLL

| REF. | SIZE | | QTY |
|---------|----------------------------------|--|-----|
| 815000* | 45 x 15 cm | | 1 |
| REF. | ACCESSORY | | QTY |
| 910000* | Soft rubber blower for inflation | | 1 |

INFLATABLE CUSHION
made of soft rubber

- red
- round



INFLATABLE CUSHION

| REF. | ORDER SIZE (Ø) CM | | | QTY |
|---------|----------------------------------|----|--|-----|
| 802000* | 42.5 | 45 | | 1 |
| REF. | ACCESSORY | | | QTY |
| 910000* | Soft rubber blower for inflation | | | 1 |

DOUBLE BULB
made of soft rubber

- with net



DOUBLE BULB

| REF. | ORDER SIZE | BALLOON VOLUME | QTY |
|---------|------------|----------------|-----|
| 900000* | 2 | 51 ml | 1 |
| | 3 | 60 ml | |
| | 4 | 77 ml | |

SINGLE BULB
made of soft rubber

- with cylindrical tube



SINGLE BULB

| REF. | ORDER SIZE | BALLOON VOLUME | QTY |
|---------|------------|----------------|-----|
| 901000* | 3 | 60 ml | 1 |
| | 4 | 77 ml | |

PESSARIES

PESSARY

| REF. | ORDER SIZE (Ø) MM | | | | | | QTY |
|---------|-------------------|----|----|----|-----|----|-----|
| 627000* | 50 | 55 | 60 | 65 | 70 | 75 | 10 |
| | 80 | 85 | 90 | 95 | 100 | | |



PESSARY
made of soft rubber

- red

PESSARY

| REF. | ORDER SIZE (Ø) MM | | | | | | QTY |
|--------|-------------------|----|----|----|-----|----|-----|
| 627001 | 50 | 55 | 60 | 65 | 70 | 75 | 10 |
| | 80 | 85 | 90 | 95 | 100 | | |



PESSARY
made of silicone

- blue
- latex-free

* These products contain natural rubber latex which may cause allergic reactions.

Teleflex is a leading global provider of specialty medical devices used for diagnostic and therapeutic procedures in critical care, urology and surgery. Our mission is to provide solutions that enable healthcare providers to improve outcomes and enhance patient and provider safety.

We specialise in devices for general and regional anaesthesia, cardiac care, respiratory care, urology, vascular access and surgery and we serve healthcare providers in more than 130 countries. Teleflex also provides specialty products for medical device manufacturers.

Our well known brands include ARROW®, DEKNATEL®, GIBECK®, HUDSON RCI®, KMEDIC®, PILLING®, PLEUR-EVAC®, RÜSCH®, SHERIDAN®, TAUT®, TFX OEM®, VASONOVA™ and WECK®, all of which are trademarks or registered trademarks of Teleflex Incorporated.

Teleflex global operations: Austria, Belgium, Canada, China, Czech Republic, France, Germany, Greece, India, Ireland, Italy, Japan, Mexico, Netherlands, Portugal, Singapore, Slovak Republic, South Africa, Spain, Switzerland, United Kingdom, Uruguay and USA.

YOUR CONTACTS FOR EUROPE, THE MIDDLE EAST AND AFRICA (EMEA):

TELEFLEX HEADQUARTERS EMEA, IRELAND

Teleflex Medical Europe Ltd., IDA Business Park, Athlone, Co. Westmeath
Phone +353 (0)9 06 46 08 00 · Fax +353 (0)14 37 07 73
orders.intl@teleflex.com

AUSTRIA +43 (0)1 402 47 72

BELGIUM +32 (0)2 333 24 60

CZECH REPUBLIC +420 (0)495 759 111

FRANCE +33 (0)5 62 18 79 40

GERMANY +49 (0)7151 406 0

GREECE +30 210 67 77 717

ITALY +39 0362 58 911

NETHERLANDS +31 (0)88 00 215 00

PORTUGAL +351 22 541 90 85

SLOVAK REPUBLIC +421 (0)3377 254 28

SOUTH AFRICA +27 (0)11 807 4887

SPAIN +34 918 300 451

SWITZERLAND +41 (0)31 818 40 90

UNITED KINGDOM +44 (0)1494 53 27 61

For detailed information see www.teleflex.com

The products in this catalogue are only available in the EMEA (Europe, Middle East, Africa) region. For further information contact your local representative. All data current at time of printing (08/2012). Subject to technical changes without further notice.

94 04 63 - 00 00 01 · REV B · MC/SF · 08 12 01