



WECK HEM-O-LOK LIGATION SYSTEM

Safety without compromise





THE CLIP BANK

EASY LOADING WITH CONVENIENT ADHESIVE STRIP

- Our Hem-o-lok ligation system combines speed and safety in one single procedure.
- Choose from a variety of ligating clip sizes for the ligation of vessels and/or tissues between 2 and 16 mm.
- The Hem-o-lok clip bank allows for effortless and precise loading, thus making surgical interventions more efficient.
- Closure of the clip is confirmed by tactile and audible feedback.
- The colour code of the clip bank corresponds to that of the applier.
- Each cartridge contains 6 clips.
- A special adhesive strip at the bottom of each cartridge enables you to choose the best possible site for clip delivery.
- Unlike absorbable clips, Hem-o-lok ligating clips can be removed with our specially developed removers if necessary.
- The clips are inert, non-toxic and radiolucent.
- The cartridges adhere to any sterile surface.



THE APPLIERS

PERFORMANCE FEATURES AND CLEANING OPTIONS

- Thanks to the special clip design, closure of the clip can be clearly felt by the surgeon during application.
- Both curved and angled appliers are available for open surgery.
- Endoscopic Hem-o-lok appliers are available in sizes ML, L, and XL. They can be taken apart for easy cleaning, thus eliminating potential contamination.
- A 360° rotation knob facilitates access and improves visibility.
- The appliers are colour-coded to match the colours of the corresponding cartridges.

THE HEM-O-LOK LIGATION SYSTEM – PRECISION IN EVERY DETAIL

Teleflex is a leading manufacturer of metal and polymer clips. Our Weck Hem-o-lok ligation system forms a line of instruments perfectly adapted to each other, meeting all the special needs of surgeons, nurses and operating theatre personnel alike.



THE CLIP

TOP-CLASS FUNCTIONALITY

| CLIP SIZE | LIGATION OF VESSELS AND TISSUES |
|----------------|---------------------------------|
| ● medium | 2 – 7 mm |
| ● medium-large | 3 – 10 mm |
| ● large | 5 – 13 mm |
| ● extra-large | 7 – 16 mm |

- Precision bosses engineered for maximum security.
- Penetrating lock: 1 x 2 teeth to provide effective closure.
- Bowed design enables added tissue to be encompassed.
- Tactile and audible response confirms that the Hem-o-lok ligating clip has locked onto vessel.
- Integrated teeth help prevent slippage in any direction.
- The surface properties of the non-absorbable polymer prevent infection.
- Living hinge provides flexibility when opening and closing.

LIGATING CLIPS



OPEN SURGERY



All Hem-o-lok ligating clips are sterile packaged. Ligating clip systems differ in their clip-closure characteristics. It is the responsibility of the user to confirm security of clips on vessels or other tissue. Not intended for contraceptive tubal occlusion.

HEM-O-LOK LIGATING CLIPS

WECK

| REF. | CLIP SIZE | LIGATION OF VESSELS AND TISSUE | CLIPS/CARTRIDGE | CARTRIDGES/BOX | CLIPS/BOX |
|----------|-----------|--------------------------------|-----------------|----------------|-----------|
| WK544220 | ● M | 2 – 7 mm | 6 | 14 | 84 |
| WK544230 | ● ML | 3 – 10 mm | 6 | 14 | 84 |
| WK544240 | ● L | 5 – 13 mm | 6 | 14 | 84 |
| WK544250 | ● XL | 7 – 16 mm | 6 | 14 | 84 |

APPLIERS

WECK

| REF. | CLIP SIZE | WORKING LENGTH | QTY |
|----------|-----------|----------------|-----|
| curved | | | |
| WK544113 | ● M | 20 cm | 1 |
| WK544170 | ● ML | 20 cm | 1 |
| WK544180 | ● L | 20 cm | 1 |

APPLIERS, LONG

WECK

| REF. | CLIP SIZE | WORKING LENGTH | QTY |
|----------|-----------|----------------|-----|
| curved | | | |
| WK544115 | ● M | 27 cm | 1 |
| WK544171 | ● ML | 27 cm | 1 |
| WK544181 | ● L | 27 cm | 1 |
| WK544191 | ● XL | 27 cm | 1 |

APPLIERS, LONG, 70° ANGLED

WECK

| REF. | CLIP SIZE | WORKING LENGTH | QTY |
|----------|-----------|----------------|-----|
| WK544114 | ● M | 27 cm | 1 |
| WK544172 | ● ML | 27 cm | 1 |
| WK544179 | ● L | 27 cm | 1 |
| WK544192 | ● XL | 27 cm | 1 |

CLIP REMOVERS

WECK

| REF. | CLIP SIZE | WORKING LENGTH | QTY |
|----------|--------------------|----------------|-----|
| WK544123 | ● M ● ML ● L | 27 cm | 1 |
| WK544124 | ● L ● XL | 27 cm | 1 |

● M medium ● ML medium-large ● L large ● XL extra-large

ENDOSCOPIC SURGERY



HEM-O-LOK TAKE-APART APPLIERS

WECK

| REF. | CLIP SIZE | Ø APPLIER | Ø TROCAR | WORKING LENGTH | QTY |
|-----------|-----------|-----------|-------------|----------------|-----|
| straight | | | | | |
| WK544965T | ● ML | 5 mm | 5 mm | 33 cm | 1 |
| curved | | | | | |
| WK544995T | ● L | 10 mm | 10 mm | 33 cm | 1 |
| WK544990T | ● XL | 10 mm | 10 or 12 mm | 33 cm | 1 |

HEM-O-LOK TAKE-APART APPLIERS, 20° ANGLED

WECK

| REF. | CLIP SIZE | Ø APPLIER | Ø TROCAR | Ø SILS | WORKING LENGTH | QTY |
|-------------|-----------|-----------|----------|--------|----------------|-----|
| WK544965T20 | ● ML | 5 mm | 10 mm | 5 mm | 33 cm | 1 |
| WK544995T20 | ● L | 10 mm | 12 mm | 10 mm | 33 cm | 1 |
| WK544990T20 | ● XL | 10 mm | 15 mm | 12 mm | 33 cm | 1 |

HEM-O-LOK NON-TAKE-APART APPLIER, LONG

WECK

| REF. | CLIP SIZE | Ø APPLIER | Ø TROCAR | WORKING LENGTH | QTY |
|-----------|-----------|-----------|-------------|----------------|-----|
| straight | | | | | |
| WK544965L | ● ML | 5 mm | 5 mm | 45 cm | 1 |
| curved | | | | | |
| WK544995L | ● L | 10 mm | 10 mm | 45 cm | 1 |
| WK544990L | ● XL | 10 mm | 10 or 12 mm | 45 cm | 1 |

HEM-O-LOK AUTO ENDO 5 APPLIERS

WECK

| REF. | CLIP SIZE | Ø APPLIER | Ø TROCAR | WORKING LENGTH | QTY |
|----------|-----------|-----------|----------|----------------|-----|
| straight | | | | | |
| WK543965 | ● ML | 5 mm | 5 mm | 33 cm | 3 |

CLIP REMOVERS

WECK

| REF. | CLIP SIZE | Ø REMOVER | Ø TROCAR | Ø SILS | WORKING LENGTH | QTY |
|-------------|-----------|-----------|-------------|--------|----------------|-----|
| straight | | | | | | |
| WK544121T | ● ML | 5 mm | 5 mm | - | 33 cm | 1 |
| WK544121L | ● ML | 5 mm | 5 mm | - | 45 cm | 1 |
| curved | | | | | | |
| WK544130T | ● L/XL | 10 mm | 10 or 12 mm | - | 33 cm | 1 |
| WK544130L | ● L/XL | 10 mm | 10 or 12 mm | - | 45 cm | 1 |
| 20° angled | | | | | | |
| WK544130T20 | ● L/XL | 10 mm | 12 or 15 mm | 12 mm | 33 cm | 1 |



ENDOSCOPIC TAKE-APART APPLIER

Our take-apart endoscopic applicators can be detached for easy cleaning. The 360° knob facilitates access and improves visibility.



Teleflex is a leading global provider of specialty medical devices used for diagnostic and therapeutic procedures in critical care, urology and surgery. Our mission is to provide solutions that enable healthcare providers to improve outcomes and enhance patient and provider safety.

We specialise in devices for general and regional anaesthesia, cardiac care, respiratory care, urology, vascular access and surgery and we serve healthcare providers in more than 140 countries. Teleflex also provides specialty products for medical device manufacturers.

Our well known brands include ARROW®, DEKNATEL®, GIBECK®, HUDSON RCI®, KMEDIC®, LMA™, PILLING®, PLEUR-EVAC®, RÜSCH®, SHERIDAN®, TAUT®, TFX OEM®, VASONOVA™ and WECK®, all of which are trademarks or registered trademarks of Teleflex Incorporated.

Teleflex global operations: Australia, Austria, Belgium, Canada, China, Czech Republic, France, Germany, Greece, India, Ireland, Italy, Japan, Malaysia, Mexico, Netherlands, New Zealand, Portugal, Singapore, Slovak Republic, South Africa, Spain, Switzerland, United Kingdom, Uruguay and USA.

YOUR INTERNATIONAL CONTACTS:

TELEFLEX HEADQUARTER INTERNATIONAL, IRELAND

Teleflex Medical Europe Ltd., IDA Business and Technology Park,
Dublin Road, Athlone, Co Westmeath
Phone +353 (0)9 06 46 08 00 · Fax +353 (0)14 37 07 73
orders.intl@teleflex.com

AUSTRALIA/NEW ZEALAND +61 (0)3 84 12 03 21

AUSTRIA +43 (0)1 402 47 72

BELGIUM +32 (0)2 333 24 60

CHINA (SHANGHAI) +86 (0)21 6381 6399

CHINA (BEIJING) +86 (0)10 6418 5699

CZECH REPUBLIC +420 (0)495 759 111

FRANCE +33 (0)5 62 18 79 40

GERMANY +49 (0)7151 406 0

GREECE +30 210 67 77 717

INDIA +91 (0)44-2836 5040

ITALY +39 0362 58 911

JAPAN +81 (0)3 3379 1511

NETHERLANDS +31 (0)88 00 215 00

PORTUGAL +351 22 541 90 85

SINGAPORE +65 6447 5022

SLOVAK REPUBLIC +421 (0)3377 254 28

SOUTH AFRICA +27 (0)11 807 4887

SPAIN +34 918 300 451

SWITZERLAND +41 (0)31 818 40 90

UNITED KINGDOM +44 (0)1494 53 27 61

For detailed information see www.teleflex.com

The products in this catalogue may not be available in all countries.
Please contact your local representative. All data current at time of printing (05/2013).
Subject to technical changes without further notice.

94 05 70 - 00 00 01 · REV A · MC / WM · 05 13 02