ETHICON

Coated VICRYL* *rapide*

(Polyglactin 910) Braided Suture



Coated **VICRYL*** rapide (Polyglactin 910) Braided Suture

Braided

For easy handling and secure knot tying.

Coated

For smooth passage through tissue and easy knot tie down.

Synthetic

For minimal tissue reaction.

Absorbable

Rapid absorption by hydrolysis. Complete by 42 days. Always predictable and reliable.

Tensile Strength

Post implantation Approximate % original strength remaining 5 days 50%

All of the original tensile strength is lost by approximately 10-14 days post implantation.

VICRYL* rapide Suture typically falls off 7-10 days post-operative or can be wiped off subsequently with sterile gauze. Normally the removal of the suture is not required.

Colour

Violet or undyed.

Range

Gauge sizes 8-0 to 1 (U.S.P.). Supplied as needled sutures and ligatures.

Sterilisation

Undyed by Irradiation. Violet by Ethylene Oxide.

These sutures comply with the requirements of U.S. Pharmacopoeia (U.S.P.) and The European Pharmacopoeia (PhEUR).

Coated VICRYL**rapide* (Polyglactin 910) Braided Suture

Needle length and descrip	tion	Length	Colour				Gauge			
⊖ 1/2 circle Taper Point plus				5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4	2 5
RB-1 plus	17mm	75cm	undyed	W9969	W9970					
SH-1 plus	22mm	75cm	undyed		W9971	W9974				
MH-1 plus	31mm	75cm	undyed				W9975			
СТ	40mm	90cm	undyed					W9977		
СТХ	48mm	90cm	undyed				W9978	W9979		

1/2 circle TAPERCUT*				5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4	2 5
V-7	26mm	90cm	undyed				V1002H			
V-34	36mm	90cm 120cm	undyed undyed				W9962 W9948	W9963	W9964	
V-37	40mm	90cm	undyed				W9946			
V-39	45mm	90cm	undyed				W9945	W9947		

▼ 1/2 circle Reverse Cutting				5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4	2 5
CP-2 26	mm	75cm	undyed				W9936			

Coated VICRYL* *rapide* (Polyglactin 910) Braided Suture

Needle length and description	Length	Colour			Gauge		
▼ 3/8 circle Reverse Cutting MultiPass* PRIME* Needle			6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
P-1 11mn	45cm 45cm 70cm	violet undyed undyed	MPVR32H MPVR9913H MPVR489H	MPVR9915H			
P-3 13mn	45cm	undyed		MPVR4930H	MPVR4940H		
PS-2 19mn	45cm 70cm	undyed undyed			MPVR4960H MPVR4260H	MPVR4970H MPVR4270H	
PS-1 24mn	45cm 70cm	undyed undyed				MPVR6830H MPVR9360H	MPVR9380H
PS 26mn	70cm	undyed				MPVR9932H	MPVR9933H

▲ 3/8 circle Conventional Cutting MultiPass* PRIME* Needle				6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
PC-3	16mm	70cm	undyed			MPVR9918H		

Coated VICRYL* *rapide* (Polyglactin 910) Braided Suture

Needle length and description	Length	Colour				Gauge			
▼ 3/8 circle Reverse Cutting PRIME* Needle			6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4
P-1 11mm	45cm	undyed	W9913	W9915					
PS-2 19mm	75cm	undyed			W9922	W9923			
PS 26mm	75cm	undyed				W9932	W9933		
PSLX 36mm	75cm	undyed				W9937	W9938		

▼ 3/8 circle Reverse Cutting				6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4
C-3	13mm	75cm	undyed			VR2290				
FS-3	16mm	75cm	undyed		VR2296					
FS-2	19mm	45cm 75cm	undyed undyed		V2910H	W9930	W9931	VR2299		
FS	26mm	75cm	undyed				W9940	W9941		

▲ 3/8 circle Conventional Cutting PRIME* Needle			6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4
PC-3 16mm	75cm	undyed			W9918	W9919			

▲ 3/8 circle Conventional Cutting			6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4
FS-3 CONV 16mm	75cm	undyed			W9924	W9925			

▲ 1/2 circle Conventional Cutting PRIME* Needle			6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4
PC-22 22mm	45cm	undyed			W9926	W9927			

▲ 1/2 circle Conventional Cutting			6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4
X-1 CONV 22mm	45cm	undyed				W9935			

5/8 circle Conventional Cutting			6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4
UX-25 26mm	75cm	undyed				W9929			

Coated VICRYL* *rapide* (Polyglactin 910) Braided Suture

Needle length and description			Colour	Gauge				
Compound Curved Reverse Cutting PRIME* Needle				6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
PS-4C	16mm	45cm	undyed		W9910	W9911		

▼ Straight Cutting				6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
<u>SC-1</u>	13mm	45cm	undyed			W9909		
KS-45	45mm	45cm	undyed			W9939		
KS	60mm	75cm	undyed				W9943	W9944

			6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
UV-17 17mm	75cm	undyed			W9928		

Needle length and description			Colour	Gauge				
- 3/8 circle MICRO-POINT spatula				8-0 0.4	7-0 0.5	6-0 0.7	5-0 1	4-0 1.5
GS-9	6.5mm	30cm	violet	V1210G		EC1534G		
GS-9	6.5mm	30cm 45cm	violet violet	V39G	V38G			
TG140-8	6.5mm	30cm	violet		W9914			

● 3/8 circle Taper Point				8-0 0.4	7-0 0.5	6-0 0.7	5-0 1	4-0 1.5
BV-1 9	9.3mm	30cm	violet	W9916				

Coated VICRYL**rapide* (Polyglactin 910) Braided Suture

Ligatures and Non-needled Sutures								
Ga	uge	undyed 1.5m lengths 12 sachets per box						
3-0	2	W9006						
2-0	3	W9007						