


ETHICON

MONOCRYL* Plus

Antibacterial (Poliglecaprone 25)

Monofilament Suture



MONOCRYL* Plus

Antibacterial (Poliglecaprone 25)

Monofilament Suture

All the benefits of MONOCRYL* (Poliglecaprone 25) Monofilament Suture: pliability, smoothness and ease of passage through tissue, with the additional benefit of an antibacterial.

Ideal for -

- Skin closure
- Subcutaneous fat
- Peritoneum, stomach and bowel
- Vagina and uterus
- Bladder and ureter

MONOCRYL* PLUS Sutures are effective against: ^{4, 7}

<i>Staphylococcus aureus</i>	<i>Staphylococcus epidermidis</i>
Methicillin-resistant <i>Staphylococcus aureus</i> (MRSA)	Methicillin-resistant <i>Staphylococcus epidermidis</i> (MRSE)
<i>Escherichia coli</i>	<i>Klebsiella pneumoniae</i>

References:

- 4 Rothenburger S, Spangler D, Bhende S, Burkley D. In vitro antibacterial evaluation of Coated VICRYL* Plus Antibacterial Suture (Coated Polyglactin 910 with Triclosan) using zone of inhibition assays. Surg Infect. 2002;3 (suppl 1):S79-S87.
- 7 Ming, X Rothenburger S, and Dachuan Y. In vitro antibacterial efficiency of MONOCRYL* Plus Antibacterial Suture (Poliglecaprone 25 with Triclosan). Surgical Infections Vol. 8, No. 2: 2007.

STRENGTH RETENTION

Violet suture 60 % AT 7 DAYS 30% AT 14 DAYS	Undyed suture 50 % AT 7 DAYS 20% AT 14 DAYS
-----------------------------------------------------------------	-----------------------------------------------------------------










Absorption completed at 90 - 120 days

Providing -

- **Smoothness** of passage through tissue
- Excellent **cosmetic** results
- Superior **strength**
- Additional **protection**








MONOCRYL* Plus Antibacterial (Poliglecaprone 25) Monofilament Suture

Needle length and description	Length	Colour	Gauge				
⊙ 3/8 circle Taper Point			4-0 1.5	3-0 2	2-0 3	0 3.5	1 4
 BB 17mm	70cm	violet	MCP3548H				

Needle length and description	Length	Colour	Gauge						
⊖ 1/2 circle Taper Point plus			6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4
 RB-2 13mm	45cm	violet	MCP3224G						
 RB-1 17mm	70cm	violet			MCP3435G	MCP215H			
 SH-1 plus 22mm	70cm	violet			MCP218H	MCP219G	MCP220H		
 JB-1 VISI-BLACK* 22mm	70cm	violet				MCP3110H			
 JB VISI-BLACK* 26mm	70cm	violet				MCP3160H	MCP3170H		
 SH plus 26mm	70cm	violet				MCP4160H	MCP4170H		
 SH plus 26mm	70cm	violet				MCP3627G	MCP3628G		
 MH-1 plus 31mm	70cm	violet				MCP3200H	MCP3441G	MCP247H	MCP3755G
 MH plus 36mm	70cm	violet					MCP227H	MCP228H	MCP229H







MONOCRYL* Plus




Antibacterial (Poliglecaprone 25) Monofilament Suture


Needle length and description	Length	Colour	Gauge							
☉ 1/2 circle Taper Point (heavier wire diameters)			6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4	
 CT	40mm	70cm	violet						MCP3758H	MCP3759G
☉ 5/8 circle Taper Point			6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4	
 UR-6	26mm	70cm	violet			MCP3625H	MCP6050H			
☑ 1/2 circle TAPERCUT*			6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4	
 V-7	26mm	70cm	violet				MCP3440G			
 V-34	36mm	70cm 90cm	violet violet				MCP3488G	MCP3489G MCP918H	MCP3490G	
▼ 1/2 circle Reverse Cutting			6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4	
 CPA	30mm	70cm	undyed				MCP1086H			
 CP-1	36mm	70cm	violet				MCP3600H	MCP4267G		
 CP	40mm	90cm	violet						MCP3770G	

MONOCRYL* Plus

Antibacterial (Poliglecaprone 25) Monofilament Suture

Needle length and description	Length	Colour	Gauge				
▼ 3/8 circle Reverse Cutting - MultiPass* PRIME* Needle			6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
P-1 	11mm	70cm	undyed	MCP489H	MCP490H		
P-3 	13mm	70cm	undyed	MCP492H	MCP493H	MCP494H	
PS-3 	16mm	70cm	undyed		MCP500H	MCP3205G	
PS-2 	19mm	70cm 45cm 70cm	violet undyed undyed		MCP495H	MCP498H MCP496H MCP4260H	MCP497H MCP4271H
PS-1 	24mm	45cm 70cm	undyed undyed			MCP6820H MCP9350H	MCP683H MCP684H
PS 	26mm	70cm	undyed			MCP3212H	MCP3213H

▼ 3/8 circle Reverse Cutting			6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
FS-3 	16mm	45cm	undyed		MCP3221G		
FS-2 	19mm	45cm 70cm	undyed undyed		MCP3209G	MCP292H	MCP4423H
FS 	26mm	70cm	undyed			MCP3326G	MCP3327G

▼ Straight Cutting			6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
KS 	60mm	70cm	undyed			MCP3649G	MCP3650G MCP3651G

