

# **BD** Medical Medical Surgical Systems Catalogue Canadian Version



Helping all people live healthy lives

Needles	
and	
Syringes	

Infusion Therapy Devices

Surgical Products Anesthesia

Sharps Disposal Products

Products

## **Table of Contents**

Click the section you'd like to jump to

### Needles and Syringes

#### 3-16

Needles and Synnyes	5-10
Safety needles and Syringes	3-12
<ul> <li>BD SafetyGlide needles and syringes</li> </ul>	3
<ul> <li>BD SafetyGlide needles and syringes for Insulin, TB and Allergy</li> </ul>	5
<ul> <li>BD Integra needles and syringes</li> </ul>	7
<ul> <li>BD Eclipse needles and syringes</li> </ul>	9
<ul> <li>BD Safety-Lok syringes</li> </ul>	11
<ul> <li>BD<sup>™</sup> Syringes and Needles</li> </ul>	13
Blunt Fill Needles	13
Pharmacy Products	14
Diabetes Products	15
<ul> <li>BD<sup>™</sup> Reusable Products</li> </ul>	16
Infusion Therapy Devices	17-26
Safety IV Catheters	18
BD Nexiva Closed IV Catheter System	18
BD Insyte AutoGuard Shielded IV Catheters	20
• BD <sup>™</sup> IV Catheters	22
<ul> <li>BD PosiFlush<sup>™</sup> Prefilled Saline Syringes</li> </ul>	23
• BD Q-Syte™	25
<ul> <li>Interlink<sup>®</sup> System</li> </ul>	26
Other Infusion Therapy Devices	27
Surgical Products	28-32
<ul> <li>BD Bard-Parker<sup>™</sup> Protected Blade System</li> </ul>	29
<ul> <li>BD Bard-Parker<sup>™</sup> Protected Disposable Scalpel</li> </ul>	30
• BD Bard-Parker <sup>™</sup> Blades and Scalpels	31
Antiseptic Surgical Scrubs	32
Heimlich Chest Drain Valve	32
Anesthesia Products	33-43

BD Spinal and Epidural Needles

## Sharps Disposal Products

44-50

34

## BD SafetyGlide<sup>™</sup>Needles and Syringes

# Safety Without Compromise

Setting the standard in safety for a wide variety of applications.

## The most versatile safety injection technology available today.

- Features single-handed activation that allows for ease of use and minimal training.
- Low profile safety mechanism remains discretely out of the way before and during injection and requires minimal change in technique.
- Activation-Assist<sup>™</sup> technology minimizes finger motion required for activation and keeps fingers behind the needle tip.
- Stainless-steel housing securely shields the needle tip after activation.
- Bevel-up design facilitates ease of injection. Safety-arm is positioned for optimum visualization of injection site.
- Available in a full array of sizes.
- Compatible with any size conventional syringe.

# **BD SafetyGlide<sup>™</sup>Needles and Syringes**

#### **Prior to activation**



#### After activation

Stainless steel latch securely shields needle tip after activation.



#### **BD SafetyGlide<sup>™</sup> Needles**

PROD #	DESCRIPTION	PACKAGING
305921	27 GA x 5/8 in. Needle	50/box, 500/case
305901	25 GA x 5/8 in. Needle	50/box, 500/case
305916	25 GA x 1 in. Needle	50/box, 500/case
305902	23 GA x 1 in. Needle	50/box, 500/case
305900	22 GA x 1 1/2 in. Needle	50/box, 500/case
305915	21 GA x 1 in. Needle	50/box, 500/case
305917	21 GA x 1 1/2 in. Needle	50/box, 500/case
305918	18 GA x 1 1/2 in. Needle	50/box, 400/case

#### 1 mL Syringes with BD SafetyGlide<sup>™</sup> Needles

PROD #	DESCRIPTION	PACKAGING
305927	27 GA x 5/8 in. Needle	50/box, 400/case
305903	25 GA x 5/8 in. Needle	50/box, 400/case

#### 3 mL Syringes with BD SafetyGlide<sup>™</sup> Needles

PROD #	DESCRIPTION	PACKAGING
305924	25 GA x 1 in. Needle	50/Box, 400/Case
305904	25 GA x 5/8 in. Needle	50/Box, 400/Case
305905	23 GA x 1 in. Needle	50/Box, 400/Case
305906	22 GA x 1 1/2 in. Needle	50/Box, 400/Case
305909	21 GA x 1 1/2 in. Needle	50/Box, 400/Case

#### 5 mL Syringes with BD SafetyGlide<sup>™</sup> Needles

PROD #	DESCRIPTION	PACKAGING
305907	22 GA x 1 1/2 in.	50 Box, 400/Case

#### 10 mL Syringes with BD SafetyGlide<sup>™</sup> Needles

PROD #	DESCRIPTION	PACKAGING
305908	22 GA x 1 1/2 in.	50 Box, 400/Case

# **BD SafetyGlide<sup>™</sup> Syringe**

# Safety In Small Doses

Completes the SafetyGlide line to meet your Subcutaneous and Intradermal injection needs.


## Available in three small syringe sizes–1 mL, 1/2 mL and 3/10 mL.

- 26 GA to 31 GA needle sizes perform a variety of subcutaneous and intradermal applications.
- 1 mL, 1/2 mL and 3/10 mL insulin syringe sizes facilitate small dose injections.
- Permanently attached needle means minimal waste space:
  - A necessity for insulin injection.
  - A potential medication cost savings if used instead of a detachable needle and syringe.
- Features single-handed activation that allows for ease of use and minimal training.
- Rotating safety arm eliminates interference with scale during drawing up, and facilitates low-angle injections.
- Cap is designed to stand on end for true one-handed passive recapping of clean needles for transport.

# **BD SafetyGlide<sup>™</sup> Syringe**

for Insulin, TB or Allergy

#### **Prior to activation**



#### BD SafetyGlide Needles with Syringes for Insulin, TB or Allergy

PROD #	DESCRIPTION	PACKAGING
305930	1 mL Insulin Syringe with 29 GA x 1/2 in. Needle	100/Box 400/Case
305932	1/2 mL Insulin Syringe with 29 GA x 1/2 in. Needle	100/Box 400/Case
305934	1/2 mL Insulin Syringe with 30 GA x 5/16 in. Needle - Thinwall	100/Box 400/Case
305935	3/10 mL Insulin Syringe with 29 GA x 1/2 in. Needle	100/Box 400/Case
305937	3/10 mL Insulin Syringe with 31 GA x 5/16 in. Needle - Thinwa	ll 100/Box 400/Case
305945	1 mL Tuberculin Syringe with 27 GA x 1/2 in. Needle (RB)	100/Box 400/Case
305946	1 mL Tuberculin Syringe with 26 GA x 3/8 in. Needle (ID)	100/Box 400/Case
305950	Allergist Tray: 1 mL Syringe with 27 GA x 1/2 in. Needle (RB)	25/Tray 1000/Case
305951	Allergist Tray: 1 mL Syringe with 26 GA x 3/8 in. Needle (ID)	25/Tray 1000/Case

(RB) = Regular bevel (ID) = Intradermal bevel

# **BD Integra**<sup>™</sup>

# Shaping the Future of Safety

The latest innovation in retracting technology.

## The only retracting syringe that truly meets your clinical needs.

- It's the first spring-based retracting syringe with a detachable needle! Allows the clinician to change needles for aspiration and administration.
- The BD Integra<sup>™</sup> Syringe has extremely low waste space. Low waste space reduces medication waste and cost.<sup>1</sup>



#### Red area = waste space

The BD Integra Syringe has up to 86% less waste space than other manufacturers' syringe designs.

- The clinician decides when to activate the retracting mechanism either before or after needle is withdrawn from patient.
- The low waste space design means there is minimal splatter after activation.
- Available in a wide range of needle gauges for use with 3mL syringes for Intramusular, Subcutaneous and Intradermal injections.

#### Visit www.bd.com/Integra for more information.

1 Data on file. Texas childrens white paper: Recent studies show that Integra can obtain more doses from multi-dose vials of flu vaccine

Customer Service: 1.866.979.9408

BD-Canada

**BD** Integra<sup>™</sup> With Retracting BD PrecisionGlide<sup>™</sup> Needle **Prior to activation BD** PrecisionGlide<sup>™</sup> Clear Barrel **Needle Technology** Wide Flanges For dosing accuracy "World Class Sharpness" For comfort and stability **Single-Handed Activation Detachable Needle** Tru-Lok Technology Ensures low waste space and a reliable connection

**After activation** 



#### Needle Retracts into the Syringe!

- Eliminates exposure to the needle after activation

## BD Integra<sup>™</sup> 3 mL Syringe with Detachable Needle

The BD Integra<sup>™</sup> Syringe with detachable needle provides the benefit of low waste space, which is important in reducing medication waste and costs.

PROD #	DESCRIPTION	PACKAGING
305269	3 mL Syringe with Detachable 25 GA x 5/8 in. Needle	100/box, 400/case
305270	3 mL Syringe with Detachable 25 GA x 1 in. Needle	100/box, 400/case
305271	3 mL Syringe with Detachable 23 GA x 1 in. Needle	100/box, 400/case
305272	3 mL Syringe with Detachable 22 GA x 1 1/2 in. Needle	100/box, 400/case
305273	3 mL Syringe with Detachable 21 GA x 1 in. Needle	100/box, 400/case
305274	3 mL Syringe with Detachable 21 GA x 1 1/2 in. Needle	100/box, 400/case

#### **BD** Integra<sup>™</sup> Needles

PROD #	DESCRIPTION	PACKAGING
305310	25 GA x 5/8 in. Retracting Needle	100/box, 1000/case
305311	25 GA x 1 in. Retracting Needle	100/box, 1000/case
305312	23 GA x 1 in. Retracting Needle	100/box, 1000/case
305313	22 GA x 1 1/2 in. Retracting Needle	100/box, 1000/case
305275 <sup>1</sup>	18 GA x 1 1/2 in. Blunt Fill Needle	100/box, 1000/case
305833 <sup>1</sup>	18 GA x 1 1/2 in. Blunt Fill Needle with Filter	50/box, 500/case

1 – 18 GA Blunt Fill and Blunt Filter Needles have nonretracting needles.

#### BD Integra<sup>™</sup> 3 mL Syringe Only

PROD #	DESCRIPTION	PACKAGING
305283	3 mL Syringe Only	100/box, 400/case

# **BD Eclipse<sup>™</sup> Injection Needle**

# Safety Made Simple


Intuitive and easy to use a highly economical approach to safety.

## The same technology for both skin injection and blood collection.

- Activates with a single finger stroke immediately after injection–no hard surface is required.
- Bevel-up design facilitates ease of injection. Safety-arm is positioned for optimum visualization of injection site.
- Wide, textured finger pad provides safe activation area for thumb or forefinger.

# **BD Eclipse**<sup>™</sup>

#### **Prior to activation**



Needle remains locked inside activated cover

#### **BD Eclipse<sup>™</sup> Needles**

With this single-handed, needle-based safety-engineered injection device, BD offers one technology for both injection and blood-collection applications. The needle is bevel-up oriented to the safety cover to facilitate low-angle injections.

PROD #	DESCRIPTION	PACKAGING
305757	30 GA x 1/2 in. Needle	100/box, 1200/case
305758	27 GA x 1/2 in. Needle	100/box, 1200/case
305767	25 GA x 1 1/2 in. Needle	100/box, 1200/case
305761	25 GA x 1 in. Needle	100/box, 1200/case
305759	25 GA x 5/8 in. Needle	100/box, 1200/case
305762	23 GA x 1 in. Needle	100/box, 1200/case
305763	22 GA x 1 1/2 in Needle	100/box, 1200/case
305768	22 GA x 1 in. Needle	100/box, 1200/case
305765	21 GA x 1 1/2 in. Needle	100/box, 1200/case
305764	21 GA x 1 in. Needle	100/box, 1200/case
305766	18 GA x 1 1/2 in. Needle	100/box, 1200/case

#### **BD Eclipse<sup>™</sup> Needles with SmartSlip<sup>™</sup> Technology\***

PROD #	DESCRIPTION	PACKAGING
305771	30 G x 1/2 in.	100/box, 1200/case
305770	27 G x 1/2 in.	100/box, 1200/case
305760	25 G x 5/8 in.	100/box, 1200/case

## 1 mL BD Luer-Lok<sup>™</sup> Syringe with Detachable Needle

PROD #	DESCRIPTION	PACKAGING
305778	30 GA x 1/2 in.	50/box, 300/case
305789	27 GA x 1/2 in.	50/box, 300/case
305780	25 GA x 5/8 in.	50/box, 300/case

\*Designed for use with Luer Slip Syringes

## 3 mL BD Luer-Lok<sup>™</sup> Syringe with Detachable Needle

PROD #	DESCRIPTION	PACKAGING
305781	25 GA x 5/8 in.	50/box, 300/case
305787	25 GA x 1 in.	50/box, 300/case
305782	23 GA x 1 in.	50/box, 300/case
305788	22 GA x 1 in	50/box, 300/case
305783	22 GA x 1 1/2 in.	50/box, 300/case
305779	21 GA x 1 in.	50/box, 300/case
305784	21 GA x 1 1/2 in.	50/box, 300/case

#### 5 mL BD Luer-Lok<sup>™</sup> Syringe with Detachable Needle

PROD #	DESCRIPTION	PACKAGING
305785	22 GA x 1 1/2 in.	50/box, 300/case

#### 10 mL BD Luer-Lok<sup>™</sup> Syringe with Detachable Needle

PROD #	DESCRIPTION	PACKAGING
305786	22 GA x 1 1/2 in.	50/box, 300/case

# **BD Safety-Lok<sup>™</sup> Syringe**



# **First In Safety**

The world's first safety-engineered syringe.

### Proven performance and dependability.

- The only device of its kind to provide color-coded visual confirmation of the safety lock.
- Protective sleeve locks into position with one motion to cover the needle after injection.
- Crystal-clear barrel and sleeve allows for dosing visualization.
- Large flanges and tapered plunger rod increase control for easier aspiration and injection.
- Interim "park" position of safety shield complements injection practice in the OR.
- Can be used with any needle 1-1/2 in. or shorter providing safety for many sizes that do not have other alternatives.
- Available in a wide array of sizes-1 mL, 3 mL, 5 mL and 10 mL.

# **BD Safety-Lok<sup>™</sup> Syringe**

### **Prior to activation**

Large Flanges and Tapered Plunger Rod increase control for easier aspiration and injection Crystal Clear Barrel and Sleeve make syringe contents easier to see

Bold Scale Markings facilitate accurate dose measurement

World Class Sharpness for patient comfort

#### After activation

Visual Lock Indicator Lock confirmed when green band fully covers red collar (Also hear and feel the "lock")

#### **BD Safety-Lok<sup>™</sup> Syringe**

The only safety-engineered syringe of its type to offer clear visual confirmation of the lock position.

PROD #	DESCRIPTION	PACKAGING
329464	1 mL Insulin Syringe with Permanently Attached 29 GA x 1/2 in. BD Ultra-Fine™ Needle (U-100)	100/box, 500/case
305553	1 mL Tuberculin Syringe with Permanently Attached 27 GA x 1/2 in. Needle	100/box, 500/case
305554	1 mL Tuberculin Syringe with Permanently Attached 25 GA x 5/8 in. Needle	100/box, 500/case
309606	3 mL Syringe only	100/box, 800/case
309592	3 mL Syringe with 25 GA x 5/8 in. Needle	100/box, 800/case
309594	3 mL Syringe with 23 GA x 1 in. Needle	100/box, 800/case
309593	3 mL Syringe with 22 GA x 1 1/2 in. Needle	100/box, 800/case
309596	3 mL Syringe with 22 GA x 1 in. Needle	100/box, 800/case
309595	3 mL Syringe with 21 GA x 1 1/2 in. Needle	100/box, 800/case
305558	5 mL Syringe only	50/box, 400/case
305561	5 mL Syringe with 21 GA x 1 1/2 in. Needle	50/box, 400/case
305559	10 mL Syringe only	50/box, 400/case
305564	10 mL Syringe with 21 GA x 1 1/2 in. Needle	50/box, 400/case

## **BD<sup>™</sup> Syringes and Needles**

#### 1 mL Tuberculin Syringe

PROD #	DESCRIPTION	PACKAGING
309602	TB Syringe Without Needle, Slip Tip	100/shelf pack, 800/case

#### 1 mL Luer-Lok<sup>™</sup> Syringe Only

PROD #	DESCRIPTION	PACKAGING
309628	BD Luer-Lok™ 1/100 mL	100/box, 800/case
309629	U-100 Insulin	100/box, 800/case

#### 3 mL Syringe Only

PROD #	DESCRIPTION	PACKAGING
309586	3 mL Syringe, Slip Tip	100/shelf pack, 800/case
309585	3 mL Syringe, BD Luer-Lok™ Tip	100/shelf pack, 800/case

#### 5 mL Syringe Only

PROD #	DESCRIPTION	PACKAGING
301603	5 mL Syringe, Slip Tip	100/shelf pack, 400/case
309603	5 mL Syringe, BD Luer-Lok™ Tip	100/shelf pack, 400/case

#### 10 mL Syringe Only

PROD #	DESCRIPTION	PACKAGING
301604	Slip Tip	100/box, 400/case
309604	BD Luer-Lok™ Tip 1/5 mL	100/box, 400/case
305482	Eccentric Tip	100/box, 400/case
309695	BD Luer-Lok™ Tip, Control Syringe	25/box, 100/case
305462	Lab Syringe, Eccentric Tip, Non-Sterile	100/box, 400/case

#### 20 mL Syringe Only

PROD #	DESCRIPTION	PACKAGING
301625	20 mL Syringe, Slip Tip	40/shelf pack, 160/case
309661	20 mL Syringe, BD Luer-Lok™ Tip	40/shelf pack, 160/case
300613	20 mL Syringe, Eccentric Tip	120/shelf pack, 480/case

#### 30 mL Syringe Only

PROD #	DESCRIPTION	PACKAGING
309651	30 mL Syringe, Slip Tip	40/shelf pack, 160/case
309650	30 mL Syringe, BD Luer-Lok <sup>™</sup> Tip	40/shelf pack, 160/case

#### 60 mL Syringe Only

PROD #	DESCRIPTION	PACKAGING
309654	60 mL Syringe, Slip Tip	40/shelf pack, 160/case
309653	60 mL Syringe, BD Luer-Lok™ Tip	40/shelf pack, 160/case
300866	60 mL Syringe, Eccentric Tip	60/shelf pack, 240/case

#### 2 oz Catheter Tip Syringe Only (with Tip Shield)

PROD #	DESCRIPTION	PACKAGING
309620	Catheter Tip	40/box, 160/case

#### **BD PrecisionGlide<sup>™</sup> Needles Regular Bevel**

PROD #	DESCRIPTION	PACKAGING
305106	30 GA x 1/2 in.	100/shelf pack, 1000/case
305128	30 GA x 1 in.	100/shelf pack, 1000/case
301629	27 GA x 1 1/2 in.	100/shelf pack, 1000/case
305109	27 GA x 1/2 in	100/shelf pack, 1000/case
305136	27 GA x 1 1/4 in.	100/shelf pack, 1000/case
305110	26 GA x 3/8 in.	100/shelf pack, 1000/case
305111	26 GA x 1/2 in.	100/shelf pack, 1000/case
305115	26 GA x 5/8 in.	100/shelf pack, 1000/case
305124	25 GA x 7/8 in.	100/shelf pack, 1000/case
305122	25 GA x 5/8 in.	100/shelf pack, 1000/case
305125	25 GA x 1 in.	100/shelf pack, 1000/case
305127	25 GA x 1 1/2 in.	100/shelf pack, 1000/case
305145	23 GA x 1 in.	100/shelf pack, 1000/case
305143	23 GA x 3/4 in.	100/shelf pack, 1000/case
305120	23 GA x 1 1/4 in.	100/shelf pack, 1000/case
		continued on top of next column>

305193	23 GA x 1 in. Thin Wall	100/shelf pack, 1000/case
305194	23 GA x 1 1/2 in. Thin Wall	100/shelf pack, 1000/case
305155	22 GA x 1 in.	100/shelf pack, 1000/case
305156	22 G x 1 1/2 in.	100/shelf pack, 1000/case
305165	21 GA x 1 in.	100/shelf pack, 1000/case
305167	21 GA x 1 1/2 in.	100/shelf pack, 1000/case
305129	21 GA x 2 in.	100/shelf pack, 1000/case
305177	21 GA x 1 in. Thin Wall	100/shelf pack, 1000/case
305175	20 GA x 1 in.	100/shelf pack, 1000/case
305176	20 GA x 1 1/2 in.	100/shelf pack, 1000/case
305186	19 GA x 1 in. Thin Wall	100/shelf pack, 1000/case
305187	19 GA x 1 1/2 in. Thin Wall	100/shelf pack, 1000/case
305195	18 GA x 1 in.	100/shelf pack, 1000/case
305196	18 GA x 1 1/2 in.	100/shelf pack, 1000/case
305185	18 GA x 1 1/2 in. Thin Wall	100/shelf pack, 1000/case
305197	16 GA x 1 in.	100/shelf pack, 1000/case
305198	16 GA x 1 1/2 in.	100/shelf pack, 1000/case

#### **BD PrecisionGlide<sup>™</sup> Needles Short Bevel**

PROD #	DESCRIPTION	PACKAGING
305159	22 GA x 1 1/2 in.	100/shelf pack, 1000/case
305178	20 GA x 1in.	100/shelf pack, 1000/case
305179	20 GA x 1 1/2 in.	100/shelf pack, 1000/case
305188	19 GA x 1 in. Thin Wall	100/shelf pack, 1000/case
305189	19 GA x 1 1/2 in. Thin Wall	100/shelf pack, 1000/case
305199	18 GA x 1 1/2 in.	100/shelf pack, 1000/case

#### 1 mL Tuberculin Syringe with Needle

PROD #	DESCRIPTION	PACKAGING
309623	TB Syringe with 27 GA x 1/2 in. Needle	100/shelf pack, 800/case
309625	TB Syringe with 26 GA x 3/8 in. Needle	100/shelf pack, 800/case
309626	TB Syringe with 25 GA x 5/8 in. Needle	100/shelf pack, 800/case
309624	TB Syringe with 21 GA x 1 in. Needle	100/shelf pack, 800/case

#### 3 mL Syringe with Detachable Needle

PROD #	DESCRIPTION	PACKAGING
309582	3 mL Syringe, 25 GA x 1 1/2 in. Needle	100/shelf pack, 800/case
309581	3 mL Syringe, 25 GA x 1 in. Needle	100/shelf pack, 800/case
309570	3 mL Syringe, 25 GA x 5/8 in. Needle	100/shelf pack, 800/case
309571	3 mL Syringe, 23 GA x 1 in. Needle	100/shelf pack, 800/case
309574	3 mL Syringe, 22 GA x 1 1/2 in. Needle	100/shelf pack, 800/case
309573	3 mL Syringe, 22 GA x 1 1/4 in. Needle	100/shelf pack, 800/case
309569	3 mL Syringe, 22 GA x 3/4 in. Needle	100/shelf pack, 800/case
309572	3 mL Syringe, 22 GA x 1 in. Needle	100/shelf pack, 800/case
309575	3 mL Syringe, 21 GA x 1 in. Needle	100/shelf pack, 800/case
309578	3 mL Syringe, 20 GA x 1 in. Needle	100/shelf pack, 800/case

### **BD<sup>™</sup> Blunt Fill Needles**

#### **Blunt Fill Needles**

PROD #	DESCRIPTION	PACKAGING
305180	BD™ Blunt Fill Needle 18 GA x 1 1/2 in.	100/box, 1000/case
305211	BD <sup>™</sup> Blunt Filter Needle 18 GA x 1 1/2 in.	100/box, 1000/case

#### 3 mL Syringe with BD Blunt Fill Needle

PROD #	DESCRIPTION	PACKAGING	
305060	BD Luer-Lok™ Tip 3 mL with 18 GA x 1 1/2 in. Blunt Fill Needle	100/box, 800/case	

#### 5 mL Syringe with BD Blunt Fill Needle

PROD #	DESCRIPTION	PACKAGING
305062	BD Luer-Lok™ Tip 5 mL with 18 GA x 1 1/2 in.	100/box, 400/case
	Blunt Fill Needle	

#### 10 mL Syringe with BD Blunt Fill Needle

PROD #	DESCRIPTION	PACKAGING
305064	BD Luer-Lok <sup>™</sup> Tip 10 mL with 18 GA x 1 1/2 in. Blunt Fill Needle	100/box, 400/case

#### **BD<sup>™</sup> Filter Needles**

BD Filter Needles feature a molded-in 5 Micron filter membrane to reduce the risk of aspirating glass particles when preparing medications from ampoules.

PROD #	DESCRIPTION	PACKAGING
305200	19 GA x 1 1/2 in., Thin Wall 5 Micron filter	100/box, 1000/case
305201	18 GA x 1 1/2 in., Thin Wall 5 Micron filter with BD Nokor™ point	100/box, 1000/case

#### **BD Nokor<sup>™</sup> Vented Needles**

The BD Nokor needle design features a cannula side port to place fluid on the side wall of a container during reconstitution. Some drug manufacturers recommend this to help reduce foaming or spray of powdered medication.

- BD Nokor design eliminates clogged cannula, retaining venting capability
- Vented needles available, designed to reduce internal container pressure

PROD #	DESCRIPTION	PACKAGING
305213	16 GA x 1 in., Thin Wall Non-Coring Vented Needle	100/box, 1000/case
305214	18 GA x 1 in., Thin Wall Non-Coring Vented Needle	100/box, 1000/case

#### **BD Nokor<sup>™</sup> Admix Needles**

**Efficiency for batch filling.** With features like a stopper at the bottom-

most position to reduce manipulation,

BD Sterile Syringe Convenience Trays help pharmacy staff save time and

effort when batch-preparing syringes.

PROD #	DESCRIPTION	PACKAGING
305216	16 GA x 1 in. Admix non-coring Needle with thin wall	100/box, 1000/case

#### BD<sup>™</sup> Sterile Syringe Convenience Trays



PROD #	DESCRIPTION	PACKAGING
309701	1 mL Syringe, BD Luer-Slip <sup>™</sup> Tip	25/tray, 12 trays/case of 300
309702	3 mL Syringe, BD Luer-Lok™ Tip	25/tray, 12 trays/case of 300
309703	5 mL Syringe, BD Luer-Lok™ Tip	25/tray, 12 trays/case of 300
309605	10 mL Syringe, BD Luer-Lok™ Tip	20/tray, 12 trays/case of 240
305617	20 mL Syringe, BD Luer-Lok™ Tip	10/tray, 12 trays/case of 120
305618	30 mL Syringe, BD Luer-Lok™ Tip	10/tray, 12 trays/case of 120
309680	60 mL Syringe, BD Luer-Lok™ Tip	20/tray, 6 trays/case of 120

#### **BD Tip Caps**

PROD #	DESCRIPTION	PACKAGING
308341	Luer Tip Cap Trays	10/tray, 200 trays, 2000/case
305819	Syringe Tip Cap	200/box, 2000/case

#### **BD Oral Syringe System**

Recommended for oral medications for patient safety. BD offers a broad range of syringes and accessories for oral medication administration. Both clear and amber syringes feature mL and teaspoon scale markings for accuracy. The amber syringes meet USP requirements for photo-sensitive medications.



PROD #	DESCRIPTION	PACKAGING
305217	1 mL Syringe, Clear with Tip Cap	100/box, 500/case
305220	3 mL Syringe, Clear with Tip Cap	100/box, 500/case
305218	5 mL Syringe, Clear with Tip Cap	100/box, 500/case
305219	10 mL Syringe, Clear with Tip Cap	100/box, 500/case
305207	1 mL Syringe, Amber with Tip Cap	100/box, 500/case
305210	3 mL Syringe, Amber with Tip Cap	100/box, 500/case
305208	5 mL Syringe, Amber with Tip Cap	100/box, 500/case
305209	10 mL Syringe, Amber with Tip Cap	100/box, 500/case

#### **BD** Tamper-Tuf<sup>™</sup> Medication Containers

BD Tamper-Tuf containers are offered in two sizes (15 mL or 30 mL) providing extra convenience with a large filling port, snap-top and volume markings. The containers are available unassembled for non-liquid and liquid oral medication or assembled for liquid medication.

PROD #	DESCRIPTION	PACKAGING
305203	15 mL Container, Assembled	250/case
305204	30 mL Container, Assembled	250/case
305205	15 mL Container, Unassembled	250/case
305206	30 mL Container, Unassembled	250/case

#### BD Cornwall<sup>™</sup> Disposable Syringe System

Comfortable, convenient fluid dispensing. The BD Cornwall Disposable Syringe System offers a user-friendly design for fluid transfer – without the capital investment required for filling pumps. Optional syringe-filling connectors allow fluid transfer to a wide variety of reservoirs, including vials, syringes and oral medication containers.



PROD #	DESCRIPTION	PACKAGING
305224	10 mL Fluid BD Cornwall Disposable Syringe	10/case
305225	Syringe Tip Connector	50/case
305223	Oral Syringe Tip Connector	50/case

#### **BD<sup>™</sup> Insulin Syringes**

## 3/10 mL BD Lo-Dose<sup>™</sup> Insulin Syringe with Permanently Attached Needle

PROD #	DESCRIPTION	PACKAGING
309300	28 GA x 1/2 in U-100 BD Micro-Fine™ IV	Self-Contained, 10ea/bag, 200/box, 1,600/case
309301	29 GA x 1/2 in U-100 BD Ultra-Fine™	Self-Contained 10ea/bag, 200/box, 1,600/case

#### 1/2 mL BD Lo-Dose Insulin Syringe with Permanently Attached Needle

PROD #	DESCRIPTION	PACKAGING
329461	28 GA x 1/2 in U-100 BD Micro-Fine IV (Orange)	Blister, 100/box, 500/case
309305	28 GA x 1/2 in U-100 BD Micro-Fine IV (Orange)	200/box, 1,600/case
309306	29 GA x 1/2 in U-100 BD Ultra-Fine	200/box, 1,600/case

#### 1 mL BD Insulin Syringe with Permanently Attached Needle

PROD #	DESCRIPTION	PACKAGING
309309	28 GA x 1/2 in U-100 BD Micro-Fine IV (Orange)	Self-Contained, 200/box, 1,600/case
309311	29 GA x 1/2 in U-100 BD Ultra-Fine	Self-Contained, 200/box, 1,600/case
329424	28 GA x 1/2 in U-100 BD Micro-Fine IV (Orange)	Blister, 100/box, 500/case
329420	28 GA x 1/2 in U-100 BD Micro-Fine IV Single unit Gradation	Blister, 100/box, 500/case

#### 1 mL BD Insulin Syringe with Detachable Needle

PROD #	DESCRIPTION	PACKAGING
329651	25 GA x 5/8 in U-100 (Orange/Blue)	Blister, 100/box, 800/case
329652	26 GA x 1/2 in U-100 (Orange/Tan)	Blister, 100/box, 800/case
329622	25 GA x 1 in U-100 (Orange/Blue)	Blister, 100/box, 800/case

#### 1 mL BD Insulin Syringe Only

PROD #	DESCRIPTION	PACKAGING
329650	U-100 Slip Tip (Orange)	Blister, 100/box, 800/case
309629	U-100 Luer Lock	Blister, 100/box, 800/case



2 2

いいたちをあるゆ

10876554320

and and and and and and and and

## **BD<sup>™</sup> Reusable Products**

#### **BD Regular Bevel Metal Hub - Nickel Chrome Plated**

PROD #	DESCRIPTION	PACKAGING
511252	30 GA 1/2 in.	12/box
511258	30 GA 1-1/2 in.	12/box
511005	27 GA 1 in.	12/box
511012	26 GA 1/2 in.	12/box
511023	25 GA 5/8 in.	12/box
511028	25 GA 2 in.	12/box
511047	23 GA 1-1/2 in.	12/box
511055	22 GA 1 in.	12/box
511057	22 GA 1-1/2 in.	12/box
511059	22 GA 3 in.	12/box
511067	21 GA 1-1/2 in.	12/box
511077	20 GA 1-1/2 in.	12/box
511078	20 GA 2 in.	12/box
511089	19 GA 3 in.	12/box
511097	18 GA 1-1/2 in.	12/box
511098	18 GA 2 in.	12/box
511108	15 GA 3-1/2 in.	12/box
511110	13 GA 3-1/2 in.	12/box

#### **BD Short Bevel Metal Hub - Nickel Chrome Plated**

PROD #	DESCRIPTION	PACKAGING
511104	16 GA 2 in.	12/box
511109	13 GA 2 in.	12/box

#### BD Tuberculin Syringe, Precision with Glass Luer Tip

PROD #	DESCRIPTION	PACKAGING
512001	1/4 Yale Tuberculin 1/4 mL in. 1/100 mL and 1/2 min.	1 ea
512002	1/2 Yale Tuberculin 1/2 mL in. 1/100 mL and 1/2 min.	1 ea
512004	1 Yale Tuberculin 1 mL in. 1/100 mL and 1/2 min.	6 ea
512027	1YTL 1 mL BD Luer-Lok™ Tip	1 ea
512019	2YB 2 mL in. 1/20 mL	1 ea

#### **BD Multifit<sup>™</sup> Reuseable Syringes**

PROD #	DESCRIPTION	PACKAGING
512140	2 mL BD Multifit Reusable Syringe	1 ea
512146	3 mL BD Multifit Reusable Syringe	1 ea
512141	5 mL BD Multifit Reusable Syringe	1 ea
512142	10 mL BD Multifit Reusable Syringe	1 ea

#### BD Cornwall<sup>™</sup> Control Syringe with BD Luer Lok Tip<sup>™</sup>

PROD #	DESCRIPTION	PACKAGING
512190	1240S 1 mL	1 ea
512191	1250S 2 mL	1 ea
512192	1260S 5 mL	1 ea
512193	1270S 10 mL	1 ea
512907	10 mL BD Multifit™ glass control syringe with BD Luer-Lok™ Tip	1 ea

Legend: Syringe Tip Abbreviations		
MULTIFIT	YALE	
mLLLuer-Lock Tip	YLLuer-Lock Tip	
MMMetal Luer Tip	YGlass Luer Tip	
MGlass Luer Tip	YMMetal Luer Tip	
MLCLuer-Lock, Control	YEGlass Luer Eccentric Tip	
	YLELuer-Lock Eccentric Tip	

#### **BD Reusable Syringe Individual Packages**

PROD #	DESCRIPTION	PACKAGING
512300	2 mL BD Yale Glass Syringe with Glass Luer Tip	1 ea
512310	2 mL BD Yale Glass Syringe with Luer-Lock™ Tip	1 ea
512306	3 mL BD Yale Glass Syringe with Glass Luer Tip	1 ea
512316	3 mL BD Yale Glass Syringe with Luer-Lock™ Tip	1 ea
512301	5 mL BD Yale Glass Syringe with Glass Luer Tip	1 ea
512311	5 mL BD Yale Glass Syringe with Luer-Lock™ Tip	1 ea
512302	10 mL BD Yale Glass Syringe with Glass Luer Tip	1 ea
512312	10 mL BD Yale Glass Syringe with Luer-Lock™ Tip	1 ea
512303	20 mL BD Yale Glass Syringe with Glass Luer Tip	1 ea
512313	20 mL BD Yale Glass Syringe with Luer-Lock™ Tip	1 ea
512314	30 mL BD Yale Glass Syringe with Luer-Lock™ Tip	1 ea
512315	50 mL BD Yale Glass Syringe with Luer-Lock™ Tip	1 ea
512317	100 mL BD Yale Glass Syringe with Luer-Lock™ tip	1 ea

#### BD Multifit<sup>™</sup> Syringe Bulk Package

PROD #	DESCRIPTION	PACKAGING
512440	2 mL glass syringe	36/box
512441	3 mL glass syringe	36/box
512442	10 mL glass syringe	36/box
512443	20 mL glass syringe	36/box
512444	30 mL glass syringe	36/box
512445	50 mL glass syringe	36/box

#### **Adaptors**

PROD #	DESCRIPTION	PACKAGING
513095	470V-1 modified two-way automatic valve-AT shaped with female BD Luer-Lok <sup>™</sup> and female BD Luer-Lok <sup>™</sup> side connection (non-return type)	1 ea
513113	mLL/ mLL adapter male Luer to male Luer	1 ea
513328	470F Fistula Tip	1 ea
513397	470W-Black Buna-N rubber washer for (470V and 3410V)	12/box
513576	470RVS Set-Four Buna-N rubber valve tacks and 2 stainless steel springs	1 ea
519067	Male BD Luer-Lok™ adapter	1 ea
513560	1220FO filling outfit, includes 470v valve, 20 in. 470t tubing, 330S metal sinker	1 ea

#### **Reuseable Anesthesia Needles**

PROD #	DESCRIPTION	PACKAGING
511118	LLX extension Needle	1 ea
511140	20 GA x 4 in. BD Yale <sup>™</sup> Quincke spinal Needle	1 ea
511168	23 GA x 1/2 in. curved tonsil Needle	1 ea
511169	23 GA x 3/4 in. straight tonsil Needle	1 ea
511219	17 GA x 3 in. Tuohy epidural Needle, thin wall	1 ea
511285	17 GA x 7 in. Tuohy epidural Needle, thin wall	1 ea
511317	20 GA x 5 in. Pudendal and Parasympathetic block Needle	1 ea
511318	19 GA x 5 in. Pudendal and Parasympathetic block Needle	1 ea
511319	20 GA x 6 in. Pudendal Parasympathetic nerve block Needle	1 ea
511320	18 GA x 6 in. Pudendal and Parasympathetic block Needle	1 ea
511364	20 GA x 6 in. BD <sup>™</sup> local anesthesia security bead Needle	1 ea

#### **Reusable Biopsy Needles**

PROD #	DESCRIPTION	PACKAGING
511348	18 GA x 1 in. bone marrow aspiration and transfusion Needle	1 ea
511346	16 GA x 1 5/16 in. bone marrow aspiration and transfusion Needl	e 1 ea
511403	16 GA x 4 in. Liver, kidney and soft tissue	1 ea
511547	14 GA x 4 1/2 in. Liver, kidney and soft tissue	1 ea
511548	14 GA x 3 3/8 in. Liver, kidney and soft tissue	1 ea



## Infusion Therapy Devices

# **BD Nexiva<sup>™</sup> Closed IV Catheter System**

Taking subtle cues from nature, we've developed the profoundly simple, yet radically advanced BD Nexiva<sup>™</sup> Closed IV Catheter System. The first of its kind, this all-in-one system is designed to reduce insertion attempts increasing first-stick success, as well as reduce your exposure to blood through the innovative blood-containment system. Having an all-in-one system decreases your exposure to blood and the potential risk of contamination.

#### - BD Q-Syte<sup>™</sup> Luer Access Split Septum

A testament to the power of simplicity. Its minimal design offers maximum benefits, including high flow rates, exceptional reliability, clear visibility and cleanability.

#### Dual Access

A second port allows for simultaneous administration of compatible fluids and medications near the insertion site. A second BD Q-Syte device is included in the package.

#### All-in-One Closed System

The pre-assembled system creates a closed single lumen fluid path, designed to minimize blood leakage from the catheter hub, and the potential for contamination and exposure to blood.

#### Built-in Stabilization Device

A soft, flexible stabilization platform is integrated into the catheter system. It is designed for stability and comfort to minimize catheter movement in the vessel, as well as irritation to the vessel.

#### ⊕ BD Instaflash<sup>™</sup> Needle Technology

Designed to reduce hit-and-miss insertion by confirming the vessel entry, enabling you to stay focused on the insertion site.

## **BD Nexiva<sup>™</sup> Closed IV Catheter System**

PROD #	COLOUR CODE	GAUGE SIZE	CATHETER LENGTH	CATHETER O.D. (mm)	CATHETER I.D. (mm)	EXTENSION TUBE I.D. (mm)	UNITS BOX/CASE	FLOW RATE H <sub>2</sub> 0 (mL/hr)
383530	Yellow	24	0.56	0.71	0.36	1.22	20/80	840
383531	Yellow	24	0.75	0.71	0.36	1.22	20/80	780
383532	Blue	22	1.00	0.90	0.52	1.22	20/80	1620
383536	Pink	20	1.00	1.10	0.64	1.65	20/80	3240
383537	Pink	20	1.25	1.10	0.64	1.65	20/80	3120
383538	Pink	20	1.75	1.10	0.64	1.65	20/80	2880
383539	Green	18	1.25	1.31	0.81	1.65	20/80	5040
383540	Green	18	1.75	1.31	0.81	1.65	20/80	4680

## BD Nexiva<sup>™</sup> Closed IV Catheter System—Single Port

#### **BD Instaflash<sup>™</sup> Needle Technology** Designed to reduce hit-and-miss insertion

by confirming vessel entry, enabling you

to stay focused on the insertion site.

#### Built-in Stabilization Platform

The built-in stabilization platform is soft and flexible, and is designed to help minimize catheter movement in the vessel, which may improve catheter dwell time.

#### All-in-One Closed System

The pre-assembled system creates a closed fluid path, designed to minimize blood leakage from the catheter hub, and minimize the potential for contamination and exposure to blood. The system is also capable of withstanding power injections for contrast-enhanced CT scans.<sup>1</sup>

#### BD Vialon<sup>™</sup> Biomaterial

Clinically-proven BD Vialon biomaterial has longer indwell times and softens up to 70% in the vessel, reducing the chance of mechanical phlebitis by up to 50%.<sup>1,2</sup>

<sup>1</sup> Maki D, Ringer, M. Risk factors for infusion-related phlebitis with small peripheral venous catheters. *Annals of Internal Medicine*. 1991;114:845-854. <sup>2</sup> Data on file.

## **BD Nexiva<sup>™</sup> Closed IV Catheter System—Single Port**

PROD #	COLOUR CODE	GAUGE SIZE	CATHETER LENGTH (in)	CATHETER O.D. (mm)*	CATHETER I.D. (mm)*	EXTENSION TUBE I.D. (mm)	UNITS BOX/CASE	FLOW RATE H <sub>2</sub> 0 (mL/hr)
383510	Yellow	24	0.56	0.71	0.53	1.22	20/80	960
383511	Yellow	24	0.75	0.71	0.53	1.22	20/80	900
383512	Blue	22	1.00	0.90	0.67	1.22	20/80	1800
383516	Pink	20	1.00	1.10	0.83	1.65	20/80	3240
383517	Pink	20	1.25	1.10	0.83	1.65	20/80	3120
383518	Pink	20	1.75	1.10	0.83	1.65	20/80	2880
383519	Green	18	1.25	1.31	0.98	1.65	20/80	4860
383520	Green	18	1.75	1.31	0.98	1.65	20/80	4560

For power injection use, refer to Instructions For Use. \*Average dimension.

# BD Insyte<sup>™</sup> Autoguard<sup>™</sup> Shielded IV Catheters

## If I can get the first stick, everything else will follow.



BD Insyte Autoguard shielded IV catheters mean added safety for you. The BD Insyte Autoguard push-button needle shielding technology was proven to effectively safeguard healthcare workers, with a 95% reduction in needlestick injuries.<sup>1</sup>

#### **Protection Is In Your Control**

Push-button needle protection gives clinicians more flexibility with the catheter, providing control of when to activate the safety mechanism. According to a clinical study, this push-button needle shielding technology has significantly reduced needlestick injuries to healthcare workers.<sup>1</sup>



BD Instaflash<sup>™</sup> Needle Technology Designed to reduce hit-and-miss insertion.



1 "Mendelcon MH, Chen LBY, Finkselstein LE, Bailey E, Kogan G, Mt Sinai medical center, New York, NY. Preliminary Evaluation of a Safety IV catheter (Insyte Autoguard, Becton Dickinson) using the Centers for Disease Control and Prevention (CDC) National Surveillance System for Hospital Healthcare workers Database. Abstr. Centers for Disease Control and Prevention - 4th Decennial International Conference on Nosocomial and Healthcare-Associated Infections, March 2000."

#### BD Insyte<sup>™</sup> Autoguard<sup>™</sup> Shielded IV Catheter

BD Insyte<sup>™</sup> Autoguard<sup>™</sup> Shielded IV Catheters mean added safety for you. A study showed that one hospital reduced its needlestick injuries by 95% by switching to BD Autoguard<sup>™</sup> Shielding Technology. The unique BD Instaflash<sup>™</sup> Needle Technology in the 20, 22, and 24 gauge sizes provides immediate confirmation of vessel entry



along the catheter shaft. Made of clinically proven BD Vialon™ Biomaterial.

PROD #	DESCRIPTION	PACKAGING
381412	24 GAx 0.75 in. (0.7 mm x 19 mm) Yellow	50/box, 200/case
381423	22 GA x 1.00 in. (0.9 mm x 25 mm) Blue	50/box, 200/case
381433	20 GA x 1.00 in. (1.1 mm x 25 mm) Pink	50/box, 200/case
381434	20 GA x 1.16 in. (1.1 mm x 30 mm) Pink	50/box, 200/case
381437	20 GA x 1.88 in. (1.1 mm x 48 mm) Pink	50/box, 200/case
381444	18 GA x 1.16 in. (1.3 mm x 30 mm) Green	50/box, 200/case
381447	18 GA x 1.88 in. (1.3 mm x 48 mm) Green	50/box, 200/case
381454	16 GA x 1.16 in. (1.7 mm x 30 mm) Gray	50/box, 200/case
381457	16 GA x 1.77 in. (1.7 mm x 45 mm) Gray	50/box, 200/case
381467	14 GA x 1.75 in. (2.1 mm x 45 mm) Orange	50/box, 200/case

#### With Wings

PROD #	DESCRIPTION	PACKAGING
381512	24 GA x 0.75 in. (0.7 mm x 19 mm) Yellow	50/box, 200/case
381523	22 GA x 1.00 in. (0.9 mm x 25 mm) Blue	50/box, 200/case
381533	20 GA x 1.00 in. (1.1 mm x 25 mm) Pink	50/box, 200/case
381534	20 GA x 1.16 in. (1.1 mm x 30 mm) Pink	50/box, 200/case
381537	20 GA x 1.88 in. (1.1 mm x 48 mm) Pink	50/box, 200/case
381544	18 GA x 1.16 in. (1.3 mm x 30 mm) Green	50/box, 200/case
381547	18 GA x 1.88 in. (1.3 mm x 48 mm) Green	50/box, 200/case
381554	16 GA x 1.16 in. (1.7 mm x 30 mm) Gray	50/box, 200/case
381557	16 GA x 1.77 in. (1.7 mm x 45 mm) Gray	50/box, 200/case

#### BD Insyte-N<sup>™</sup> Autoguard<sup>™</sup> Shielded IV Catheter

BD Insyte-N<sup>™</sup> Autoguard<sup>™</sup> Shielded IV Catheters were designed specifically for neonates and incorporates the patented BD Autoguard<sup>™</sup> Shielding Technology. The unique BD Instaflash<sup>™</sup> Needle Technology provides immediate confirmation of vessel entry along the catheter shaft. Made of clinically proven BD Vialon<sup>™</sup> Biomaterial.



PROD #	DESCRIPTION	PACKAGING
381411	24 GA x 0.56 in. (0.7 mm x 14 mm) Yellow	50/box, 200/case
381511	24 GA x 0.56 in. (0.7 mm x 14 mm) with Wings, Yellow	50/box, 200/case

#### BD Angiocath<sup>™</sup> Autoguard<sup>™</sup> Shielded IV Catheter

BD Angiocath<sup>™</sup> Autoguard<sup>™</sup> Shielded IV Catheters combine the proven shielding technology with the name you trust in IV catheters, BD Angiocath.<sup>™</sup> The FEP polymer catheter material features a thin-wall design and a tapered tip that results in smooth insertion and optimal flow rates. The unique BD Instaflash<sup>™</sup>



Needle Technology in the 20, 22, and 24 gauge sizes provides immediate confirmation of vessel entry along the catheter shaft.

PROD #	DESCRIPTION	PACKAGING
381700	24 GA x 0.75 in. (0.7 mm x 19 mm) Yellow	50/box, 200/case
381701	22 GA x 1.00 in. (0.9 mm x 25 mm) Blue	50/box, 200/case
381702	20 GA x 1.00 in. (1.1 mm x 25 mm) Pink	50/box, 200/case
381703	20 GA x 1.16 in. (1.1 mm x 30 mm) Pink	50/box, 200/case
381704	20 GA x 1.88 in. (1.1 mm x 48 mm) Pink	50/box, 200/case
381705	18 GA x 1.16 in. (1.3 mm x 30 mm) Green	50/box, 200/case
381706	18 GA x 1.88 in. (1.3 mm x 48 mm) Green	50/box, 200/case
381707	16 GA x 1.16 in. (1.7 mm x 30 mm) Gray	50/box, 200/case
381708	16 GA x 1.77 in. (1.7 mm x 45 mm) Gray	50/box, 200/case
381709	14 GA x 1.75 in. (2.1 mm x 45 mm) Orange	50/box, 200/case

#### BD Angiocath-N<sup>™</sup> Autoguard<sup>™</sup> Shielded IV Catheter

The BD Angiocath-N<sup>™</sup> Autoguard<sup>™</sup> Shielded IV Catheter was designed specifically for neonates and incorporates the patented BD Autoguard<sup>™</sup> Shielding Technology. The FEP polymer catheter material features a thin-wall design and a tapered tip that results in smooth insertion and optimal flow rates. The unique BD Instaflash<sup>™</sup> Needle Technology provides immediate confirmation of vessel entry along the catheter shaft.



 PROD #
 DESCRIPTION
 PACKAGING

 381720
 24 GA x 0.56 in. (0.7 mm x 14 mm) Yellow
 50/box, 200/case

#### BD Saf-T-Intima™ Closed IV Catheter System

The BD Saf-T-Intima<sup>™</sup> Closed IV Catheter System combines needlestick safety with a closed system design. The result is a virtually bloodless procedure and greater safety for you, all from the first closed IV catheter system available for vascular access. Made of clinically proven BD Vialon<sup>™</sup> Biomaterial.



#### With Removable PRN Adapter

PROD #	DESCRIPTION	PACKAGING
383312	24 GA x 0.75 in. (0.7 mm x 19 mm) Yellow	25/box, 200/case
383322	22 GA x 0.75 in. (0.9 mm x 19 mm) Blue	25/box, 200/case
383335	20 GA x 1.00 in. (1.1 mm x 25 mm) Pink	25/box, 200/case

#### With Y-Port and Removable PRN Adapter

PROD #	DESCRIPTION	PACKAGING
383313	24 GA x 0.75 in. (0.7 mm x 19 mm) Yellow	25/box, 200/case
383323	22 GA x 0.75 in. (0.9 mm x 19 mm) Blue	25/box, 200/case
383336	20 GA x 1.00 in. (1.1 mm x 25 mm) Pink	25/box, 200/case
383346	18 GA x 1.00 in. (1.3 mm x 25 mm) Green	25/box, 200/case

**Infusion Therapy Devices** 

### **Infusion Therapy Devices**

All BD<sup>™</sup> IV Catheters incorporate the sharp BD<sup>™</sup> Needle Technology and a tapered catheter tip to facilitate penetration and insertion. Gauge sizes are colour coded for easy identification. BD Instaflash<sup>™</sup> Needle Technology provides immediate visual confirmation of vessel entry and is designed to enhance first-stick success. Available in 24, 22, and 20 gauge sizes.

#### **BD Vialon<sup>™</sup> Biomaterial**

#### **Clinically Superior**

BD Insyte<sup>™</sup> and BD Insyte-N<sup>™</sup> IV Catheters are made with BD Vialon<sup>™</sup> Biomaterial, clinically proven to<sup>1</sup>:

- Increase the ease of penetration
- Reduce splitting and peelback
- Resist kinking
- Reduce the risk of phlebitis by up to 50%
- Reduce the risk of catheter-related complications
- Increase indwelling times

BD Insyte-N<sup>™</sup> IV Catheters are recommended for access to small or compromised veins. The broad line of BD Angiocath<sup>™</sup> IV Catheters features inert FEP polymer material with high lubricity for easy insertion.

#### **BD Insyte<sup>™</sup> IV Catheters**

PROD #	DESCRIPTION	PACKAGING
381212	24 GA x .75 in. (Yellow)	50/shelf pack, 200/case
381223	22 GA x 1.0 in. (Blue)	50/shelf pack, 200/case
381233	20 GA x 1.0 in. (Pink)	50/shelf pack, 200/case
381234	20 GA x 1.16 in. (Pink)	50/shelf pack, 200/case
381237	20 GA x 1.88 in. (Pink)	50/shelf pack, 200/case
381244	18 GA x 1.16 in. (Green)	50/shelf pack, 200/case
381247	18 GA x 1.88 in. (Green)	50/shelf pack, 200/case
381254	16 GA x 1.16 in. (Gray)	50/shelf pack, 200/case
381257	16 GA x 1.77 in. (Gray)	50/shelf pack, 200/case
381267	14 GA x 1.75 in. (Orange)	50/shelf pack, 200/case

#### **BD Insyte<sup>™</sup> IV Catheters with Wings**

PROD #	DESCRIPTION	PACKAGING
381312	24 GA x .75 in. (Yellow)	50/shelf pack, 200/case
381323	22 GA x 1.0 in. (Blue)	50/shelf pack, 200/case
381333	20 GA x 1.0 in. (Pink)	50/shelf pack, 200/case
381334	20 GA x 1.16 in. (Pink)	50/shelf pack, 200/case
381337	20 GA x 1.88 in. (Pink)	50/shelf pack, 200/case
381344	18 GA x 1.16 in. (Green)	50/shelf pack, 200/case
381347	18 GA x 1.88 in. (Green)	50/shelf pack, 200/case
381354	16 GA x 1.16 in. (Gray)	50/shelf pack, 200/case
381357	16 GA x 1.77 in. (Gray)	50/shelf pack, 200/case

<sup>1</sup> "Mendelcon MH, Chen LBY, Finkselstein LE, Bailey E, Kogan G, Mt Sinai medical center, New York, NY. Preliminary Evaluation of a Safety IV catheter (Insyte Autoguard, Becton Dickinson) using the Centers for Disease Control and Prevention (CDC) National Surveillance System for Hospital Healthcare workers Database. Abstr. Centers for Disease Control and Prevention - 4th Decennial International Conference on Nosocomial and Healthcare-Associated Infections, March 2000."

#### **BD Insyte-N<sup>™</sup> IV Catheters**

PROD #	DESCRIPTION	PACKAGING
381211	24 GA x .56 in. (Yellow)	50/shelf pack, 200/case
381311	24 GA x .56 in with Wings (Yellow)	50/shelf pack, 200/case

#### **BD Angiocath<sup>™</sup> IV Catheters**

PROD #	DESCRIPTION	PACKAGING
381112	24 GA x .75 in. (Yellow)	50/shelf pack, 200/case
381123	22 GA x 1.0 in. (Blue)	50/shelf pack, 200/case
381134	20 GA x 1.16 in. (Pink)	50/shelf pack, 200/case
381137	20 GA x 1.88 in. (Pink)	50/shelf pack, 200/case
381144	18 GA x 1.16 in. (Green)	50/shelf pack, 200/case
381147	18 GA x 1.88 in. (Green)	50/shelf pack, 200/case
381154	16 GA x 1.16 in. (Gray)	50/shelf pack, 200/case
381157	16 GA x 1.88 in. (Gray)	50/shelf pack, 200/case
381164	14 GA x 1.16 in. (Orange)	50/shelf pack, 200/case
381167	14 GA x 1.88 in. (Orange)	50/shelf pack, 200/case

#### **BD Angiocath<sup>™</sup> Catheters for Special Placement**

PROD #	DESCRIPTION	PACKAGING
382258	16 GA x 3.25 in. (Gray)	10/shelf pack, 50/case
382259	16 GA x 5.25 in. (Gray)	10/shelf pack, 50/case
382268	14 GA x 3.25 in. (Orange)	10/shelf pack, 50/case
382269	14 GA x 5.25 in. (Orange)	10/shelf pack, 50/case
382277	12 GA x 3.0 in. (Light Blue)	10/shelf pack, 50/case
382287	10 GA x 3.0 in. (Tan)	10/shelf pack, 50/case

# **BD Posiflush<sup>™</sup> Syringes** Designed To Meet Your Needs

## BD PosiFlush<sup>™</sup> Syringes are designed to:

## Improve patient outcomes

Unique Syringe Design

- Designed to eliminate syringe-induced blood reflux\*
- Designed to prevent solution from entering a non-sterile area of the syringe

# Reduce the risk of medication errors

Clear Labeling

- Greater visibility of syringe contents
- Bold print for clarity
- Colour and bar coded for easy identification and verification



## **Reduce waste and costs**

Stubby Syringe Profile

Selecting smaller size syringes (3, 5 mL) for peripheral lines:

- Reduces storage and disposal costs
- Minimizes environmental waste



This graphic demonstrates the number of syringes that it takes to fill a sharps container. Data shows that you can fit more than twice as many BD PosiFlush 3 mL size syringes into a 5.1 L In-Room Sharps Collector as you can by using a standard 12 mL syringe!\*\*

# Reduce the risk of catheter damage

Standard 10 mL Syringe Diameter

- Generates significantly lower pressure (PSI) compared to standard 3 mL syringes
- All sizes comply with PICC manufacturer flushing recommendations

\*Average reflux as measured in 4Fr PICC; data on file at BD.

\*\*Data on file at BD.



#### BD PosiFlush<sup>™</sup> syringes containing 0.9% sodium chloride.

PROD #	DESCRIPTION	PACKAGING		
306570	Sterile Field (XS) 3 mL	30 Units/Box	8 Boxes/Case	240 Units/Case
306571	Sterile Field (XS) 5 mL	30 Units/Box	8 Boxes/Case	240 Units/Case
306572	Sterile Field (XS) 10 mL	30 Units/Box	8 Boxes/Case	240 Units/Case
306573	Sterile Pathway (SP) 3 mL	30 Units/Box	16 Boxes/Case	480 Units/Case
306574	Sterile Pathway (SP) 5 mL	30 Units/Box	16 Boxes/Case	480 Units/Case
306575	Sterile Pathway (SP) 10 mL	30 Units/Box	16 Boxes/Case	480 Units/Case

All products are Latex-Free, Preservative-Free

# **BD Q-Syte<sup>™</sup> Luer Access Split Septum**









#### Advanced design plus proven material create the ultimate performance in closed luer access

The BD Q-Syte<sup>™</sup> Closed Luer Access Device is a testament to the power of simplicity. Its minimal design offers maximum benefits, including the highest flow rates, exceptional reliability, and evident cleanability.

- Extremely reliable
- Unrestricted flow
- Crystal clear visibility
- Simple to clean, access and connect

#### **BD Q-Syte Closed Luer Access Split Septum**

PROD #	UNITS BOX/CASE	FLOW RATE <sup>1</sup> H₂0 (ML/PER MIN)	FLOW RATE <sup>1</sup> H₂0 (L/PER HR)	PRIMING VOL. (ML)
385100	50/200	525	32	0.10

#### **BD Q-Syte Extension Sets**

PROD #	DESCRIPTION	UNITS BOX/CASE	FLOW RATE <sup>1</sup> H₂0 (ML/MIN)	FLOW RATE <sup>1</sup> H <sub>2</sub> 0 (L/PER HR)	PRIMING VOL. (ML)
385101	15 cm Standard Bore	25/200	446	27	1.14
385102	15 cm Micro Bore	50/200	49	3	0.21

#### **BD Q-Syte Accessory**

PROD #	DESCRIPTION	PACKAGING
385108	Vial Access Adapter 25/Shelf Pack	100/Case

#### Acacia Extension Sets with BD Q-Syte Luer Access Split Septum

PROD #	DESCRIPTION	PACKAGING	FLOW RATE <sup>1</sup> H₂0 (ML/MIN)	PRIMING VOL. (ML)
385150	15 cm (6 in.) BD RightBore™ 18 Ext. Set, spin nut	50/Case	19	0.60
385151	15 cm (6 in.) Small Bore Ext. Set, spin nut	50/Case	2	0.25
385152	15 cm (6 in.) Macro Bore Ext. Set, spin nut	50/Case	22	1.00
385153	36 cm (14 in.) Small Bore Ext. Set, spin nut	50/Case	3	0.55
385161	15 cm (6 in.) Macro Bore Bi-Ext. Set, spin nut	50/Case	20	1.60
385162	15 cm (6 in.) Macro Bore Tri-Ext. Set, spin nut	50/Case	20	2.25
385163	15 cm (6 in.) Small Bore Bi-Ext. Set, spin nut	50/Case	2	0.45
385164	15 cm (6 in.) Small Bore Tri-Ext. Set, spin nut	50/Case	2	0.80
385165	20 cm (8 in.) "Y" Macro Bore Ext. Set, spin nut	50/Case	22	1.40
385166	50 cm (20 in.) "2Y" Macro Bore Ext. Set, spin nut	50/Case	18	3.20

1 Test Method; ISO 10555-5 Annex B (max flow @ 525 ± 25mL/min)

# Interlink<sup>®</sup> System and BD<sup>™</sup> Cannulas

The BD<sup>™</sup> Blunt Plastic Cannula is specifically designed for use with the Interlink<sup>®</sup> injection site, distinctively marked with a coloured alert ring. Also compatible with Abbott LifeShield<sup>®</sup> and Braun Safeline<sup>®</sup> sites.



PROD #	DESCRIPTION	PACKAGING
303345	BD <sup>™</sup> Blunt Plastic Cannula	100/box, 1000/case
303346	3 mL Syringe with BD <sup>™</sup> Blunt Plastic Cannula	100/box, 800/case
303347	5 mL Syringe with BD <sup>™</sup> Blunt Plastic Cannula	100/box, 400/case
303348	10 mL Syringe with BD <sup>™</sup> Blunt Plastic Cannula	100/box, 400/case
303367	Vial Access Cannula	100/box, 1000/case
303401	3 mL Syringe with Vial Access Cannula	100/box, 800/case
303403	5 mL Syringe with Vial Access Cannula	100/box, 400/case
303405	10 mL Syringe with Vial Access Cannula	100/box, 400/case
303369	Threaded Lock Cannula	100/box, 1000/case
303370	Lever Lock Cannula	100/box, 1000/case
303365	Y-Lock Cannula	100/box, 1000/case
303380	BD <sup>™</sup> Blood Collection Assembly	25/box, 200/case

#### **BD Twinpak<sup>™</sup> Dual Cannula Device**

The BD Twinpak<sup>™</sup> Dual Cannula is a truly universal device that simplifies needleless IV access. It offers twice the performance of traditional methods. The BD<sup>™</sup> Twinpak Device features two cannulas in one shield: one red hub, steel cannula for syringe filling and one clear, blunt plastic cannula for safe access into a split septum. The BD Twinpak<sup>™</sup> Device helps to



streamline IV work-flow practices by reducing the number of components required to perform IV push and flush protocols, thereby saving time and reducing costs.

PROD #	DESCRIPTION	PACKAGING
303390	BD Twinpak <sup>™</sup> Cannula	100/box, 1000/case
303391	3 mL Syringe with BD Twinpak™ Cannula	100/box, 800/case
303392	5 mL Syringe with BD Twinpak™ Cannula	100/box, 400/case
303393	10 mL Syringe with BD Twinpak™ Cannula	100/box, 400/case



#### **BD Posiflow<sup>™</sup> Positive Pressure Valve**



PROD #	DESCRIPTION	PACKAGING
385300	Luer Activated Valve	50/box, 200/case
385301	6 in. BD Luer-Lok™ Extension Set	50/box, 200/case

## BD Persist<sup>™</sup> Plus Protective Barrier with Chlorhexidine in Alcohol

BD Persist<sup>™</sup> Plus is a 1% chlorhexidine gluconate preparation combined with 75% ethanol to give an effective antimicrobial barrier. It provides a fast, one-step application to immediately kill all susceptible microbes.

PROD #	DESCRIPTION	PACKAGING
386401	Single Swab Pack	50/box, 500/case
386402	Triple Swab Pack	25/box, 250/case

#### BD<sup>™</sup> J Loop

PROD #	DESCRIPTION	PACKAGING
385121	Connector Loop with Flexible Extension Tubing	50/box, 200/case

#### BD<sup>™</sup> IV Loop

PROD #	DESCRIPTION	PACKAGING
385120	Connector Loop	50/box, 200/case

#### **BD<sup>™</sup> PRN Adapters**

PROD #	DESCRIPTION	PACKAGING
385111	BD Luer-Lok <sup>™</sup> Adapter with Injection Site-Short 3/4 in.	50/box, 200/case
385112	BD Luer-Lok™ Adapter with Injection Site-Mid Size 1 1/16 in.	50/box, 200/case
385113	Male Luer Adapter with Injection Site-Long 1 5/8 in.	50/box, 200/case
385114	BD Luer-Lok <sup>™</sup> Adapter with Injection Site-Long 1 5/8 in.	50/box, 200/case





## **Surgical Products**

# BD Bard-Parker<sup>™</sup> Protected Blade System

#### **Protection Without Compromise**

The BD Bard-Parker<sup>™</sup> Protected Blade System has been designed with the weight and feel that clinicians prefer, with convenient and easy one-handed activation of the protective shield. The blade remains covered during both assembly and disposal to help prevent accidental contact with the blade.

#### Get a lock on safety

The system offers a unique locking shield that provides easy activation:

- Clear shield locks before and after use for protection
- Top activation button enables left and right-handed use
- Audible button click and tactile sensation confirm lock is in place
- Blade assembly and disassembly are quick and easy, no extra steps

#### No compromise to technique

The system provides all the quality you've come to expect from BD Bard-Parker<sup>™</sup>:

- BD Bard-Parker<sup>™</sup> Stainless-Steel Blade
- Stainless-steel handle with conventional weight and feel
- Durable handle design for extended life


#### **BD Bard-Parker<sup>™</sup> Protected Blade System**

#### Protected Blade System Metal Handles

PROD #	FOR USE WITH BLADE SIZE	HANDLE SIZE	PACKAGING
374030	10, 11, 15, 15C	3	5/box
374031	10, 11, 15, 15C	3L	5/box
374040	20-25	4	5/box
374041	20-25	4L	5/box

#### Protected Surgical Blades

PROD #	BLADE SIZE	PACKAGING
373910	10	50/box, 3/case
373911	11	50/box, 3/case
373915	15	50/box, 3/case
373920	20	50/box, 3/case
373921	21	50/box, 3/case
373922	22	50/box, 3/case
373923	23	50/box, 3/case
374305	15C	50/box, 3/case



## Protect yourself. Protect your staff... in the OR.

## BD Bard-Parker<sup>™</sup> Protected Disposable Scalpel

Featuring Locking Retractable Shield

**Protect Yourself. Protect Your Staff.** Designed to minimize the risk of surgical blade injuries, the BD Bard-Parker<sup>™</sup> Protected Disposable Scalpel gives clinicians, throughout the hospital...

## Protection Without Compromise

#### • In Procedure:

Permitting clinicians to focus first on the surgical site, rather than on the instrument.

#### • In Handling:

Providing instrument characteristics of weight, handling, and feel that clinicians prefer.

#### • In Technique:

Allowing clinicians to retract the protective shield, easily and safely with one hand.

#### • In Reliability:

Assuring quality and performance you have come to expect from the BD Bard-Parker<sup>™</sup> brand.

PROD #	BLADE SIZE	PACKAGING	
372610	#10	10/box, 10/case	
372611	#11	10/box, 10/case	
372612	#12	10/box, 10/case	
372613	#12B	10/box, 10/case	
372615	#15	10/box, 10/case	
372620	#20	10/box, 10/case	
372621	#21	10/box, 10/case	
372622	#22	10/box, 10/case	
372623	#23	10/box, 10/case	



BD-Canada

**Surgical Products** 

## **Surgical Products**

#### BD Bard-Parker<sup>™</sup> Safetylock<sup>™</sup> Blades – Sterile BD Rib-Back<sup>™</sup> Carbon Steel Blades, Individually Packaged

PROD #	BLADE SIZE	PACKAGING
371150	10	50/shelf pack, 150/case
371151	11	50/shelf pack, 150/case
371152	12	50/shelf pack, 150/case
371153	15	50/shelf pack, 150/case
371154	20	50/shelf pack, 150/case
371155	21	50/shelf pack, 150/case
371156	22	50/shelf pack, 150/case
371157	23	50/shelf pack, 150/case

#### BD Bard-Parker<sup>™</sup> Sterile Rib-Back<sup>™</sup> Carbon Steel Blades, Individually Packaged

PROD #	BLADE SIZE	PACKAGING
371110	10	50/shelf pack, 150/case
371111	11	50/shelf pack, 150/case
371112	12	50/shelf pack, 150/case
371115	15	50/shelf pack, 150/case
371120	20	50/shelf pack, 150/case
371121	21	50/shelf pack, 150/case
371122	22	50/shelf pack, 150/case
371123	23	50/shelf pack, 150/case

#### BD Bard-Parker Surgical Scalpels Sterile, Disposable Scalpels with Stainless Steel Blade

PROD #	BLADE SIZE	PACKAGING
371610	10	10/box, 100/case
371611	11	10/box, 100/case
371612	12	10/box, 100/case
371615	15	10/box, 100/case
371620	20	10/box, 100/case
371621	21	10/box, 100/case
371623	23	10/box, 100/case
371613	12B	10/box, 100/case

#### BD Bard-Parker<sup>™</sup> Sterile Stainless Steel Blades, Individually Packaged

PROD #	BLADE SIZE	PACKAGING
371210	10	50/shelf pack, 150/case
371700	10A	50/shelf pack, 150/case
371211	11	50/shelf pack, 150/case
371212	12	50/shelf pack, 150/case
371712	12B	50/shelf pack, 150/case
371215	15	50/shelf pack, 150/case
371716	15C	50/shelf pack, 150/case
371220	20	50/shelf pack, 150/case
371221	21	50/shelf pack, 150/case
371222	22	50/shelf pack, 150/case
371223	23	50/shelf pack, 150/case

#### BD Bard-Parker<sup>™</sup> Rib-Back<sup>™</sup> Carbon Steel Blades Nonsterile, Individually Packaged, in 25 Strips of 6

PROD #	BLADE SIZE	PACKAGING
371310	10	150/case
371311	11	150/case
371315	15	150/case
371320	20	150/case
371321	21	150/case
371322	22	150/case
371325	25	150/case
371340	60 Autopsy (fits #8 Handle only)	150/case
371312	12	150/case
371323	23	150/case
371324	24	150/case
371325	31	150/case

#### **BD Bard-Parker<sup>™</sup> Surgical Blade Handles**

PROD #	FOR USE WITH BLADE SIZE	PACKAGING
371030	3/10 through 15	5/case
371031	3L/10 through 15	5/case
371040	4/20 through 23	5/case
371041	4L/20 through 23	5/case
371050	5/10 through 15	5/case
371060	6/20 through 23	5/case
371070	7/10 through 15	5/case
371080	8/Autopsy Handle (for Blade 60)	5/case
371090	9/10 through 15	5/case

## BD<sup>™</sup> E-Z Scrub Impregnated Preoperative Scrub Brushes and Accessories

PROD #	DESCRIPTION	PACKAGING
377479	BD E-Z Scrub with 4% CHG (brush/sponge)	10 dispensers of 30 ea.
374089	BD E-Z Scrub with BD Ultradex <sup>™</sup> Solution; 3% Chloroxylenol (PCMX) with emollients (brush/sponge)	10 dispensers of 30 ea.
371901	Bracket	1 ea.

#### **BD Chest Drainage**

Heimlich Chest Drain Valve-Sterile	
------------------------------------	--

PROD #	DESCRIPTION	PACKAGING
373460	Drainage of the Pleural Cavity	10/case



#### **BD<sup>™</sup> Whitacre Spinal Needles**

#### Standard Lengths

PROD #	DESCRIPTION	PACKAGING
405010	22 GA x 3 1/2 in.	10/box, 50/case
405133	24 GA x 3 1/2 in. High Flow	10/box, 50/case
405138	25 GA x 3 1/2 in. High Flow	10/box, 50/case
405079	27 GA x 3 1/2 in. High Flow	10/box, 50/case

#### Short Length

-		
PROD #	DESCRIPTION	PACKAGING
409442	25 GA x 4 11/16 in.	10/box, 50/case
405140	25 GA x 5 in. High Flow	10/box, 50/case
409443	27 GA x 4 11/16 in. High Flow	10/box, 50/case
405144	27 GA x 5 in. High Flow	10/box, 50/case

#### Spinal Introducers

PROD #	DESCRIPTION	PACKAGING
405160	20 GA x 1 1/4 in. (use with 25 GA & 27 GA Needle or smaller)	25/box, 100/case

#### **BD<sup>™</sup> Quincke Spinal Needles**

PROD #	DESCRIPTION	PACKAGING
405174	18 GA x 3 in.	25/box, 100/case
405184	18 GA x 3 1/2 in.	25/box, 100/case
405071	20 GA x 2 1/2 in.	25/box, 100/case
405172	20 GA x 3 in.	25/box, 100/case
405182	20 GA x 3 1/2 in.	25/box, 100/case
405074	22 GA x 2 1/2 in.	25/box, 100/case
405171	22 GA x 3 in.	25/box, 100/case
405181	22 GA x 3 1/2 in.	25/box, 100/case
400106	23 GA x 3 1/2 in.	25/box, 100/case
405170	25 GA x 3 in.	25/box, 100/case
405180	25 GA x 3 1/2 in.	25/box, 100/case
405164	26 GA x 3 1/2 in.	25/box, 100/case
405081	27 GA x 3 1/2 in.	25/box, 100/case

#### Long Length

PROD #	DESCRIPTION	PACKAGING
408360	18 GA x 6 in.	10/box, 50/case
405211	20 GA x 6 in.	10/box, 50/case
405148	22 GA x 5 in.	10/box, 50/case
405149	22 GA x 7 in.	10/box, 50/case
405234	25 GA x 4 11/16 in.	10/box, 50/case

#### Short Length

PROD #	DESCRIPTION	PACKAGING
405162	20 GA x 1 ½ in.	25/box, 100/case
405161	22 GA x 1 ½ in.	25/box, 100/case
405073	25 GA x 1 in.	10/box, 50/case
405078	25 GA x 2 in.	25/box, 100/case

#### Spinal Introducers

PROD #	DESCRIPTION	PACKAGING
405160	20 GA x 1 1/4 in. (use with 25 GA & 27 GA Needle or smaller)	25/box, 100/case



#### BD<sup>™</sup> Whitacre Pencil Point Available in gauges from 22 to 27 and in lengths from 3 1/2 in. to 5 in.



BD<sup>™</sup> Quincke-Type Point Available in gauges from 18 to 27 and in lengths from 1 in. to 7 in.

#### **BD Perisafe<sup>™</sup> Epidural Needles**

#### Weiss (Fixed Wing) Needles

PROD #	DESCRIPTION	PACKAGING
405191	17 GA x 3 1/2 in. Modified	10/box, 50/case
405193	17 GA x 5 in. Modified	10/box, 50/case

#### Tuohy (Non-Winged) Needle

PROD #	DESCRIPTION	PACKAGING
405092	17 GA x 3 1/2 in. Modified	10/box, 50/case
405018	18 GA x 3 1/2 in. Modified	10/box, 50/case
405028	20 GA x 3 1/2 in. Modified	10/box, 50/case



Tuohy Weiss

#### **BD Perisafe<sup>™</sup> Epidural Mini-Kits**

PROD #	DESCRIPTION	PACKAGING
400273	18 GA X 3 1/2 in. Weiss Epidural Needle 20 GA X 36 in. Polyamide Catheter, Closed End Three Hole with Guide 1-Catheter Connector 122 Micron Filter 10 mL BD Epilor <sup>™</sup> Loss of Resistance Syringe	10/box, 50/case

#### **BD** Perisafe<sup>™</sup> Nylon Epidural Catheters with Connectors

PROD #	DESCRIPTION	PACKAGING
400507	20 GA x 36 in. Closed End for use with 18 GA x 3 1/2 in. Needle	10/box, 50/case

#### **BD Epilor<sup>™</sup> Loss of Resistance Syringe**

PROD #	DESCRIPTION	PACKAGING
405291	7 mL BD Luer-Lok™, Plastic	10/box, 50/case
405292	7 mL Luer-Slip, Plastic	10/box, 50/case





**BD Epilor<sup>™</sup> Plastic Syringe** 

## BD<sup>™</sup> Combined Spinal/Epidural (CSE) System

#### BD Durasafe<sup>™</sup> CSE Needle Set

Provides rapid onset, extended duration, and pain management

- The only CSE matched-set system with BD<sup>™</sup> Whitacre high flow pencil point needle
- BD Durasafe<sup>™</sup> matched sets are designed to allow 12 mm extension of spinal needle when matched hubs are fully engaged

#### BD Durasafe<sup>™</sup> Needle Sets

PROD #	DESCRIPTION	PACKAGING
405139	25 GA x 4 11/16 in. Whitacre / 17 GA x 3 1/2 in. Weiss Needle	10/box, 50/case
405452	27 GA x 4 11/16 in. Whitacre / 18 GA x 3 1/2 in. Weiss Needle	10/box, 50/case



BD Durasafe<sup>™</sup> Combined Spinal/Epidural Set



## BD<sup>™</sup> Anesthesia Products

#### BD<sup>™</sup> Suprapubic Catheter, Sterile, Single Use, Reusable

PROD #	DESCRIPTION	PACKAGING
408289	Bonanno Catheter Tray 1-Teflon® Catheter, 14 GA x 11 in.; 1-Puncture Needle, 18 GA; 1-Adapter Clamp	6/case

#### **BD<sup>™</sup> Local Anesthesia Security Bead Needle, Reusable**

PROD #	DESCRIPTION	PACKAGING
511364	20 GA x 6 in.	1/box

#### **BD<sup>™</sup> Epidural Anesthesia Needle, Reusable**

PROD #	DESCRIPTION	PACKAGING
511219	17 GA 3 in. Tuohy Thin-Wall	1/box
511285	17 GA 7 in. Tuohy Thin-Wall	1/box

#### **BD<sup>™</sup> Pudendal and Parasympathetic Block Needle**

PROD #	DESCRIPTION	PACKAGING
511318	19 GA x 5 in. Sympathetic Nerve Needle	1/box
511317	20 GA x 5 in. Sympathetic Nerve Needle	1/box
511320	18 GA x 6 in. Sympathetic Nerve Needle	1/box
511319	20 GA x 6 in. Parasympathetic Nerve Needle	1/box

#### BD Yale<sup>™</sup> Spinal Needle Reusable, Quincke

PROD #	DESCRIPTION	PACKAGING
511140	20 GA x 4 in.	1/box

#### **BD<sup>™</sup> Nerve Block Needle, Security Bead, Block Bevel**

PROD #	DESCRIPTION	PACKAGING
408348	22 GA x 1 1/2 in. (Regional block with security bead)	10/box, 50/case

#### **BD Pudendal/Local Block Anesthesia Tray**

PROD #	DESCRIPTION	PACKAGING
402601	Pudendal/Local Block Tray with 10 mL Control Syringe	20/case

## Anesthesia Trays

#### **BD<sup>™</sup> Spinal Trays**

Precision Components

- Your choice of spinal needle:
- BD<sup>™</sup> Whitacre High-flow Pencil Point
- Pencil point is designed to spread the dural fibers and help reduce post-dural puncture headaches (PDPH)

#### OR


#### **BD<sup>™</sup> Quincke Spinal Needle**

- Key/slot arrangement of stylet and cannula hubs facilitates proper needle bevel orientation. Precision-formed side hole
- BD<sup>™</sup> Hypodermic Needle
- BD Luer-Lok<sup>™</sup> (latex free) Syringe

#### All Trays Contain:

- BD<sup>™</sup> Filter Needle (5 micron), blunt tip 18 GA x 1.5 in.
- Prep section including towel, swab sticks and gauze pads
- Fenestrated drape (blue paper or clear plastic)
- Sharps stick pad

PROD #	TRAY DESCRIPTION	PACKAGING
405606	Quincke Spinal Needle 22 GA x 3.5 in.	
405607	Quincke Spinal Needle 22 GA x 3.5 in.	
405608	Quincke Spinal Needle 25 GA x 3.5 in.	
405661	Whitacre Spinal Needle 25 GA x 3.5 in.	
405662	Whitacre Spinal Needle 25 GA x 3.5 in.	
405696	Whitacre Spinal Needle 27 GA x 3.5 in.	
405605	Spinal Support Tray with No Procedure Needle	
405691	Whitacre Spinal Needle 25 GA x 3.5 in.	
405660	Whitacre Spinal Needle 25 GA x 3.5 in.	



#### Legend: Syringe Tip Abbreviations

LL ......Luer-Lok Tip LS .....Luer Slip

P......Plastic (clear) hub





## Anesthesia Trays

## **BD**<sup>™</sup> Epidural Trays

#### Precision Components

- BD<sup>™</sup> Modified Tuohy Epidual Needle (non-winged)
- Your choice of glass or plastic BD Epilor<sup>™</sup> Loss of Resistance (LOR) Syringe
- BD<sup>™</sup> Hypodermic Needle (general use)
- BD Luer-Lok<sup>™</sup> Syringe (latex-free)

#### All Trays Contain:

- BD<sup>™</sup> Filter Needle, (5 micron), 18 GA x 1.5 in.
- Prep section including towel, swab sticks and gauze pads
- Fenestrated drape (blue paper or clear plastic)
- Sharps stick pad

PROD #	DESCRIPTION	PACKAGING
406078	Tuohy Epidural Needle 17 GA x 3.5 in	
406031	Tuohy Epidural Needle 17 GA x 3.5 in	
406076	Tuohy Epidural Needle 17 GA x 3.5 in	
406062	Epidural Support Tray with No Procedure Needle	
406063	Epidural Support Tray with No Procedure Needle	



BD SKU #	Tray Style	Epidu	iral N	leed	le		Cat	hete	r an	d Ac	cesso	ories			LOR S	yrin	ge		A	dditiona	al Needles				Additio Syring			A	diti	iona	l Cor	mpoi	nent	s
	SD = Single Deck		Gauge	Length in.	Hub	Depth Markings		Gauge	Length (inches)	Catheter Connector	Insertion Guide	Catheter Support Pad	Syringe & Catheter Tip Protector	Procedure Filter		mL	Type	Hypo 1 in.	Hypo 1.25 in.	Hypo 1.5 in.	Hypo 1.5 in.	19 GA x 1.5 in. Fltr Asp	Sharps Stick Pad	3 mL	5 mL	10 mL	20 mL	Adhesive Strip, 3/4 in. x 3 in.	Towel	Swab Sticks	Gauze (3 in. x 3 in.)	Tegaderm/Uniflex Dressing	Fenestrated Drape	Clear Fenestrated Drape
406078	SD	Tuohy	17	3.5	м	х	n/a								Plastic	7	LL			18 GA	25 GA	х	х	LL		LL			х	3	3		х	
406031	SD	Tuohy	17	3.5	М	Х	Polyamide	19	36	BD	х	х	х	х	Glass	5	LS	25 GA		18 GA	22 GA	х	х	LL			LL		х	3	3		х	
406076	SD	Tuohy	17	3.5	М	Х	n/a								Glass	5	LS	25 GA		18 GA	22 GA	х	х	LL			LL		х	3	3		х	
406062	SD	No Needle					n/a								Glass	10	LL		27 GA	18 GA	22 GA	х	х		LL/LS	LL			х	3	3	т		х
406063	SD	No Needle					n/a								Plastic Glass	7 5	LS LL			18 GA	22 GA/25 GA	х	х	LL			LL	х	х	3	3		х	

Legend: Syringe	Tip Abbreviations
LLLuer-Lok Tip	LSLuer Slip
LORLoss of resistance	MMetal Hub

**Anesthesia Products** 

## Anesthesia Trays

### **BD<sup>™</sup> Anesthesia Safety Trays**

#### Premium Trays

- Spinal and epidural trays with safety engineered devices
- Standard and customized trays available

#### Precision Components

- BD SafetyGlide<sup>™</sup> shielding hypodermic needle allows you to shield the needle immediately after use to reduce the potential of accidental needlesticks
- BD SafetyGlide<sup>™</sup> needle securely attaches to any standard syringe
- Point-Lok<sup>®</sup> device lets you cover needle tips after use; but, prior to disposal, requires only a single-handed, single-step operation




#### Point-Lok<sup>®</sup> sharps safety device for your safety

• Place Point-Lok® device on flat surface, insert tip of used sharps device and press down. The sharp is now locked in the device.

#### BD SafetyGlide<sup>™</sup> shielding hypodermic needle for your safety

- Activation-Assist<sup>™</sup> technology for fast, easy needle tip shielding
- Single-handed activation
- Stainless steel latch securely shields needle tip after activation to reduce the potential of accidental needlesticks



## • Prep basin, prep swabs, gauze pads,

towel, and fenestrated drape

#### All trays with catheters also contain:

 0.22 micron filter, catheter connector, catheter pad, and catheter warning labels

#### **Legend: Syringe Tip Abbreviations**

- CP ...... Cleap plastics drape LOR ..... Loss of resistance
- SD ......Single-deck tray with a removable prep tray



# BD<sup>™</sup> Radiology and Biopsy Products


#### BD<sup>™</sup> Angiography Needles for Seldinger Technique

Potts-Cournand Needle, Thin-Wall

<b>PROD #</b> 408270	DESCRIPTION 18 GA x 2 3/64 in. for .032 in038 in. Guide Wire	PACKAGING 10/box. 50/case
		· · · · · · · · · · · · · · · · · · ·
	leedle Outer Blunt Cannula	DACKACINC
Arterial N PROD # 408273	Jeedle Outer Blunt Cannula DESCRIPTION 18 GA x 2 7/8 in for 032 in - 038 in Guide Wire	PACKAGING 10/box, 50/case

#### PROD # DESCRIPTION PACKAGING

408295	18 GA x 2 3/4 in. One-Part Needle	10/box, 50/case
	with Wings for .032 in038 in. Guide Wire	



One-Part Needle with Wings



## BD<sup>™</sup> Radiology and Biopsy Products

#### **BD<sup>™</sup> Arthrogram Trays**

PROD #	DESCRIPTION	PACKAGING
408501	Arthrogram Tray, Sterile, Single Use	
	Tray consists of:	
	1 - 5 mL BD Luer-Lok <sup>™</sup> Syringe	1 - Absorbent Fenestrated Drape
	1 - 10 mL BD Luer-Lok™ Syringe	2 - Small Basins
	1 - 30 mL BD Luer-Lok™ Syringe	2 - Absorbent Towels
	1 - 18 GA x 1 1/2 in. Transfer Needle	1 - 2 oz. Medicine Cup
	1 - 20 GA x 1 1/2 in. Short Bevel Needle	1 - 10 in. Connecting Tube
	1 - 25 GA x 5/8 in. Needle	1 - Adhesive Bandage
	5 - 3 in. x 3 in. Gauze Pads	2 - Prep Sponges

#### **BD<sup>™</sup> Westcott Fine Needle Aspiration Biopsy**

PROD #	DESCRIPTION	PACKAGING
408265	22 GA x 3 1/2 in.	10/box, 50/case
408266	22 GA x 5 1/2 in.	10/box, 50/case
408267	22 GA x 7 in.	10/box, 50/case
408262	20 GA x 3 1/2 in.	10/box, 50/case
408263	20 GA x 5 1/2 in.	10/box, 50/case
408264	20 GA x 7 in.	10/box, 50/case

#### BD Chiba Fine Needle Aspiration Biopsy/ Transhepatic Cholangiography Needles

PROD #	DESCRIPTION	PACKAGING
408268	22 GA x 6 in. Chiba Type, Thin-Wall	10/box, 50/case
408269	22 GA x 8 in. Chiba Type, Thin-Wall	10/box, 50/case

#### **BD<sup>™</sup> Breast Lesion Localization Needles**

PROD #	DESCRIPTION	PACKAGING
409407	BD PrecisionGlide <sup>™</sup> 20 GA x 3 1/2 in. (8.9 cm) with Hook Wire	10/box, 50/case
409408	BD PrecisionGlide <sup>™</sup> 20 GA x 2 1/2 in. (6.4 cm) with Hook Wire	10/box, 50/case
409409	BD PrecisionGlide <sup>™</sup> 20 GA x 1 1/2 in. (3.8 cm) with Hook Wire	10/box, 50/case

#### **BD<sup>™</sup> Echogenic Procedure Needles** for Use With Ultrasound

PROD #	DESCRIPTION	PACKAGING
400727	BD Ultra-vue <sup>™</sup> 22 GA x 3 1/2 in. Spinal Type Point	10/box, 50/case
400730	BD Ultra-vue <sup>™</sup> 20 GA x 3 1/2 in. Spinal Type Point	10/box, 50/case

#### BD Liver, Kidney and Soft Tissue Biopsy, Reusable

PROD #	DESCRIPTION	PACKAGING
511548	14 GA x 3 3/8 in. Silverman – BD, Modified Franklin, TW	1/box
511547	14 GA x 4 ½ in. Silverman – BD, Modified Franklin, TW	1/box
511403	16 GA x 4 in. Klatskin, TW	1/box

# BD Bone Marrow, Aspiration and Transfusion Needle, Reusable


PROD #	DESCRIPTION	PACKAGING
511346	16 GA x 1 5/16 in. Rosenthal	1/box
511348	18 GA x 1 in. Osgood	1/box





Westcott Fine Needle Aspiration Biopsy

Breast Lesion Localization Needles



Echogenic Procedure Needles for use with Ultrasound

## BD<sup>™</sup> Anesthesia Products

#### **BD<sup>™</sup>** Needle Tubing Adapter, Sterile, Single Use

PROD #	DESCRIPTION	PACKAGING
408205	15 GA x 1/2" .066" I.D071" I.D.	25/box, 100/case
408206	16 GA x 1/2" .058" I.D064" I.D.	25/box, 100/case
408208	18 GA x 1/2" .042" I.D049" I.D.	25/box, 100/case
408210	20 GA x 1/2" .030" I.D035" I.D.	25/box, 100/case
408213	23 GA x 1/2" .020" I.D025" I.D.	25/box, 100/case
408214	24 GA x 1/2" .018" I.D022" I.D.	25/box, 100/case

#### **BD<sup>™</sup>** Catheter Adapter, Sterile, Single Use

PROD #	DESCRIPTION	PACKAGING
403250	Catheter Adapter	100/box, 1000/case

#### **BD<sup>™</sup> Adapters, Sterile**

BD Luer-Lok<sup>™</sup> Adapter, Sterile, Single Use

PROD #	DESCRIPTION	PACKAGING
408299	LL/MLL Adapter Female to Male	50/box, 500/case
BD LUEF-L	ok™ Caps, Sterile, Single Use	PACKAGING
	,	PACKAGING 50/case

#### **BD<sup>™</sup> Miscellaneous Adapters**

PROD #	DESCRIPTION	PACKAGING
408541	BD Luer-Lok™ Female Cap, Non-Sterile	1000/case





Needle Tubing Adapter, Sterile, Single Use

Catheter Adapter, Sterile

#### Adapters, Sterile





408299

408531

408530





## **Sharps Collectors**

## **Sharps Collectors – Chemotherapy**

## Chemotherapy

Chemotherapy Sharps Collectors feature a wide entry port design for disposal of bulky chemotherapy waste. In addition, these collectors are designed to offer increased leak resistance.

## Features:

- Point-first disposal
- Wide entry port for disposal of bulky chemotherapy waste items
- Exclusive safety plug for temporary closure
- Single-piece molding no assembly required
- Sturdy puncture-resistant construction
- Attached locking lid for secure closure and handles for safe transport
- Not autoclaveable

#### Chemotherapy

SIZE (L)	DESCRIPTION	PKG. QTY. 3/cs
6.6	Plug cap	12
10.3	X-large plug cap	12
18.7	X-large plug cap	8
	6.6 10.3	6.6Plug cap10.3X-large plug cap





X-large plug cap



#### Accessories

ACCESSORY #	DESCRIPTION	PKG. QTY.
305497	Quick release bracket	12
305497	Quick release bracket	12
305495	X-large security bracket	12
305495	X-large security bracket	12
	305497 305497 305495	305497Quick release bracket305495X-large security bracket

## Sharps Collectors – Extra Large

## Extra Large

These stackable and sturdy sharps collectors feature large openings for bulky medical waste.

#### Features:

- Large access openings for up to 10 inch diameter medical waste
- Temporary and permanent closures
- Sturdy puncture-resistant construction
- Nesting for efficient storage
- Built in handles for easy handling and disposal

#### Extra Large

PRODUCT #	SIZE (L)	DESCRIPTION	PKG. QTY.
305579	34	10" opening	8
305647	72	10" opening	5





10" opening

#### Accessories

PRODUCT #	ACCESSORY #	DESCRIPTION	PKG. QTY.
305579	305091	Transport trolley, foot operated	1
305647	305093	Transport trolley, foot operated	1

## **Sharps Collectors – Phlebotomy**

## Phlebotomy

Clear top provides visualization of collector fill-levels. Increased fill-level identification helps prevent collector OVER-FILL...a leading cause of disposal injuries.

#### Features:

- Strong and durable puncture resistant construction
- Two-way safety slot securely grips needle hub and permits smooth, one-handed needle removal
- Temporary and permanent closures allow for safer transport and reduced risk of tampering and accidental needlesticks
- Larger all purpose opening permits disposal of bulkier medical wastes, such as swabs and needle holders

#### Phlebotomy

PRODUCT #	SIZE (L)	DESCRIPTION	PKG. QTY.
367216	1.0	Universal Unwinder	60





**Universal Unwinder** 

#### Accessories

PRODUCT #	ACCESSORY #	DESCRIPTION	PKG. QTY.
367216	305970	Tray bracket	12

**Sharps Disposal Products** 

## Sharps Collectors – Point-of-use

## Point-of-use

One-hand, one-step technique... they feature a counter-balanced door. Container provides security against tampering.<sup>1</sup>

#### Features:

- Counterbalanced door returns automatically for one hand, one step disposal
- Counterbalanced door stays shut to deter overfilling
- Safety transport handle
- Sturdy puncture-resistant construction

#### **Point-of-use**

PRODUCT #	SIZE (L)	DESCRIPTION	PKG. QTY.
300475	5.1	Side entry	12
300974	5.1	Horizontal entry	20
303050	7.6	Horizontal entry	10
303051	11.3	Horizontal entry	10



#### Accessories

PRODUCT #	ACCESSORY #	DESCRIPTION	PKG. QTY.
300475	305447	Locking wall bracket	10
300974	305447	Locking wall bracket	10
	305017	Patient room cabinet	1
303050	305409	Bracket	10
303051	305409	Bracket	10

<sup>1</sup> Selecting, Evaluating, and Using Sharps Disposal Containers CDC/NIOSH Publication 1997-111

**Sharps Disposal Products** 

## **Sharps Collectors – One-Piece**

## **One-Piece**

Convenient, vertical "point first" drop, these robust sharps collectors are preassembled and offer a patented one-way funnel valve to minimize needlesticks and overfilling.

#### **Features**

- Point-first disposal
- One-way tamper resistant entry
- Secure cap closure
- Single piece molding no assembly required
- Sturdy puncture-resistant construction
- Attached locking lid for secure closure and safe transport





**One-Piece** 

SIZE (L)

1.4

1.4

3.1

10.3

DESCRIPTION

X-large plug cap

Funnel cap

Funnel cap

Plug cap

PRODUCT #

300460

300462

300466

300452



PKG. QTY.

36

36

24

12

Plug cap

X-large plug cap

# 

**Accessories** 

PRODUCT #	ACCESSORY #	DESCRIPTION	PKG. QTY.
300460	305451	Tray bracket	12
300462	305451	Tray bracket	12
300466	305497	Quick release bracket	12
300452	305497	Quick release bracket	12

**BD-Canada** 

## **Sharps Collectors – Nestable**

## Nestable

With an exclusive funnel-shaped opening and safety neck, these sharps collectors are nestable for easy storage. They allow point-first drop and reduce the risk of tampering and spilling.

#### **Features**

- Clear top helps deter overfilling
- Point-first disposal
- One-way tamper resistant entry and safety neck
- Nestable for efficient storage
- Sturdy carrying handles
- Sturdy puncture-resistant construction
- Secure cap closure for safe transport



**Funnel cap** 



SIZE (L)

3.1

7.6

13.2

22.7

22.7

DESCRIPTION

Point-first vertical entry

X-large funnel cap

Funnel cap

Funnel cap

Open cap

Nestable PRODUCT #

300450

300439

300443

300434

300483

**Open cap** 



PKG. QTY. 3/cs

36

24

20

12

12

X-large funnel cap



#### Accessories

PRODUCT #	ACCESSORY #	DESCRIPTION	PKG. QTY.
300450	305494	Security padlock	12
	305475	Wall cabinet	12
300439	305494	Security padlock	12
	305485	Bracket for nestable containers	12
300443	305494	Security padlock	12
	305485	Bracket for nestable containers	12
300434	305494	Security padlock	12
	305485	Bracket for nestable containers	12

#### Accessories

PRODUCT #	ACCESSORY #	DESCRIPTION	PKG. QTY.
300483	305494	Security padlock	12
	305485	Bracket for nestable containers	12
300988	305494	Security padlock	12
	305485	Bracket for nestable containers	12

**BD-Canada** 

**Sharps Disposal Products** 



LifeShield® is a trademark of Abbot. © 2010 Safeline® is a trademark of Braun. © 2010 BD, BD Logo and all other trademarks are property of Becton, Dickinson and Company. © 2010 BD

2100 Derry Road W.,Suite 100 Mississauga, Ontario L5N 0B3 www.bd.com