

Urology

**tyco**

Healthcare



**KENDALL**

- Specialty Catheters
  - CURITY Suprapubic Trocar Catheters
  - CURITY Temperature Monitoring Silicone Catheters
- Catheters
  - Foley Catheters
    - DOVER 100% Silicone
    - ULTRAMER
    - CURITY Silicone Coated
    - DOVER T.U.R.P. Silicone Elastomer Coated
  - Urethral Catheters
    - CURITY
    - Nelaton
- Urethral Catheterisation Trays
  - DOVER (Open)
  - DOVER (Closed)
- Foley Trays
  - Catheter Insertion Trays
    - DOVER
  - Catheterisation Trays
    - DOVER 100% Silicone
    - DOVER Silicone Coated
    - CURITY ULTRAMER
- Drainage Bags
  - CURITY MONO-FLO
  - DOVER
  - KENGUARD
  - CURITY
  - CURITY Ureteral
  - CURITY T.U.R.P.
- Drainage Bag Foley Trays
  - CURITY MONO-FLO
  - DOVER KENGUARD

- Urine Meters
  - URIPROGRAM 500
  - PRECISION 400
  - DOVER MEASURIT
  - CURITY 200ml
- Urine Meter Foley Trays
  - PRECISION 400
  - DOVER MEASURIT
- Paediatric Urology Products
  - CURITY Urodynamic Catheters
  - PRECISION 150
  - CURITY KIDS
- Urinary Collection
  - Leg Bags
    - MONO-FLO TAG & NACHT
    - CURITY
  - Male External Catheters
    - TEXAS
    - DOVER Catheter Straps
  - URI-DRAIN Extension Tubing
- Bladder Management
  - Vaginal Irrigation (Douche) Bags
  - Enema Bags
  - Flatus Bag
  - Rectal Tube
  - Enema Bucket
- Specimen Collection
  - Commode Specimen Collection
  - SAGE Midstream Catch Kits
  - EASY CATCH Midstream Series
  - CURITY Clean Catch Kit

CURITY, KENGUARD, MONO-FLO, PRECISION, SAFEGUARD, SPECIPAN, SPLASHGUARD, and ULTRAMER are trademarks of Tyco Healthcare Group LP.


DOVER, MEASURIT, TEXAS CATHETER, and URI-DRAIN are trademarks of Sherwood Services AG. EASY-CATCH is a registered trademark of Superior Healthcare Group, Inc.

URIPROGRAM is a registered trademark of Proclinics S.A.

VELCRO is a registered trademark of Velcro Industries B.V.

# KENDALL Catheters

Offered with the following features:




### Smooth Reinforced Tips

— Provides added strength and ease of insertion.



### Large Smooth Drainage Eyes

— Provides easy insertion and uninhibited drainage flow.



### Reinforced Ribbed Balloons

— Provide superior strength and proper positioning within the bladder.




### Large Round Drainage Lumens

— Provide maximum draining and reduces the risk of clotting.




### An Assortment of Materials

- Natural Latex
- 100% Silicone
- Clear Vinyl
- Red Rubber



## Catheters

## Specialty Catheters



### CURITY Suprapubic Trocar Catheters

- A sharp edged smooth semi-open trocar tip with a clearly defined small cutting angle.
- Little or no bleeding in puncture channel; no stab incision of the cutis with scalpel.
- Minimum exertion required to pierce tissue.
- No tearing, punching and crushing of the tissue.
- 'Part open' design ensures ease of separation between the trocar and catheter.
- Catheter easily threads down the trocar shaft.
- Radiopaque tip to verify placement.
- Sterile.

Code	French Size	Ship Case
------	-------------	-----------

#### *Silicone Catheter & Trocar*

◆ 5040-02	10	20
◆ 5046-02	12	20

#### *Silicone Catheter Exchange Set*

◆ 5042-02	10	20
◆ 5044-02	12	20
◆ 5048-02	14	20

#### *Temperature Monitoring Catheter Set*


◆ 5041-02	10	20
◆ 5047-02	12	20

#### *Temperature Monitoring Catheter Exchange Set*

◆ 5043-02	10	20
◆ 5045-02	12	20
◆ 5049-02	14	20

#### *Male Connector Cables*

◆ 2001-02	Cable for SIEMENS Monitors	1
◆ 2002-02	Cable for HP Monitors	1



### CURITY Temperature Monitoring Silicone Catheters

- Allows urine drainage and temperature monitoring at the same time.
- Precision temperature probe located at tip of catheter.
- Reliable continuous measurement of core temperature (Closely parallels arterial blood temperature).
- Eliminates additional strain on patient of alternative temperature measuring techniques.
- Measurement confined to one place only (Bladder).
- Other means of temperature measurement (rectal, oesophageal, tympanic) are no longer required.
- Compatible with virtually all conventional patient monitors.
- Sterile.

Code	French Size	Balloon	Ship Case
------	-------------	---------	-----------

◆ 90053 T	8	3ml	12
◆ 90054 T	10	3ml	12
◆ 90055 T	12	5ml	12
◆ 90056 T	14	5ml	12
◆ 90051 T	16	5ml	12
◆ 90052 T	18	5ml	12

Code	Description	Ship Case
------	-------------	-----------

<i>THERMISTOR 400 - Female Connector Cables</i>		
◆ 502-0400	5m, with 1/4" jack plug	1
◆ 502-0411	3m, for use with HP Monitors	1
◆ 502-0401	5m, for use with HP Monitors	1
◆ 502-0415	3m, for use with SIEMENS SC Monitors	1
◆ 502-0405	5m, for use with SIEMENS SC Monitors	1
◆ 502-0410	3m, for use with SIEMENS Monitors	1

◆ Latex Free.



### DOVER 100% Silicone Foley Catheters

- Larger lumen than ordinary latex foley catheters.
- Reinforced, ribbed balloon.
- Non-reactive, non-adhering.
- Colour-coded inflation valve accepts luer lock and luer slip syringes.
- Packaged in an outer peel-back pouch and an easy-to-open inner sheath.
- Sterile.

Code 8887-	French Size	Length	Ship Case
<i>MALE Catheters</i>			
<i>10ml Balloon, 2-way</i>			
♦ 605122	12	43cm	10
♦ 605148	14	43cm	10
♦ 605163	16	43cm	10
♦ 605189	18	43cm	10
♦ 605205	20	43cm	10
♦ 605221	22	43cm	10
♦ 605247	24	43cm	10
♦ 605262	26	43cm	10
<i>30ml Balloon, 2-way</i>			
♦ 630161	16	43cm	10
♦ 630187	18	43cm	10
♦ 630203	20	43cm	10
♦ 630229	22	43cm	10
♦ 630245	24	43cm	10
♦ 630260	26	43cm	10
<i>30ml Balloon, 3-way</i>			
♦ 665183	18	43cm	10
♦ 665209	20	43cm	10
♦ 665225	22	43cm	10
♦ 665241	24	43cm	10
<i>FEMALE Catheters</i>			
<i>5ml Balloon, 2-way</i>			
♦ 615121	12	23cm	10
♦ 615147	14	23cm	10
♦ 615162	16	23cm	10
♦ 615188	18	23cm	10
<i>PAEDIATRIC Catheters</i>			
<i>3ml, 2-way</i>			
♦ 603085	8	28cm	10
♦ 603101	10	28cm	10

♦ Latex Free.

## Catheters

# KENDALL



### ULTRAMER Foley Catheters

- Latex foley catheter with ULTRAMER coating.
- ULTRAMER coating provides superior lubricity, increasing biocompatibility and potentially decreasing the incidence of catheter-related UTI.
- 4.2 times more lubricious than latex catheters and 4.6 times more lubricious than typical 100% silicone catheters.
- Maintains lubricity by absorbing mucosal fluid, providing easy insertion and removal with enhanced patient protection.
- Greater patient comfort.
- Ideal for long-term use.
- Sterile.

## Catheters

Code	French Size	Length	Ship Case
------	-------------	--------	-----------

#### MALE Catheters

##### 10ml 2-Way

1612	12	40cm	12
1614	14	40cm	12
1616	16	40cm	12
1619	18	40cm	12
1621	20	40cm	12
1622	22	40cm	12
1624	24	40cm	12
1626	26	40cm	12
1628	28	40cm	12
1631	30	40cm	12

##### 30ml 2-Way

1416	16	40cm	12
1419	18	40cm	12
1420	20	40cm	12
1422	22	40cm	12
1424	24	40cm	12
1426	26	40cm	12
1427	28	40cm	12
1431	30	40cm	12

##### 10ml 3-Way

2718	18	40cm	12
2720	20	40cm	12
2722	22	40cm	12
2723	24	40cm	12

##### 30ml 3-Way

2819	18	40cm	12
2821	20	40cm	12
2822	22	40cm	12
2824	24	40cm	12
2826	26	40cm	12

#### FEMALE Catheters

##### 10ml 2-Way

16128	12	23cm	12
16148	14	23cm	12
16168	16	23cm	12
16198	18	23cm	12
16218	20	23cm	12
16228	22	23cm	12

##### 30ml 2-Way

14168	16	23cm	12
14198	18	23cm	12
14208	20	23cm	12
14228	22	23cm	12

#### ULTRAMER Polyurethane Paediatric Catheter 3ml 2-Way with Luer Lock Adapter

◆ 6006*	6	28cm	12
---------	---	------	----

◆ Latex Free.

\* Non-CE marked at time of print.



**CURITY Silicone Coated Latex Foley Catheters**

- Silicone coating provides a barrier between the latex catheter shaft and the patients urethra increasing biocompatibility.
- Large eyes holes and drainage lumen facilitates urine flow.
- Thick funnel helps prevent collapse when irrigating.
- Reinforced tip protects catheter from stylet puncture.
- Strong symmetric balloons enhances confidence in catheter reliability.
- For medium/long-term use.
- Sterile.

Code	French Size	Length	Ship Case
<i>3ml Balloon, 2-way</i>			
2908-02	8	28cm	12
2910-02	10	28cm	12
<i>5ml Balloon, 2-way</i>			
1125-02	12	40cm	12
1185-02	14	40cm	12
1213-02	16	40cm	12
1394-02	18	40cm	12
1440-02	20	40cm	12
1575-02	22	40cm	12
1620-02	24	40cm	12
1746-02	26	40cm	12
<i>5ml Balloon, 2-way, TIEMANN type</i>			
9812-02	12	40cm	12
9814-02	14	40cm	12
9816-02	16	40cm	12
9818-02	18	40cm	12
9820-02	20	40cm	12
9822-02	22	40cm	12
9824-02	24	40cm	12
<i>30ml Balloon, 3-way</i>			
9718-02	18	40cm	12
9720-02	20	40cm	12
9722-02	22	40cm	12
9724-02	24	40cm	12



**DOVER T.U.R.P. Silicone Elastomer Coated Latex Foley Catheters**


- Silicone Elastomer-coated latex with nylon spiral embedded in catheter wall.
- Extra strong 30-50ml reinforced ribbed balloon can be safely inflated up to 90ml.
- Four large eyes.
- Colour-coded inflation valve accepts luer lock and luer slip syringes.
- Packaged in an outer peel-back pouch and an easy-to-open inner sheath.
- Sterile.

Code 8887-	French Size	Balloon Size	Type	Ship Case
648221*	22	30-50ml	2-way	10
648247*	24	30-50ml	2-way	10
649229*	22	30-50ml	3-way	10
649245*	24	30-50ml	3-way	10

Catheters

\* Non-CE marked at time of print.


## Catheters & Urethral Catheter Trays



### CURITY Intermittent Urethral Catheters

- Straight, thermosensitive PVC catheter.
- Soft, rounded end.
- Staggered oval eyes x 2.
- With colour-coded connector.
- Sterile.

Code	French Size	Length	Carton	Ship Case
<i>CURITY Male Catheter</i>				
♦ NT0620	6	40cm	25	675
♦ NT0820	8	40cm	25	675
♦ NT1020	10	40cm	25	675
♦ NT1220	12	40cm	25	675
♦ NT1420	14	40cm	25	675
♦ NT1620	16	40cm	25	675
♦ NT1820	18	40cm	25	675
♦ NT2020	20	40cm	25	675



### Nelaton Intermittent Urethral Catheters

- Straight, thermosensitive PVC catheter.
- Soft, rounded end.
- Staggered oval eyes x 2.
- Funnel connector.
- Sterile.

Code 1180-	French Size	Length	Ship Case
<i>Nelaton Male Catheter</i>			
♦ 836080	8	40cm	100
♦ 836106	10	40cm	100
♦ 836122	12	40cm	100
♦ 836148	14	40cm	100
♦ 836163	16	40cm	100
♦ 836189	18	40cm	100
♦ 836205	20	40cm	100
<i>Nelaton Female Catheter</i>			
♦ 835124	12	20cm	100
♦ 835140	14	20cm	100
♦ 835165	16	20cm	100



### DOVER Urethral Catheterisation Trays (Open)

#### Sterile contents:

- Underpad, waterproof drape, lubricating jelly, PVI solution, prepping balls, forceps, latex gloves and specimen container. (CSR wrap).

Code 8887-	Description	Ship Case
600065*	No Catheter	20
600073*	14Fr Clear Vinyl Catheter	20
600057*	14Fr Red Rubber Catheter	20



### DOVER Urethral Catheterisation Trays (Closed)

#### Sterile contents:

- Underpad, waterproof drape, lubricating jelly, PVI solution, prepping balls, forceps, latex gloves, specimen container and 1500ml collection bag. (CSR wrap).

Code 8887-	Description	Ship Case
<i>Tray with Catheter and Preconnected Drain Bag</i>		
600016*	14Fr Clear Vinyl Catheter	20
600032*	14Fr Red Rubber Catheter	20

♦ Latex Free.

\* Non-CE marked at time of print.



**DOVER Foley Catheter Insertion Trays**

**Sterile Contents:**

- Underpad, latex gloves, fenestrated drape, PVI solution (10%), prepping balls, forceps, lubricating jelly, specimen container, graduated basin. (CSR wrap).

Code 8887-	Description	Ship Case
601220*	10ml Prefilled Syringe	20
601261*	30ml Prefilled Syringe - Lidded tray	20



**DOVER Foley Catheterisation Trays**

**Sterile Contents:**

- Underpad, latex gloves, fenestrated drape, PVI solution (10%), prepping balls, forceps, lubricating jelly, specimen container, graduated basin, 10ml prefilled syringe. (CSR wrap).

Code 8887-	Catheter Size	Ship Case
<i>DOVER 100% Silicone Catheter</i>		
600776*	16Fr, 5ml	20
600784*	18Fr, 5ml	20
<i>DOVER Silicone Coated Latex Catheter</i>		
607169*	16Fr, 5ml	20
607185*	18Fr, 5ml	20

Foley Trays



**CURITY ULTRAMER Foley Catheterisation Trays**

**Sterile Contents:**

- Underpad, latex gloves, fenestrated drape, PVI solution (10%), prepping balls, forceps, lubricating jelly, specimen container, graduated basin, 10ml prefilled syringe. (CSR wrap).


Code	Catheter Size	Ship Case
2101*	16Fr, 5ml	20
3044*	18Fr, 5ml	20

\* Non-CE marked at time of print.

# The 4 Critical Sites for Bacteria to Enter the Urinary Tract<sup>1</sup>

Kendall offers 4 products/features to reduce the incidence of infection introduction:

1. Catheter/Meatal Junction
  - ULTRAMER Foley Catheters
2. Catheter/Tubing Junction
  - TAMPER EVIDENT SEAL
3. Tubing/Drainage Bag Junction
  - MONO-FLO Anti-Reflux Device
4. Drain Spout
  - SPLASHGUARD




**1**  
ULTRAMER  
Foley Catheters



**2**  
TAMPER EVIDENT  
SEAL



**4**  
SPLASHGUARD



**3**  
MONO-FLO  
Anti-Reflux Device



“Wash Your Hands”  
Logo

<sup>1</sup> Calvin M. Kunin, *Detection, Prevention and Management of Urinary Tract Infections, Edition 4, Lea & Febiger, Philadelphia, 1987, pp 245-297.*

# KENDALL Drainage Bags, Meters & Trays

Offered with the following features:



## MONO-FLO Anti-Reflux Device (Latex Free)

- Helps prevent urine backflow.
- Unique one-way flap valve.
- Prevents urine reflux even if bag is accidentally laid on its side or tipped upside down.



## SPLASHGUARD II (Latex Free)

- Helps prevent backsplash of urine.
- Protects provider and patient.
- Prevents drain tube from direct contact with the ground.



## TAMPER EVIDENT SEAL

- Helps secure Foley catheter to drainage tubing junction.
- Ensures integrity of closed system.



## SAFEGUARD Needleless Sampling Port (Latex Free)

- Accessed with standard luer slip syringe or blunt cannula.
- Pre-slit "latex free" elastomer rubber.



## Standard Prep Tray (Latex Free)

- Sterile contents: Fenestrated drape, underpad, **latex free** gloves, forceps, prepping balls, 10ml prefilled sterile water inflation syringe, unique bellows lubricant syringe, Povidone-Iodine solution, specimen container and label. Catheter pre-connected to drain bag or urine meter. (CSR wrap).

## Drainage Bags



### CURITY MONO-FLO Premium Drainage Bags

- CLOSED SYSTEM.
- MONO-FLO Anti-Reflux valve automatically closes when bag is tilted, lifted, or laid on side.
- Large diameter drip chamber provides dry barrier against bacteria, facilitates passage of clots.
- Universal Hanger.
- Large diameter non-kink tube.
- 2000ml capacity.
- With needle sampling port as standard on all bags.
- SPLASHGUARD II prevents backsplash of urine.
- SAFEGUARD Needleless Sampling Port prevents needle stick injury.
- Sterile.



Code	Description	Ship Case
6208-02	With 90cm Drainage Tubing	30
6209-02	With 120cm Drainage Tubing	30
6208	With SPLASHGUARD II drain spout & 90cm Drainage Tubing	20
6308	With SPLASHGUARD II drain spout, SAFEGUARD Needleless Sampling Port & 120cm Drainage Tubing	20



**DOVER Anti-Reflux Drainage Bags**

- Compact teardrop shape; 2000ml capacity.
- Centered inlet and exit drains.
- Exclusive FLO-CHEK top line anti-reflux valve at catheter connection.
- 135cm drainage tubing.
- Vented catheter connector.
- KRATON sample port at top of line.
- One-piece hinged hook and a string/capstan hanger for variable hanging heights.
- Sterile.

Code	Description	Ship Case
8887-		
601709	With 7.5mm I.D. tubing, drip chamber, and Forman bottom drain	40
601600	With 7.0mm I.D. tubing and latex bottom drain tube with Aggogle Clamp	40
600602	With 7.0mm I.D. tubing, vented connector and sample port with Aggogle Clamp	40

**KENGUARD Drain Bags**

- Contoured design for more even filling and complete emptying.
- 8.5mm drainage tube provides quick drainage, eliminates standing column of urine.
- Universal hanger fits securely on bedrails and wheelchairs.
- Convenient carrying handle for ambulatory patients.
- Sterile.

Code	Description	Ship Case
3512	2000ml bag with Anti-Reflux chamber	20
3502	2000ml bag without Anti-Reflux chamber	20

Drainage Bags

## Drainage Bags



UB3990



UB3994/UB3996



UB3010



UB3020/UB3030



UB3101/  
UB3102/UB3103



UB6207

### CURITY Economy Drain Bags

- 2000ml capacity.
- Terraced catheter connector.
- Anti-reflux flap valve.
- Rectangular shaped bag (except UB6207).

Code	Description	Ship Case
<i>Sterile</i>		
◆ UB3030	Pull/push bottom drain, 90cm drain tube	160
◆ UB3101*	Pull/push bottom drain, rigid plastic hanger bar & string loop, 90cm drain tube	200
◆ UB3103*	Pull/push bottom drain, rigid plastic hanger bar & string loop, 120cm drain tube	200
◆ UB3102*	Pull/push bottom drain, rigid plastic hanger bar & plastic loops, 90cm drain tube	200
◆ UB3010	Pull/push bottom drain, plastic loop hangers, 90cm drain tube	160
◆ UB3020	Pull/push bottom drain, plastic loop hangers, 120cm drain tube	150
◆ UB3996*	Aspirating port & twist bottom drain, no hangers, 120cm drain tube, white backing	200
UB6207	Tear Shaped bag with pull/push bottom drain, latex sample port, rigid plastic hanger bar & adjustable plastic strap, 120cm wide-bore drain tube and partial white backing	100
<i>Non-sterile bulk</i>		
◆ UB3990*	No bottom drain, 90cm drain tube, white backing	250
◆ UB3994*	Aspirating port, twist bottom drain, 120cm drain tube, white backing	250

◆ Latex Free.

\* Non-CE marked at time of print.



**CURITY Ureteral Drainage Bags**

- Integrated collection system. Double-bag keeps ureteral and bladder urine separate.
- Ureteral catheter adapter accepts 4Fr - 8Fr catheters, I.V. and feeding tubes.
- Breakaway safety feature disconnects ureteral catheter when pulled.
- Closed-system aspirating port permits repeated fresh urine ureteral catheter irrigation port may also be used for specimen collection.
- Universal hanging device fits easily over bedrail, serves as carrying handle.
- Sterile.

Code	Description	Ship Case
6251	Ureteral Drainage Set	10



**CURITY T.U.R.P. Drainage Bag**

- CLOSED SYSTEM.
- MONO-FLO Anti-Reflux valve automatically closes when bag is tilted, lifted, or laid on side.
- Large diameter drip chamber provides dry barrier against bacteria, facilitates passage of clots.
- Strong Hanger.
- Large diameter non-kink tubing.
- 4000ml capacity.
- 100ml graduation marks.
- Sterile.

Code	Description	Ship Case
10082	T.U.R.P. Drainage Bag	25

Drainage Bags



CURITY Foley Tray with 6308 Drain Bag

**Sterile Contents:**

- Standard Prep Tray.
- 2000ml MONO-FLO Drain Bag with SAFEGUARD Needleless Sampling Port & SPLASHGUARD II Drain Spout.
- Catheter preconnected to drain tube.

Code	Description	Ship Case
6165*	16Fr ULTRAMER Catheter	10
6166*	18Fr ULTRAMER Catheter	10
◆ 6146*	16Fr 100% Silicone Catheter	10
◆ 6148*	18Fr 100% Silicone Catheter	10



CURITY Foley Tray with 6209 Drain Bag

**Sterile Contents:**

- Standard Prep Tray.
- 2000ml MONO-FLO Drain Bag with MONO-FLO Anti-reflux device.
- Catheter preconnected to drain tube.

Code	Description	Ship Case
6014*	16Fr ULTRAMER Catheter	10
6080*	18Fr ULTRAMER Catheter	10
6140*	16Fr 100% Silicone Catheter	10
6210*	18Fr 100% Silicone Catheter	10

## Drainage Bag Foley Trays



DOVER Foley Tray with 3512 Drain Bag

Sterile Contents:

- Standard Prep Tray.
- 2000ml KENGUARD Drain Bag with Anti-reflux valve.
- Catheter preconnected to drain tube.

Code 8887-	Description	Ship Case
670662*	16Fr Silicone Coated Catheter	20
670688*	18Fr Silicone Coated Catheter	20
640665*	16Fr 100% Silicone Catheter	20
640681*	18Fr 100% Silicone Catheter	20



DOVER Foley Tray with Drain Bag


Sterile Contents:

- Standard Prep Tray.
- 2000ml DOVER Drain Bag.
- Catheter preconnected to drain tube.

Code 8887-	Description	Ship Case
<i>With DOVER Drain Bag 8887-601600</i>		
671660*	16Fr Silicone Coated Catheter	20
671686*	18Fr Silicone Coated Catheter	20
641663*	16Fr 100% Silicone Catheter	20
641689*	18Fr 100% Silicone Catheter	20
<i>With DOVER Drain Bag 8887-601709</i>		
671769*	16Fr Silicone Coated Catheter	20
671785*	18Fr Silicone Coated Catheter	20
641762*	16Fr 100% Silicone Catheter	20
641788*	18Fr 100% Silicone Catheter	20

## Drainage Bag Foley Trays

\* Non-CE marked at time of print.



#### URIPROGRAM 500 Urine Meter

- Transparent triple chamber urine meter: 100ml, 100ml to 300ml, 300ml to 500ml.
- Sequential drainage tap.
- PVC drainage tubing with clamp.
- Catheter connector with sampling port.
- Anti-reflux valve.
- Hydrophobic air filter.
- Easy-adjustable bed rail straps.
- 1500ml capacity available with or without bottom drain.
- Sterile.

Code 1180-	Description	Ship Case
<i>With 115cm Drainage Tubing</i>		
400101	With Closed Bag	10
400143	With Drainage Bag (attached) and Sequential Drainage Tap	10
400119	With Drainage Bag and Sequential Drainage Tap	10
<i>With 90cm Drainage Tubing</i>		
400200	With Closed Bag	10
400218	With Drainage Bag	10
<i>Replacement Closed Bags</i>		
400846	Sterile	100
400861	With Sequential Drainage Tap	100
400854	Non-Sterile	100
400879	Non-Sterile with Sequential Drainage Tap	100



#### PRECISION 400 Urine Meters

- SPLASHGUARD II Drain Spout prevents splash-back of urine.
- SAFEGUARD Needleless Sampling Port prevents needle stick injury.
- Spring-loaded meter sample port.
- Small burette chamber.
- White meter background.
- Extra-long 150cm tubing.
- Frosted writing area.
- Teardrop-shaped bag.
- Sterile.

Code	Description	Ship Case
2000	With SPLASHGUARD II	10
♦ 2000LF	Latex Free, with SPLASHGUARD II	10
7000	With SPLASHGUARD II and SAFEGUARD Needleless Sampling Port	10

## Urine Meters



**DOVER MEASURIT Urine Meters**

- Transparent dual-chamber urine meter (30ml burette chamber).
- Rotating valve on burette chamber.
- 1800ml capacity drainage bag.
- Vented catheter connector with plastic-backed latex sample port and 7.5mm I.D. tubing.
- 510cm Forman bottom drain.
- One-piece hinged hook with string for variable hanging heights.
- Sterile.

Code 8887-	Description	Ship Case
600800	200ml Urine Meter	10
600891	400ml Urine Meter	10

**CURITY 200 Urine Meters**

- Economic Hygiene Concept.
- 2000ml capacity drainage bag.
- 200ml capacity meter.
- Clear plastic meter calibrated in 1ml increments (2-40ml) and 5ml increments (40-2000ml) for easy reading.
- Transfer tube with anti-reflux flap between meter and bag.
- Anti-reflux valve at drain tube.
- Kink resistant drain tubing.
- Sterile.

Code	Description	Ship Case
8000	With spring loaded sampling port	10
6500	With 120cm drainage tubing	10
6502	With 150cm drainage tubing	10
400010*	Without Anti-reflux valve at drain tube	10

Urine Meters

\* Non-CE marked at time of print.



PRECISION 400 Urine Meter Foley Trays

**Sterile Contents:**

- Standard Prep Tray.
- PRECISION 400ml Meter with SAFEGUARD Needleless Sampling Port & SPLASHGUARD II Drain Spout.
- Catheter preconnected to drain tube.

Code	Description	Ship Case
♦ 2006LF*	16Fr 100% Silicone Catheter	10
♦ 2008LF*	18Fr 100% Silicone Catheter	10




DOVER MEASURIT Urine Meter Foley Trays

**Sterile Contents:**

- Standard Prep Tray.
- 400ml MEASURIT Meter.
- Catheter preconnected to drain tube.

Code	Description	Ship Case
8887-	400ml DOVER Meter	
640962*	16Fr 100% Silicone Catheter	10
640988*	18Fr 100% Silicone Catheter	10

## Urine Meter Foley Trays



**CURITY Urodynamic Catheter**

- Available with double or triple lumens.
- FILL line, CMG line and UPP line have unique colour-coded luer-lock hubs.
- French sizes are clearly marked.
- Flexible polyurethane becomes more pliable as it warms to body temperature.
- Rounded tip and flexible polyurethane construction allow for smoother insertion.
- Radiopaque markers are uniquely positioned inside the catheter.
- Catheter is conveniently marked in centimetres.

Code	French Size	Ship Case
<i>Dual Lumen Catheters</i>		
◆ 6427*	7	12
◆ 6429*	10	12
<i>Triple Lumen Catheters</i>		
◆ 6437*	7	12
◆ 6439*	10	12



**CURITY KIDS 150ml Urine Meter**

- Pull-push emptying and sample ports.
- Injection port for fresh sampling.
- Universal hanger.
- Anti-reflux valve.
- Narrow-bore tubing.
- Luer-lock or luer-slip connection.

Code	Description	Ship Case
6600	PRECISION 150ml Urine Meter with Anti-Reflux Device, Luer-Lock Catheter Connector, Top Line Filter for Drainage and Universal Hanging Device	10

Paediatric  
Urology  
Products



**CURITY KIDS Drainage Bag**

- CLOSED SYSTEM.
- MONO-FLO Anti-reflux valve.
- Large diameter drip chamber.
- 1000 ml capacity.
- Luer lock catheter connector.
- Top Line Filter for drainage.

Code	Description	Ship Case
6510	With Universal Hanging Device	20
<i>ULTRAMER Polyurethane Paediatric Catheter</i>		
◆ 6006*	6Fr 3ml 2-way with Luer Lock adapter	12

◆ Latex Free.

\* Non-CE marked at time of print.

## Urinary Collection



### MONO-FLO TAG & NACHT Leg Bags

- CLOSED SYSTEM.
- Drainage Bag for the mobile catheterised patient.
- Bag can be continuously used for up to 14 days with minimal risk of UTI.
- Soft viscose backing for patient comfort.
- Discreet, even with MONO-FLO Anti-reflux dry drip chamber.
- 1000ml capacity.
- Sterile.

Code	Description	Ship Case
5079	With 30cm drainage tube	20
5081	With 60cm drainage tube	20
5083	With 90cm drainage tube	20
<i>Adjustable Leg Straps - required for TAG &amp; NACHT</i>		
31377	With Velcro tags	10



### CURITY Urine Leg Bags

- Anti-Reflux valve.
- Longer catheter adapter for greater security with all types of catheters.
- Twist valve for easier emptying and less spill; no bottom cap to lose or misplace.
- Soft, cloth straps.
- Sterile.

Code	Size	Capacity	Ship Case
3433	Large, with 30cm extension tube	1000ml	50
145514	Medium	500ml	50
145516	Large	1000ml	50



**TEXAS External Male Catheters**

- Soft, lightweight latex.
- Drainage connector permits attachment to connective tubing or to a collection bag.
- Optional self-adhering elastic foam strap for holding catheter in position.
- Size – standard.

Code 8884-	Description	Ship Case
<i>Sterile</i>		
730300	Individually packaged catheter with 25mm DOVER adhesive elastic foam strap	144
<i>Non-Sterile</i>		
730200	Bulk packaged catheters – method of attachment required	144
♦ 730600	Bulk packaged 25mm DOVER adhesive elastic foam straps	500

**DOVER Catheter Straps**

- Holds catheter securely in place without stressing or irritating skin.

Code 8887-	Description	Ship Case
♦ 600149	DOVER catheter leg strap	10

**URI-DRAIN Urinary Extension Tubing with Connector**

- For use with all external devices and leg bags.
- PVC tubing can be cut for use between catheter and bag.
- Permanently attached connector.

Code 8884-	Description	Ship Case
731900*	URI-DRAIN extension tubing with connector, 8.5mm x 45cm	24

Urinary  
Collection

♦ Latex Free.

\* Non-CE marked at time of print.



### Vaginal Irrigation (Douche) Bag

- Easy to hang and easy to read.
- Polyethylene.
- Sterile.

Code	Description	Ship Case
145617*	1500ml Bag with Flip Top slide lock closure, attached 150cm tube with shut-off clamp, douche tip and douche powder	50



### Enema Bags

- 1500ml Polyethylene Bag with Flip Top Slide Lock Closure.

Code	Description	Ship Case
<i>Non-Sterile</i>		
♦ 145540*	With attached 24Fr, 150cm enema tube, depth mark 13cm from pre-lubricated tip, castile soap and underpad	50

## Bladder Management



### Flatus Bag, Rectal Tube

- Smooth, rounded tips for easy insertion.

Code	Description	Ship Case
<i>Non-Sterile Vinyl Flatus Bag</i>		
♦ 145524	Complete with 24Fr, 48cm length rectal tube with pre-lubricated tip	50
<i>Rectal Tube</i>		
♦ 155731	24Fr, 50cm length, pre-lubricated tip	50



### Enema Bucket

Code	Description	Ship Case
<i>Non-Sterile</i>		
♦ 145546	1400ml plastic bucket, 24Fr, 150cm enema tube, depth mark 13cm from pre-lubricated tip, castile soap and underpad	50



Commode Specimen Collection Systems

- Ideal collection unit for urine and stool samples or kidney stones.
- Wide lip and firm base helps prevent splashing, spilling, and tipping.
- Clearly marked gradients are easy to read.

Code	Description	Ship Case
◆ 4014	CURITY SPECIPAN with clearly marked gradients	100
◆ 2400SA	Commode Specimen Measuring System with Lid	60
◆ 2500SA	Commode Specimen Measuring System	100



SAGE Midstream Catch Kits

- Sterile midstream collection kits feature a moisture resistant plastisol liner that minimizes chance of leakage.
- Clear packaging for easy visualisation.
- Sequentially packed.
- Pneumatic capable.

Code	Description	Ship Case
◆ 2090SA	Midstream Kit with Castile Soap Towelettes	100
◆ 4090SA	Midstream Kit with Protective Ring and Castile Soap Towelettes	100
◆ 2040SA	Midstream Kit with Funnel and Castile Soap Towelettes	100
◆ 2110SA	Calculi Strainer	100

Specimen Collection



EASY-CATCH Midstream Series

- Features:
  - 133ml graduated Specimen Container with rim protector.
  - Perfect-seal Screw Cap.
  - Protective Cover.
  - Patient Label.
- Sterile fluid path.

Code	Description	Ship Case
◆ 25000	With BZK towelettes	100
◆ 25015	With Castile Soap towelettes	100

◆ Latex Free.



CURITY Clean Catch Kit

- 4.5oz capacity.
- Screw cap.

Code	Description	Ship Case
◆ 5205	Castile Soap Towelettes	48
◆ 5204	BZK Towelettes	48

Code	Product	Catalogue	Page
10082	T.U.R.P. Drainage Bag w/Vented Connector, 4000ml	UR	13
1180-400101	URIPROGRAM Urine Meter w/Closed Bag & 115cm Tubing	UR	16
1180-400119	URIPROGRAM Urine Meter w/Drainage Bag & 115cm Tubing	UR	16
1180-400143	URIPROGRAM Urine Meter w/Drainage Bag & 115cm Tubing	UR	16
1180-400200	URIPROGRAM Urine Meter w/Closed Bag & 90cm Tubing	UR	16
1180-400218	URIPROGRAM Urine Meter w/Drainage Bag & 90cm Tubing	UR	16
1180-400846	URIPROGRAM Replacement Closed Bag, Sterile	UR	16
1180-400854	URIPROGRAM Replacement Closed Bag, Non-Sterile	UR	16
1180-400861	URIPROGRAM Replacement Closed Bag, Sterile	UR	16
1180-400879	URIPROGRAM Replacement Closed Bag, Non-Sterile	UR	16
1180-835124	Nelaton Urethral Cath. w/Funnel, PVC, 12Fr, 20cm	UR	6
1180-835140	Nelaton Urethral Cath. w/Funnel, PVC, 14Fr, 20cm	UR	6
1180-835165	Nelaton Urethral Cath. w/Funnel, PVC, 16Fr, 20cm	UR	6
1180-836064	Nelaton Urethral Cath., PVC, 6Fr, 40cm	UR	6
1180-836080	Nelaton Urethral Cath., PVC, 8Fr, 40cm	UR	6
1180-836106	Nelaton Urethral Cath., PVC, 10Fr, 40cm	UR	6
1180-836122	Nelaton Urethral Cath., PVC, 12Fr, 40cm	UR	6
1180-836148	Nelaton Urethral Cath., PVC, 14Fr, 40cm	UR	6
1180-836163	Nelaton Urethral Cath., PVC, 16Fr, 40cm	UR	6
1180-836189	Nelaton Urethral Cath., PVC, 18Fr, 40cm	UR	6
1180-836205	Nelaton Urethral Cath., PVC, 20Fr, 40cm	UR	6
1185-02	CURITY Foley Cath., Sil. Coated Latex, 5ml, 2-way, 14Fr, 40cm	UR	5
1213-02	CURITY Foley Cath., Sil. Coated Latex, 5ml, 2-way, 16Fr, 40cm	UR	5
1216	ULTRAMER Foley Cath., Female, 30ml, 3-way, 16Fr, 23cm	UR	4
1218	ULTRAMER Foley Cath., Female, 30ml, 3-way, 18Fr, 23cm	UR	4
1220	ULTRAMER Foley Cath., Female, 30ml, 3-way, 20Fr, 23cm	UR	4
1222	ULTRAMER Foley Cath., Female, 3-way, 22Fr, 23cm	UR	4
1394-02	CURITY Foley Cath., Sil. Coated Latex, 5ml, 2-way, 18Fr, 40cm	UR	5
1416	ULTRAMER Foley Cath., Male, 30ml, 2-way, 16Fr, 40cm	UR	4
1419	ULTRAMER Foley Cath., Male, 30ml, 2-way, 18Fr, 40cm	UR	4
1420	ULTRAMER Foley Cath., Male, 30ml, 2-way, 20Fr, 40cm	UR	4
1422	ULTRAMER Foley Cath., Male, 30ml, 2-way, 22Fr, 40cm	UR	4
1424	ULTRAMER Foley Cath., Male, 30ml, 2-way, 24Fr, 40cm	UR	4
1426	ULTRAMER Foley Cath., Male, 30ml, 2-way, 26Fr, 40cm	UR	4
1427	ULTRAMER Foley Cath., Male, 30ml, 2-way, 28Fr, 40cm	UR	4
1431	ULTRAMER Foley Cath., Male, 30ml, 2-way, 30Fr, 40cm	UR	4
1440-02	CURITY Foley Cath., Sil. Coated Latex, 5ml, 2-way, 20Fr, 40cm	UR	5
145514	CURITY Urine Leg Bag, Medium, 500ml	UR	20
145516	CURITY Urine Leg Bag, Large, 1000ml	UR	20
145524	Flatus Bag, Vinyl, w/24Fr, 48cm Pre-Lub. Rectal Tube	UR	22
145540	Enema Bag, 15000ml, w/attached 24Fr, 150cm Enema Tube, Non-Sterile	UR	22
145546	Enema Bucket, Plastic, 1400ml, w/24Fr, 150cm Enema Tube, Non-Sterile	UR	22
145617	Vaginal Irrigation (Douche) Bag, 1500ml	UR	22
155731	Rectal Tube, 24Fr, 50cm, Pre-Lub.	UR	22
1575-02	CURITY Foley Cath., Sil. Coated Latex, 5ml, 2-way, 22Fr, 40cm	UR	5
1612	ULTRAMER Foley Cath., Male, 10ml, 2-way, 12Fr, 40cm	UR	4
1614	ULTRAMER Foley Cath., Male, 10ml, 2-way, 14Fr, 40cm	UR	4
1616	ULTRAMER Foley Cath., Male, 10ml, 2-way, 16Fr, 40cm	UR	4
1619	ULTRAMER Foley Cath., Male, 10ml, 2-way, 18Fr, 40cm	UR	4
1620-02	CURITY Foley Cath., Sil. Coated Latex, 5ml, 2-way, 24Fr, 40cm	UR	5
1621	ULTRAMER Foley Cath., Male, 10ml, 2-way, 20Fr, 40cm	UR	4
1622	ULTRAMER Foley Cath., Male, 10ml, 2-way, 22Fr, 40cm	UR	4
1624	ULTRAMER Foley Cath., Male, 10ml, 2-way, 24Fr, 40cm	UR	4
1626	ULTRAMER Foley Cath., Male, 10ml, 2-way, 26Fr, 40cm	UR	4
1628	ULTRAMER Foley Cath., Male, 10ml, 2-way, 28Fr, 40cm	UR	4
1631	ULTRAMER Foley Cath., Male, 10ml, 2-way, 30Fr, 40cm	UR	4
1712	ULTRAMER Foley Cath., Female, 10ml, 2-way, 12Fr, 23cm	UR	4
1714	ULTRAMER Foley Cath., Female, 10ml, 2-way, 14Fr, 23cm	UR	4
1716	ULTRAMER Foley Cath., Female, 10ml, 2-way, 16Fr, 23cm	UR	4
1718	ULTRAMER Foley Cath., Female, 10ml, 2-way, 18Fr, 23cm	UR	4
1720	ULTRAMER Foley Cath., Female, 10ml, 2-way, 20Fr, 23cm	UR	4
1722	ULTRAMER Foley Cath., Female, 10ml, 2-way, 22Fr, 23cm	UR	4
1746-02	CURITY Foley Cath., Sil. Coated Latex, 5ml, 2-way, 26Fr, 40cm	UR	5
2000	CURITY Precision 400 Urine Meter	UR	16
2000LF	CURITY Precision 400 Urine Meter, Latex Free	UR	16

## Index by Code Number

Code	Product	Catalogue	Page
2001-02	Male Connector Cable for SIEMENS Monitors	UR	2
2002-02	Male Connector Cable for HP Monitors	UR	2
2006LF	PRECISION 400 Urine Meter Foley Tray w/100% Sil. Cath., 16Fr	UR	18
2008LF	PRECISION 400 Urine Meter Foley Tray w/100% Sil. Cath., 18Fr	UR	18
2040SA	SAGE Midstream Catch Kit	UR	23
2090SA	SAGE Midstream Catch Kit	UR	23
2101	CURITY ULTRAMER Catheterisation Tray w/ULTRAMER Cath., 16Fr	UR	7
2110SA	SAGE Calculi Strainer	UR	23
24000SA	Commode Spec. Meas. Collection w/Lid	UR	23
25000	EASY-CATCH Midstream Series 25000 Spec. Cont., 133ml	UR	23
25000SA	Commode Spec. Meas. Collection	UR	23
25015	EASY-CATCH Midstream Series 25000 Spec. Cont., 133ml	UR	23
2718	ULTRAMER Foley Cath., Male, 10ml, 3-way, 18Fr, 40cm	UR	4
2720	ULTRAMER Foley Cath., Male, 10ml, 3-way, 20Fr, 40cm	UR	4
2722	ULTRAMER Foley Cath., Male, 10ml, 3-way, 22Fr, 40cm	UR	4
2723	ULTRAMER Foley Cath., Male, 10ml, 3-way, 24Fr, 40cm	UR	4
2819	ULTRAMER Foley Cath., Male, 30ml, 3-way, 18Fr, 40cm	UR	4
2821	ULTRAMER Foley Cath., Male, 30ml, 3-way, 20Fr, 40cm	UR	4
2822	ULTRAMER Foley Cath., Male, 30ml, 3-way, 22Fr, 40cm	UR	4
2824	ULTRAMER Foley Cath., Male, 30ml, 3-way, 24Fr, 40cm	UR	4
2826	ULTRAMER Foley Cath., Male, 30ml, 3-way, 26Fr, 40cm	UR	4
2908-02	CURITY Foley Cath., Sil. Coated Latex, 3ml, 2-way, 8Fr, 28cm	UR	5
2910-02	CURITY Foley Cath., Sil. Coated Latex, 3ml, 2-way, 10Fr, 28cm	UR	5
3044	CURITY ULTRAMER Catheterisation Tray w/ULTRAMER Cath., 18Fr	UR	7
31377	Adjustable Leg Straps for TAG & NACHT Leg Bags	UR	20
3433	CURITY Urine Leg Bag w/30cm Ext. Tube, Large, 1000ml	UR	20
3502	KENGUARD Drainage Bag w/o Anti-Reflux, 2000ml	UR	11
3512	KENGUARD Drainage Bag w/Anti-Reflux, 2000ml	UR	11
400010	CURITY Urine Meter w/o Anti-Reflux, 200ml	UR	17
4014	CURITY SPECIPAN Collection Unit	UR	23
4090SA	SAGE Midstream Catch Kit	UR	23
502-0400	THERMISTOR 400 Female Connector Cable, 5m, w/¼" Jack Plug	UR	2
502-0401	THERMISTOR 400 Female Connector Cable, 5m, for HP Monitors	UR	2
502-0405	THERMISTOR 400 Female Connector Cable, 5m, for SIEMENS SC Monitors	UR	2
502-0410	THERMISTOR 400 Female Connector Cable, 3m, for SIEMENS Monitors	UR	2
502-0411	THERMISTOR 400 Female Connector Cable, 3m, for HP Monitors	UR	2
502-0415	THERMISTOR 400 Female Connector Cable, 3m, for SIEMENS SC Monitors	UR	2
5040-02	CURITY Suprapubic Trocar Cath., Silicone, 10Fr	UR	2
5041-02	CURITY Temperature Monitoring Cath. Set, 10Fr	UR	2
5042-02	CURITY Suprapubic Trocar Cath. Exchange Set, Silicone, 10Fr	UR	2
5043-02	CURITY Temperature Monitoring Cath. Exchange Set, 10Fr	UR	2
5044-02	CURITY Suprapubic Trocar Cath. Exchange Set, Silicone, 12Fr	UR	2
5045-02	CURITY Temperature Monitoring Cath. Exchange Set, 12Fr	UR	2
5046-02	CURITY Suprapubic Trocar Cath., Silicone, 12Fr	UR	2
5047-02	CURITY Temperature Monitoring Cath. Set, 12Fr	UR	2
5048-02	CURITY Suprapubic Trocar. Cath. Exchange Set, Silicone, 14Fr	UR	2
5049-02	CURITY Temperature Monitoring Cath. Exchange Set, 14Fr	UR	2
5079	MONO-FLO TAG & NACHT Leg Bag w/30cm Drain Tube	UR	20
5081	MONO-FLO TAG & NACHT Leg Bag w/60cm Drain Tube	UR	20
5083	MONO-FLO TAG & NACHT Leg Bag w/90cm Drain Tube	UR	20
5204	CURITY Clean Catch Kit, BZK Towelettes	UR	23
5205	CURITY Clean Catch Kit, Castile Soap Towelettes	UR	23
6006	ULTRAMER Paed. Foley Cath. w/L/Lock Adapter, Polyurethane, 6Fr	UR	4, 19
6014	CURITY Foley Tray w/6209 Drain Bag & ULTRAMER Cath., 16Fr	UR	14
6080	CURITY Foley Tray w/6209 Drain Bag & ULTRAMER Cath., 18Fr	UR	14
6140	CURITY Foley Tray w/6209 Drain Bag & 100% Sil. Cath., 16Fr	UR	14
6146	CURITY Foley Tray w/6308 Drain Bag & 100% Sil. Cath., 16Fr	UR	14
6148	CURITY Foley Tray w/6308 Drain Bag & 100% Sil. Cath., 18Fr	UR	14
6165	CURITY Foley Tray w/6308 Drain Bag & ULTRAMER Cath., 16Fr	UR	14
6166	CURITY Foley Tray w/6308 Drain Bag & ULTRAMER Cath., 18Fr	UR	14
6208	CURITY MONO-FLO Drain Bag w/Anti-Reflux, 2000ml	UR	10
6208-02	CURITY MONO-FLO Drain Bag w/Anti-Reflux, 2000ml	UR	10
6209-02	CURITY MONO-FLO Drain Bag w/Anti-Reflux, 2000ml	UR	10
6210	CURITY Foley Tray w/6209 Drain Bag & 100% Sil. Cath., 18Fr	UR	14
6251	CURITY Ureteral Drainage Set	UR	13

Index by  
Code Number

Code	Product	Catalogue	Page
6308	CURITY MONO-FLO Drain Bag w/Anti-Reflux & SAFEGUARD, 2000ml	UR	10
6427	CURITY Urodynamic Cath., D/Lumen, 7Fr	UR	19
6429	CURITY Urodynamic Cath., D/Lumen, 10Fr	UR	19
6437	CURITY Urodynamic Cath., T/Lumen, 7Fr	UR	19
6439	CURITY Urodynamic Cath., T/Lumen, 10Fr	UR	19
6500	CURITY Urine Meter w/120cm tubing, 200ml	UR	17
6502	CURITY Urine Meter w/150cm tubing, 200ml	UR	17
6510	CURITY Paed. Drainage Bag, 1000ml	UR	19
6600	PRECISION Paed. Urine Meter w/Anti-Reflux, 150ml	UR	19
7000	CURITY Precision 400 Urine Meter w/SAFEGUARD	UR	16
8000	CURITY Urine Meter, 200ml	UR	17
8884-730200	TEXAS Male External Cath. w/o Strap, Bulk, Non-Sterile	UR	21
8884-730300	TEXAS Male External Cath. w/25mm Foam Strap, Sterile	UR	21
8884-730600	DOVER Elastic Foam Strap, 25mm	UR	21
8884-731900	URI-DRAIN Urinary Ext. Tubing w/Connector, 8.5mm x 45cm	UR	21
8887-600016	DOVER Urethral Catheterisation Tray (Closed) w/Clear Vinyl Cath., 14Fr	UR	6
8887-600032	DOVER Urethral Catheterisation Tray (Closed) w/Red Rubber Cath., 14Fr	UR	6
8887-600057	DOVER Urethral Catheterisation Tray (Open) w/Red Rubber Cath., 14Fr	UR	6
8887-600065	DOVER Urethral Catheterisation Tray (Open) w/o Cath.	UR	6
8887-600073	DOVER Urethral Catheterisation Tray (Open) w/Vinyl Cath., 14Fr	UR	6
8887-600149	DOVER Catheter Leg Strap	UR	21
8887-600602	DOVER Drainage Bag w/Vented Connector, 2000ml	UR	11
8887-600776	DOVER Foley Catheterisation Tray w/100% Sil. Cath., 16Fr	UR	7
8887-600784	DOVER Foley Catheterisation Tray w/100% Sil. Cath., 18Fr	UR	7
8887-600800	DOVER MEASURIT Urine Meter, 200ml	UR	17
8887-600891	DOVER MEASURIT Urine Meter, 400ml	UR	17
8887-601220	DOVER Foley Cath. Insert. Tray w/10ml Syringe	UR	7
8887-601261	DOVER Foley Cath. Insert. Tray w/30ml Syringe, Lidded tray	UR	7
8887-601600	DOVER Anti-Reflux Drainage Bag w/FLO-CHEK, 2000ml	UR	11
8887-601709	DOVER Anti-Reflux Drainage Bag w/FLO-CHEK, 2000ml	UR	11
8887-603085	DOVER Foley Cath., 100% Silicone, Male, Paed. 3ml, 2-way, 8Fr, 28cm	UR	3
8887-603101	DOVER Foley Cath., 100% Silicone, Male, Paed. 3ml, 2-way, 10Fr, 28cm	UR	3
8887-605122	DOVER Foley Cath., 100% Silicone, Male, 10ml, 2-way, 12Fr, 43cm	UR	3
8887-605148	DOVER Foley Cath., 100% Silicone, Male, 10ml, 2-way, 14Fr, 43cm	UR	3
8887-605163	DOVER Foley Cath., 100% Silicone, Male, 10ml, 2-way, 16Fr, 43cm	UR	3
8887-605189	DOVER Foley Cath., 100% Silicone, Male, 10ml, 2-way, 18Fr, 43cm	UR	3
8887-605205	DOVER Foley Cath., 100% Silicone, Male, 10ml, 2-way, 20Fr, 43cm	UR	3
8887-605221	DOVER Foley Cath., 100% Silicone, Male, 10ml, 2-way, 22Fr, 43cm	UR	3
8887-605247	DOVER Foley Cath., 100% Silicone, Male, 10ml, 2-way, 24Fr, 43cm	UR	3
8887-605262	DOVER Foley Cath., 100% Silicone, Male, 10ml, 2-way, 26Fr, 43cm	UR	3
8887-607169	DOVER Foley Catheterisation Tray w/Sil. Elast. Cath., 16Fr	UR	7
8887-607185	DOVER Foley Catheterisation Tray w/Sil. Elast. Cath., 18Fr	UR	7
8887-615121	DOVER Foley Cath., 100% Silicone, Female, 5ml, 2-way, 12Fr, 23cm	UR	3
8887-615147	DOVER Foley Cath., 100% Silicone, Female, 5ml, 2-way, 14Fr, 23cm	UR	3
8887-615162	DOVER Foley Cath., 100% Silicone, Female, 5ml, 2-way, 16Fr, 23cm	UR	3
8887-615188	DOVER Foley Cath., 100% Silicone, Female, 5ml, 2-way, 18Fr, 23cm	UR	3
8887-630161	DOVER Foley Cath., 100% Silicone, Male, 30ml, 2-way, 16Fr, 43cm	UR	3
8887-630187	DOVER Foley Cath., 100% Silicone, Male, 30ml, 2-way, 18Fr, 43cm	UR	3
8887-630203	DOVER Foley Cath., 100% Silicone, Male, 30ml, 2-way, 20Fr, 43cm	UR	3
8887-630229	DOVER Foley Cath., 100% Silicone, Male, 30ml, 2-way, 22Fr, 43cm	UR	3
8887-630245	DOVER Foley Cath., 100% Silicone, Male, 30ml, 2-way, 24Fr, 43cm	UR	3
8887-630260	DOVER Foley Cath., 100% Silicone, Male, 30ml, 2-way, 26Fr, 43cm	UR	3
8887-640665	DOVER Foley Cath. Tray w/100% Sil. Cath., 16Fr	UR	15
8887-640681	DOVER Foley Cath. Tray w/100% Sil. Cath., 18Fr	UR	15
8887-640962	DOVER MEASURIT Urine Meter Foley Tray, Lge, 16Fr	UR	18
8887-640988	DOVER MEASURIT Urine Meter Foley Tray, Lge, 18Fr	UR	18
8887-641663	DOVER Foley Cath. Tray w/100% Sil. Cath., 16Fr	UR	15
8887-641689	DOVER Foley Cath. Tray w/100% Sil. Cath., 18Fr	UR	15
8887-641762	DOVER Foley Cath. Tray w/100% Sil. Cath., 16Fr	UR	15
8887-641788	DOVER Foley Cath. Tray w/100% Sil. Cath., 18Fr	UR	15
8887-648221	DOVER T.U.R.P. Balloon Cath., Sil./Elast. Latex, 2-way, 22Fr, 30-50ml	UR	5
8887-648247	DOVER T.U.R.P. Balloon Cath., Sil./Elast. Latex, 2-way, 24Fr, 30-50ml	UR	5
8887-649229	DOVER T.U.R.P. Balloon Cath., Sil./Elast. Latex, 2-way, 22Fr, 30-50ml	UR	5
8887-649245	DOVER T.U.R.P. Balloon Cath., Sil./Elast. Latex, 2-way, 24Fr, 30-50ml	UR	5
8887-665183	DOVER Foley Cath., 100% Silicone, Male, 30ml, 3-way, 18Fr, 43cm	UR	3

## Index by Code Number

Code	Product	Catalogue	Page
8887-665209	DOVER Foley Cath., 100% Silicone, Male, 30ml, 3-way, 20Fr, 43cm	UR	3
8887-665225	DOVER Foley Cath., 100% Silicone, Male, 30ml, 3-way, 22Fr, 43cm	UR	3
8887-665241	DOVER Foley Cath., 100% Silicone, 30ml, Male, 3-way, 24Fr, 43cm	UR	3
8887-670662	DOVER Foley Cath. Tray w/Sil. Elast. Coated Latex Cath., 16Fr	UR	15
8887-670688	DOVER Foley Cath. Tray w/Sil. Elast. Coated Latex Cath., 18Fr	UR	15
8887-671660	DOVER Foley Cath. Tray w/Sil. Elast. Coated Latex Cath., 16Fr	UR	15
8887-671686	DOVER Foley Cath. Tray w/Sil. Elast. Coated Latex Cath., 18Fr	UR	15
8887-671769	DOVER Foley Cath. Tray w/Sil. Elast. Coated Latex Cath., 16Fr	UR	15
8887-671785	DOVER Foley Cath. Tray w/Sil. Elast. Coated Latex Cath., 18Fr	UR	15
90051 T	CURITY Temperature Monitoring Sil. Cath., 5ml, 16Fr	UR	2
90052 T	CURITY Temperature Monitoring Sil. Cath., 5ml, 18Fr	UR	2
90053 T	CURITY Temperature Monitoring Sil. Cath., 3ml, 8Fr	UR	2
90054 T	CURITY Temperature Monitoring Sil. Cath., 3ml, 10Fr	UR	2
90055 T	CURITY Temperature Monitoring Sil. Cath., 5ml, 12Fr	UR	2
90056 T	CURITY Temperature Monitoring Sil. Cath., 5ml, 14Fr	UR	2
9718-02	CURITY Foley Cath., Sil. Coated Latex, 30ml, 3-way, 18Fr, 40cm	UR	5
9720-02	CURITY Foley Cath., Sil. Coated Latex, 30ml, 3-way, 20Fr, 40cm	UR	5
9722-02	CURITY Foley Cath., Sil. Coated Latex, 30ml, 3-way, 22Fr, 40cm	UR	5
9724-02	CURITY Foley Cath., Sil. Coated Latex, 30ml, 3-way, 24Fr, 40cm	UR	5
9812-02	CURITY Foley Cath., Sil. Coated Latex, Tiemann Type, 5ml, 2-way, 12Fr, 40cm	UR	5
9814-02	CURITY Foley Cath., Sil. Coated Latex, Tiemann Type, 5ml, 2-way, 14Fr, 40cm	UR	5
9816-02	CURITY Foley Cath., Sil. Coated Latex, Tiemann Type, 5ml, 2-way, 16Fr, 40cm	UR	5
9818-02	CURITY Foley Cath., Sil. Coated Latex, Tiemann Type, 5ml, 2-way, 18Fr, 40cm	UR	5
9820-02	CURITY Foley Cath., Sil. Coated Latex, Tiemann Type, 5ml, 2-way, 20Fr, 40cm	UR	5
9822-02	CURITY Foley Cath., Sil. Coated Latex, Tiemann Type, 5ml, 2-way, 22Fr, 40cm	UR	5
9824-02	CURITY Foley Cath., Sil. Coated Latex, Tiemann Type, 5ml, 2-way, 24Fr, 40cm	UR	5
NT0620	CURITY Intermittent Urethral Cath., PVC, 8Fr, 40cm	UR	6
NT0820	CURITY Intermittent Urethral Cath., PVC, 8Fr, 40cm	UR	6
NT1020	CURITY Intermittent Urethral Cath., PVC, 10Fr, 40cm	UR	6
NT1220	CURITY Intermittent Urethral Cath., PVC, 12Fr, 40cm	UR	6
NT1420	CURITY Intermittent Urethral Cath., PVC, 14Fr, 40cm	UR	6
NT1620	CURITY Intermittent Urethral Cath., PVC, 16Fr, 40cm	UR	6
NT1820	CURITY Intermittent Urethral Cath., PVC, 18Fr, 40cm	UR	6
NT2020	CURITY Intermittent Urethral Cath., PVC, 20Fr, 40cm	UR	6
UB3010	CURITY Economy Drain Bag w/Pull/Push Drain, 90cm Drain Tube	UR	12
UB3020	CURITY Economy Drain Bag w/Pull/Push Drain, 120cm Drain Tube	UR	12
UB3030	CURITY Economy Drain Bag w/Pull/Push Drain, 90cm Drain Tube	UR	12
UB3101	CURITY Economy Drain Bag w/Pull/Push Drain, 90cm Drain Tube	UR	12
UB3102	CURITY Economy Drain Bag w/Pull/Push Drain, 90cm Drain Tube	UR	12
UB3103	CURITY Economy Drain Bag w/Pull/Push Drain, 120cm Drain Tube	UR	12
UB3990	CURITY Economy Drain Bag w/o Drain, 90cm Drain Tube, Non-Sterile	UR	12
UB3994	CURITY Economy Drain Bag w/Twist Drain, 120cm Drain Tube, Non-Sterile	UR	12
UB3996	CURITY Economy Drain Bag w/Twist Drain, 120cm Drain Tube	UR	12
UB6207	CURITY Economy Drain Bag w/Pull/Push Drain, 120cm Drain Tube	UR	12

Index by  
Code Number

Notes





---

**tyco**

*Healthcare*



**KENDALL**

**Regional Export,**

154 Fareham Road, Gosport, Hampshire, PO13 0AS, UK.

**Tel: +44 (0)1329 224451 Fax: +44 (0)1329 224086**