



Jupiterimages/Jupiterimages/Thinkstock

Urology Products
1-800-962-9888



COVIDIEN

Urology Products

Table of Contents

Dover™ Silver Products

Foley Catheters	1
Urine Meter Foley Trays	2
Foley Trays	3
Insertion Trays	3

Dover™ 100% Silicone Products

Foley Catheters	4
Urine Meter Foley Trays	5-6
Insertion Trays	6
Foley Trays	7

Dover™ Latex Products

Foley Catheters	8-9
Urine Meter Foley Trays	10
Foley Trays	11-12
Insertion Trays	12

Catheter Accessories

Universal Trays	13
Miscellaneous Accessories	13

Dover™ Collection Products

Urine Meters	14
Drainage Bags	15
Add-A-Foley Trays	15

Dover™ Urethral Products

Latex Free Catheters	16
Latex Catheters	17
Urethral Trays	18

Ancillary Products

Male External Catheters	19
Leg Bags & Accessories	19
Irrigation Trays	20
Midstream Catch Kits	20
Specimen Containers	21
Patient Care Products (Latex Free)	21

Index

Index	22-25
-------------	-------

Dover™ Silver Foley Catheters



- 100% silicone with silver hydrogel coating
- 16" length
- Rounded tip, 2 opposing eyes
- 100% silicone Foley catheter eliminates latex allergy concerns
- Outer and inner catheter and balloon surfaces contain a patented silver and hydrogel coating which provides a consistent release of ionic silver particles
- Hydrogel coating eases insertion to prevent urethra trauma
- Larger internal lumens offer improved flow rate for better drainage

Code	Size (Fr)	Ship Case
5 cc, 2-Way		
605122IC	12 Fr	10
605148IC	14 Fr	10
605163IC	16 Fr	10
605189IC	18Fr	10
605205IC	20 Fr	10
605221IC	22 Fr	10
605247IC	24 Fr	10
605262IC	26Fr	10

30 cc, 2-Way

630187IC	18 Fr	10
630203IC	20 Fr	10
630229IC	22 Fr	10
630245IC	24 Fr	10
630260IC	26 Fr	10

5 cc, 2-Way, 400 Series Temperature Sensing

50516ICT	16 Fr	12
50518ICT	18 Fr	12

Code	Size (Fr)	Ship Case
30 cc, 3-Way		
665183IC	18 Fr	10
665209IC	20 Fr	10
665225IC	22 Fr	10
665241IC	24 Fr	10

Dover™ Silver Urine Meter Trays



P4P16XSD

- Dover™ urine meter with 100% silicone catheter with silver hydrogel coating
- 400 mL urine meter capacity
- Extra long 60" tubing
- Standard Foley prep tray (100% latex-free)
Sterile tray contains: CSR wrap, fenestrated drape, underpad, blue nitrile exam gloves, forceps, prepping balls, 10 cc prefilled sterile water inflation syringe, unique bellows lubricant syringe, povidone-iodine solution, specimen container and label



Code	Size (Fr)	Securement	Top Vent	Sample Port Luer Lock	Temperature Sensing 400	Tube Management	Sheet Clip	Meter Port		Enclosed Drain Spout	Self-leveling Bed Hook	Ship Case
								Luer Lock	Spring			
P4P14XSD	14 Fr	•	•	•		•	•	•		•	•	10
P4P16XSD	16 Fr	•	•	•		•	•	•		•	•	10
P4P18XSD	18 Fr	•	•	•		•	•	•		•	•	10
P4P14XS	14 Fr		•	•		•	•	•		•	•	10
P4P16XS	16 Fr		•	•		•	•	•		•	•	10
P4P18XS	18 Fr		•	•		•	•	•		•	•	10
P4P16XTSD	16 Fr	•	•	•	•	•	•	•		•	•	10
P4P18XTSD	18 Fr	•	•	•	•	•	•	•		•	•	10
P4P16XTS	16 Fr		•	•	•	•	•	•		•	•	10
P4P18XTS	18 Fr		•	•	•	•	•	•		•	•	10
7004ICLL	14 Fr			•			•		•	•	•	10
7006ICLL	16 Fr			•			•		•	•	•	10
7008ICLL	18 Fr			•			•		•	•	•	10

Dover™ Silver Foley Trays



PP16XSD

- Dover™ Foley trays with 100% silicone catheter with silver hydrogel coating
- 2000 mL capacity
- 48" tubing
- Standard Foley prep tray (100% latex-free)
Sterile tray contains: CSR wrap, fenestrated drape, underpad, blue nitrile exam gloves, forceps, prepping balls, 10 cc prefilled sterile water inflation syringe, unique bellows lubricant syringe, povidone-iodine solution, specimen container and label



Code	Size (Fr)	Securement	Top Vent	Sample Port Luer Lock	Tube Management	Sheet Clip	Enclosed Drain Spout	Anti-Reflux Device	Self-leveling Bed Hook	Ship Case
PP14XSD	14 Fr	•	•	•	•	•	•	•	•	10
PP16XSD	16 Fr	•	•	•	•	•	•	•	•	10
PP18XSD	18 Fr	•	•	•	•	•	•	•	•	10
PP14XS	14 Fr		•	•	•	•	•	•	•	10
PP16XS	16 Fr		•	•	•	•	•	•	•	10
PP18XS	18 Fr		•	•	•	•	•	•	•	10
6144ICLL	14 Fr			•		•	•	•	•	10
6146ICLL	16 Fr			•		•	•	•	•	10
6148ICLL	18 Fr			•		•	•	•	•	10

Dover™ Silver Insertion Trays

Dover™ 100% Silicone Foley Catheter with Silver Hydrogel Coating

Sterile tray contains: CSR wrap, underpad, fenestrated drape, PVP solution, prepping balls, forceps, lubricating jelly, 10cc prefilled water syringe, specimen container, graduated basin, blue nitrile exam gloves

Code	Size (Fr)	Ship Case
2116IC	16 Fr	20
2118IC	18 Fr	20



Dover™ Silicone Foley Catheters



8887605163

- 100% silicone Foley catheter eliminates latex allergy concern
- 16" length
- Biocompatible, inorganic, non-porus, non-toxic, non-adhering material reduces incidence of encrustation, urethritis, and stricture vs. latex catheters
- Color-coded inflation valve accepts Luer lock and Luer slip syringes
- Packaged in an outer peel-back pouch and an easy-to-open inner sheath

2-Way

Code	Size (Fr)	Ship Case
3 cc, 2-way (Pediatric)		
8887603085	8 Fr	10
8887603101	10 Fr	10
6108	8 Fr (Luer Lock)	12
6010	10 Fr (Luer Lock)	12

5 cc, 2-way

8887605122	12 Fr	10
8887605148	14 Fr	10
8887605163	16 Fr	10
8887605189	18 Fr	10
8887605205	20 Fr	10
8887605221	22 Fr	10
8887605247	24 Fr	10
8887605262	26 Fr	10
8887605282	28 Fr	10
8887605303	30 Fr	10

30 cc, 2-way

8887630142	14 Fr	10
8887630163	16 Fr	10
8887630187	18 Fr	10
8887630203	20 Fr	10
8887630229	22 Fr	10
8887630245	24 Fr	10
8887630260	26 Fr	10
8887630284	28 Fr	10
8887630305	30 Fr	10

3-Way

Code	Size (Fr)	Ship Case
5 cc, 3-way		
8887664161	16 Fr	10
8887664182	18 Fr	10
8887664203	20 Fr	10
8887664224	22 Fr	10
8887664245	24 Fr	10
8887664266	26 Fr	10
8887664287	28 Fr	10

30 cc, 3-way

8887665161	16 Fr	10
8887665183	18 Fr	10
8887665209	20 Fr	10
8887665225	22 Fr	10
8887665241	24 Fr	10
8887665262	26 Fr	10
8887665283	28 Fr	10

Specialty

Code	Size (Fr)	Ship Case
5 cc, 2-way, 400 Series Temperature Sensing		
50514T	14 Fr	12
50516T	16 Fr	12
50518T	18 Fr	12

5 cc, 2-way, Coudé tip

20512C	12 Fr	10
20514C	14 Fr	10
20516C	16 Fr	10
20518C	18 Fr	10
20520C	20 Fr	10
20522C	22 Fr	10
20524C	24 Fr	10

30 cc, 2-way, Coudé tip

23018C	18 Fr	10
23020C	20 Fr	10
23022C	22 Fr	10
23024C	24 Fr	10
23026C	26 Fr	10

5 cc, 2-way, Council tip

40516L	16 Fr	10
40518L	18 Fr	10
40520L	20 Fr	10
40522L	22 Fr	10
40524L	24 Fr	10

Dover™ Silicone Urine Meter Trays



P4P16SD

- 100% silicone Foley catheter
- Extra long 60" tubing
- 400 mL or 200 mL urine meter capacity
- Standard Foley prep tray (100% latex-free)
Sterile tray contains: CSR wrap, fenestrated drape, underpad, blue nitrile exam gloves, forceps, prepping balls, 10 cc prefilled sterile water inflation syringe, unique bellows lubricant syringe, povidone-iodine solution, specimen container and label



Code	Size (Fr)	Securement	Top Vent	Sample Port		Temperature Sensing 400	Tube Management	Sheet Clip	Luer Lock Meter Port	Enclosed Drain Spout	Self-leveling Bed Hook	Ship Case
				Luer Lock	Needleless							
400 mL												
P4P14SD	14 Fr	•	•	•			•	•	•	•	•	10
P4P16SD	16 Fr	•	•	•			•	•	•	•	•	10
P4P18SD	18 Fr	•	•	•			•	•	•	•	•	10
P4P14S	14 Fr		•	•			•	•	•	•	•	10
P4P16S	16 Fr		•	•			•	•	•	•	•	10
P4P18S	18 Fr		•	•			•	•	•	•	•	10
P4A14S	14 Fr		•		•		•	•	•	•	•	10
P4A16S	16 Fr		•		•		•	•	•	•	•	10
P4A18S	18 Fr		•		•		•	•	•	•	•	10
P4P14TSD	14 Fr	•	•	•		•	•	•	•	•	•	10
P4P16TSD	16 Fr	•	•	•		•	•	•	•	•	•	10
P4P18TSD	18 Fr	•	•	•		•	•	•	•	•	•	10
P4P14TS	14 Fr		•	•		•	•	•	•	•	•	10
P4P16TS	16 Fr		•	•		•	•	•	•	•	•	10
P4P18TS	18 Fr		•	•		•	•	•	•	•	•	10
200 mL												
P2A16S	16 Fr		•		•		•	•	•	•	•	10
P2A18S	18 Fr		•		•		•	•	•	•	•	10

Dover™ Silicone Urine Meter Trays (continued)



Code	Size (Fr)	Sample Port Luer Lock	Sheet Clip	Spring Meter Port	Drain Enclosed Spout	Self-leveling Bed Hook	Ship Case
------	-----------	-----------------------	------------	-------------------	----------------------	------------------------	-----------

400 mL

2004LL	14 Fr	•	•	•	•	•	10
7006LL	16 Fr	•	•	•	•	•	10
7008LL	18 Fr	•	•	•	•	•	10

Dover™ Silicone Insertion Trays

Dover™ 100% Silicone Foley Catheter

Tray contains: CSR wrap, underpad, fenestrated drape, PVP solution, prepping balls, forceps, lubricating jelly, 10cc prefilled water syringe, specimen container, graduated basin, blue nitrile exam gloves



8887600776

Code	Size (Fr)	Ship Case
8887600776	16 Fr	20
8887600784	18 Fr	20

Dover™ Silicone Foley Trays



PP16SD

- 100% silicone Foley catheters
- 2000 mL capacity
- 48" tubing
- Standard Foley prep tray (100% latex-free)
Sterile tray contains: CSR wrap, fenestrated drape, underpad, blue nitrile exam gloves, forceps, prepping balls, 10 cc prefilled sterile water inflation syringe, unique bellows lubricant syringe, povidone-iodine solution, specimen container and label



Code	Size (Fr)	Securement	Top Vent	Sample Port		Tube Management	Sheet Clip	Drain Enclosed Spout	Anti-Reflux Device	Self-leveling Bed Hook	Ship Case
				Luer Lock	Needleless						
PP14SD	14 Fr	•	•	•		•	•	•	•	•	10
PP16SD	16 Fr	•	•	•		•	•	•	•	•	10
PP18SD	18 Fr	•	•	•		•	•	•	•	•	10
PP14S	14 Fr		•	•		•	•	•	•	•	10
PP16S	16 Fr		•	•		•	•	•	•	•	10
PP18S	18 Fr		•	•		•	•	•	•	•	10
PA14S	14 Fr				•	•	•	•	•	•	10
PA16S	16 Fr				•	•	•	•	•	•	10
PA18S	18 Fr				•	•	•	•	•	•	10
6144LL	14 Fr			•			•	•	•	•	10
6146LL	16 Fr			•			•	•	•	•	10
6148LL	18 Fr			•			•	•	•	•	10

Dover™ Latex Foley Catheters



- 16" length
- Coatings ease insertion to help prevent urethra trauma

Hydrogel Coating

Code	Size (Fr)	Coating	Ship Case
5 cc, 2-way			
1612	12 Fr	Hydrogel	12
1614	14 Fr	Hydrogel	12
1616	16 Fr	Hydrogel	12
1619	18 Fr	Hydrogel	12
1621	20 Fr	Hydrogel	12
1622	22 Fr	Hydrogel	12
1624	24 Fr	Hydrogel	12
1626	26 Fr	Hydrogel	12
1628	28 Fr	Hydrogel	12
1631	30 Fr	Hydrogel	12

30 cc, 2-way

1416	16 Fr	Hydrogel	12
1419	18 Fr	Hydrogel	12
1420	20 Fr	Hydrogel	12
1422	22 Fr	Hydrogel	12
1424	24 Fr	Hydrogel	12
1426	26 Fr	Hydrogel	12
1427	28 Fr	Hydrogel	12
1431	30 Fr	Hydrogel	12

5 cc, 2-way Kits

Includes: latex Foley catheter with lubricating jelly and 10 cc syringe

1814	14 Fr	Hydrogel	12
1817	16 Fr	Hydrogel	12
1818	18 Fr	Hydrogel	12

5 cc, 3-way

2718	18 Fr	Hydrogel	12
2720	20 Fr	Hydrogel	12
2722	22 Fr	Hydrogel	12
2723	24 Fr	Hydrogel	12

30 cc, 3-way

2819	18 Fr	Hydrogel	12
2821	20 Fr	Hydrogel	12
2822	22 Fr	Hydrogel	12
2824	24 Fr	Hydrogel	12
2826	26 Fr	Hydrogel	12

Silicone Elastomer Coating

Code	Size (Fr)	Coating	Ship Case
3 cc, 2-way (Pediatric)			
8887654088	8 Fr	Silicone Elastomer	10
8887654101	10 Fr	Silicone Elastomer	10

5 cc, 2-way

402712	12 Fr	Silicone Elastomer	10
402714	14 Fr	Silicone Elastomer	10
402716	16 Fr	Silicone Elastomer	10
402718	18 Fr	Silicone Elastomer	10
402720	20 Fr	Silicone Elastomer	10
402722	22 Fr	Silicone Elastomer	10
402724	24 Fr	Silicone Elastomer	10
402726	26 Fr	Silicone Elastomer	10
402728	28 Fr	Silicone Elastomer	10
402730	30 Fr	Silicone Elastomer	10

30 cc, 2-way

403712	12 Fr	Silicone Elastomer	10
403714	14 Fr	Silicone Elastomer	10
403716	16 Fr	Silicone Elastomer	10
403718	18 Fr	Silicone Elastomer	10
403720	20 Fr	Silicone Elastomer	10
403722	22 Fr	Silicone Elastomer	10
403724	24 Fr	Silicone Elastomer	10
403726	26 Fr	Silicone Elastomer	10
403728	28 Fr	Silicone Elastomer	10
403730	30 Fr	Silicone Elastomer	10

30 cc, 3-way

8887689167	16 Fr	Silicone Elastomer	10
8887689183	18 Fr	Silicone Elastomer	10
8887689209	20 Fr	Silicone Elastomer	10
8887689225	22 Fr	Silicone Elastomer	10
8887689241	24 Fr	Silicone Elastomer	10
8887689266	26 Fr	Silicone Elastomer	10

5 cc, 3-way

8887688169	16 Fr	Silicone Elastomer	10
8887688185	18 Fr	Silicone Elastomer	10
8887688201	20 Fr	Silicone Elastomer	10
8887688227	22 Fr	Silicone Elastomer	10
8887688243	24 Fr	Silicone Elastomer	10
8887688268	26 Fr	Silicone Elastomer	10

Dover™ Latex Foley Catheters (continued)

Silicone Oil Coating

Code	Size (Fr)	Coating	Ship Case
5 cc, 2-way			
8887623125	12 Fr	Silicone oil	10
3558	14 Fr	Silicone oil	10
3560	16 Fr	Silicone oil	10
3563	18 Fr	Silicone oil	10
3565	20 Fr	Silicone oil	10
3567	22 Fr	Silicone oil	10
3570	24 Fr	Silicone oil	10
3573	26 Fr	Silicone oil	10
3575	28 Fr	Silicone oil	10
3577	30 Fr	Silicone oil	10

30 cc, 2-way,

8887624149	14 Fr	Silicone oil	10
3601	16 Fr	Silicone oil	10
3607	18 Fr	Silicone oil	10
3611	20 Fr	Silicone oil	10
3614	22 Fr	Silicone oil	10
3618	24 Fr	Silicone oil	10
3619	26 Fr	Silicone oil	10
3623	28 Fr	Silicone oil	10
3631	30 Fr	Silicone oil	10

Specialty Catheters

Code	Size (Fr)	Coating	Ship Case
5 cc, 2-way, Coudé tip			
1612C	12 Fr	Hydrogel	12
1614C	14 Fr	Hydrogel	12
1616C	16 Fr	Hydrogel	12
1618C	18 Fr	Hydrogel	12
1620C	20 Fr	Hydrogel	12
1622C	22 Fr	Hydrogel	12
1624C	24 Fr	Hydrogel	12

5 cc, 2-way, Red Latex, Radiopaque, Coudé tip

1512C	12 Fr	Hydrogel	12
1514C	14 Fr	Hydrogel	12
1516C	16 Fr	Hydrogel	12
1518C	18 Fr	Hydrogel	12
1520C	20 Fr	Hydrogel	12
1522C	22 Fr	Hydrogel	12
1524C	24 Fr	Hydrogel	12

30 cc, 2-way, Coudé tip

1416C	16 Fr	Hydrogel	12
1418C	18 Fr	Hydrogel	12
1420C	20 Fr	Hydrogel	12
1422C	22 Fr	Hydrogel	12
1424C	24 Fr	Hydrogel	12

T.U.R.P 30 – 50 cc, up to 90 cc, 2-way, 6 staggered eyes

8887648221	22 Fr	Silicone Elastomer	10
8887648247	24 Fr	Silicone Elastomer	10

T.U.R.P 30 – 50 cc, up to 90 cc, 3-way, 4 opposing eyes

8887649229	22 Fr	Silicone Elastomer	10
8887649245	24 Fr	Silicone Elastomer	10



1616C



1516C



888764247

Dover™ Latex Urine Meter Foley Trays



P4P16ULD

- Latex Foley catheter with hydrogel coating
- 400 mL or 200 mL urine meter capacity
- Extra long 60" tubing
- Standard Foley prep tray (100% latex-free)
Sterile tray contains: CSR wrap, fenestrated drape, underpad, blue nitrile exam gloves, forceps, prepping balls, 10 cc prefilled sterile water inflation syringe, unique bellows lubricant syringe, povidone-iodine solution, specimen container and label



Code	Size (Fr)	Securement	Top Vent	Sample Port		Tube Management	Sheet Clip	Meter Port		Enclosed Drain Spout	Self-leveling Bed Hook	Ship Case
				Luer Lock	Needless			Luer Lock	Spring			

400 mL

P4P14ULD	14 Fr	•	•	•		•	•	•		•	•	10
P4P16ULD	16 Fr	•	•	•		•	•	•		•	•	10
P4P18ULD	18 Fr	•	•	•		•	•	•		•	•	10
P4P14UL	14 Fr		•	•		•	•	•		•	•	10
P4P16UL	16 Fr		•	•		•	•	•		•	•	10
P4P18UL	18 Fr		•	•		•	•	•		•	•	10
P4A14UL	14 Fr		•		•	•	•	•		•	•	10
P4A16UL	16 Fr		•		•	•	•	•		•	•	10
P4A18UL	18 Fr		•		•	•	•	•		•	•	10

400 mL, Coudé tip

P4P16UCL	16 Fr		•	•		•	•	•		•	•	10
----------	-------	--	---	---	--	---	---	---	--	---	---	----

400 mL

7014LL	14 Fr			•		•		•	•	•	•	10
7016LL	16 Fr			•		•		•	•	•	•	10
7018LL	18 Fr			•		•		•	•	•	•	10

200 mL

P2A16UL	16 Fr		•		•	•	•			•	•	10
P2A18UL	18 Fr		•		•	•	•			•	•	10

Dover™ Latex Foley Trays



PP16ULD

- Latex Foley catheter with hydrogel coating
- 2000 mL capacity
- 48" tubing
- Standard Foley prep tray (100% latex-free)
Sterile tray contains: CSR wrap, fenestrated drape, underpad, blue nitrile exam gloves, forceps, prepping balls, 10 cc prefilled sterile water inflation syringe, unique bellows lubricant syringe, povidone-iodine solution, specimen container and label



Code	Size (Fr)	Securement	Top Vent	Sample Port		Tube Management	Sheet Clip	Enclosed Drain Spout	Anti-Reflux Device	Self-leveling Bed Hook	Ship Case
				Luer Lock	Needleless						
PP14ULD	14 Fr	•	•	•		•	•	•	•	•	10
PP16ULD	16 Fr	•	•	•		•	•	•	•	•	10
PP18ULD	18 Fr	•	•	•		•	•	•	•	•	10
PP14UL	14 Fr		•	•		•	•	•	•	•	10
PP16UL	16 Fr		•	•		•	•	•	•	•	10
PP18UL	18 Fr		•	•		•	•	•	•	•	10
PA14UL	14 Fr				•	•	•	•	•	•	10
PA16UL	16 Fr				•	•	•	•	•	•	10
PA18UL	18 Fr				•	•	•	•	•	•	10

Coudé tip

PP16UCL	16 Fr		•	•		•	•	•	•	•	10
---------	-------	--	---	---	--	---	---	---	---	---	----

Dover™ Latex Foley Trays (continued)



Code	Size (Fr)	Sample Port		Sheet Clip	Drain		Anti-Reflux Device	Self-leveling Bed Hook	Ship Case
		Luer Lock	Needle		Enclosed Spout	Tube			
6164LL	14 Fr	•		•	•		•	•	10
6165LL	16 Fr	•		•	•		•	•	10
6166LL	18 Fr	•		•	•		•	•	10

Silicone coating, 40" tubing

3716	16 Fr		•	•		•	•	•	10
3718	18 Fr		•	•		•	•	•	10

Dover™ Latex Insertion Trays

Dover™ Latex Foley catheter

Sterile tray contains: CSR wrap, underpad, fenestrated drape, PVP solution, prepping balls, forceps, lubricating jelly, 10cc prefilled water syringe, specimen container, graduated basin, blue nitrile exam gloves



Code	Size (Fr)	Ship Case
2101	16 Fr	20
3044	18 Fr	20

Catheter Accessories

Universal Trays

All trays feature (unless indicated otherwise): CSR wrap, underpad, blue nitrile exam gloves, PVP solution 10%, prepping balls, forceps, lubricating jelly (5 g).



5029

Code	Antiseptic	Syringe	Exam Gloves	Cotton balls	Forceps	Lubricating Jelly	Specimen Container	Underpad/ Drape	Lid	Ship Case
5029	PVP Solution	10 cc	•	•	•	•	•	•	•	20
7100	BZK Swabs	10 cc	•			•	•	•	•	20
5027	PVP Solution	30 cc	•	•	•	•	•	•	•	20
76000	PVP Swabs	30 cc	•			•		•	•	20
76010	PVP Swabs	10 cc	•			•		•	•	20
76012	PVP Swabs	10 cc	•			•	•	•	•	20
76020	BZK Swabs	30 cc	•			•		•	•	20
76030	BZK Swabs	10 cc	•			•		•	•	20

Miscellaneous Accessories



54456

Code	Description	Ship Case
54456	Dover™ Catheter Securement	50
1600	Catheter Plug	50
95516	Sterile Catheter Care Kit	50
8887600149	Foley Leg Strap	10
1106	Lubricating Jelly Bellows Syringe (10 g)	280
8887636002	Lubricated Jelly Pre-filled Syringe (5 cc)	50

Collection



Urine Meter

Designed for precise measurement when output is critically important. Capacity ranges from 150 mL to 400 mL.

Next Generation:

- Lowered & centered entry point
- Horizontal face plate with self leveling chamber
- Top vent promotes optimal flow within the collection system
- Luer Lock sample port and meter port helps prevent needle stick injury
- Tube management & sheet clip improve flow dynamics

Urine Meters



Code	Meter Capacity	Top Vent	Sample Port Luer Lock	Tube Management	Sheet Clip	Meter Port		Drain Enclosed Spout	Self-leveling Bed Hook	Ship Case
						Luer Lock	Spring			

Extra-long 60" tubing

P4P	400 mL	•	•	•	•	•		•	•	10
P4A	400 mL	•	•	•	•	•		•	•	10
P2A	200 mL	•	•	•	•	•		•	•	10
P2K	200 mL		•			•		•	•	10
7000LL	400 mL		•		•		•	•	•	10

Add-A-Foley 400 mL Urine Meter Trays

Standard Foley prep tray (100% latex-free) Sterile tray contains: CSR wrap, fenestrated drape, underpad, blue nitrile exam gloves, forceps, prepping balls, unique bellows lubricant syringe, povidone-iodine solution, specimen container and label

Code	Syringe	Ship Case
P4PA	10 cc	10
P4AA	10 cc	10
P2AA	10 cc	10

Collection



Drainage Bag

Typical capacity is 2000 mL, but range from 1000 mL to 4000 mL.

Next Generation:

- Updated enclosed drain spout on a new 45° pedestal
- Pivoting hanger hooks
- Top vent promotes optimal flow within the collection system
- Luer Lock sample port helps in preventing needle stick injury
- Tube management & sheet clip improve flow dynamics

Drainage Bags



Code	Size	Top Vent	Sample Port			Tube Management	Sheet Clip	Drain		Anti-reflux Device	Hanger			Ship Case
			Luer Lock	Needle-less	Needle			Enclosed Spout	Drain Tube		Bed Hook	String	Hook & Loop	
48" tubing														
PP	2000 mL	•	•			•	•	•		•	•			20
PA	2000 mL			•			•	•		•	•			20
PAB *	2000 mL			•			•	•		•	•			20
6251	4000 mL				•		•		2	2		•		15
6261	4000 mL				•		•		•	•		•		20
6308LL	2000 mL		•				•	•		•	•			20
40" tubing														
3512N *	2000 mL			•					•	•	•			20
3512V *	2000 mL				•				•	•			•	20
3512 *	2000 mL				•				•	•	•			20
3513CF	2000 mL								•	•			•	20
3502 *	2000 mL				•				•		•			20

*Indicates poly bag packaging – interior of drainage bag and fluid pathway are guaranteed sterile

Add-A-Foley 2000 mL Drainage Bag Foley Trays

Standard Foley prep tray (100% latex-free) Sterile tray contains: CSR wrap, fenestrated drape, underpad, blue nitrile exam gloves, forceps, prepping balls, unique bellows lubricant syringe, povidone-iodine solution, specimen container and label

Code	Syringe	Ship Case
PPA	10 cc	10
PAA	10 cc	10
3515	10 cc	10
3531	30 cc	10
3532	30 cc	10

Dover™ Urethral Catheters

100% Latex Free Urethral Catheters (PVC)



8888492041

Code	Size (Fr)	Ship Case
------	-----------	-----------

Rob-Nel catheters, 14" length, smooth rounded tip

8888492017	8 Fr	100
8888492025	10 Fr	100
8888492033	12 Fr	100
8888492041	14 Fr	100
8888492058	16 Fr	100
8888492066	18 Fr	100
8888492120	20 Fr	100
8888492122	22 Fr	100
8888492124	24 Fr	100
8888492126	26 Fr	100
8888492128	28 Fr	100
8888492130	30 Fr	100



400614

Robinson catheters, 14" length, smooth rounded tip, clear vinyl

400608	8 Fr	100
400610	10 Fr	100
400612	12 Fr	100
400614	14 Fr	100
400616	16 Fr	100
400618	18 Fr	100



6006

Pediatric urethral catheters, 16" length, hydrogel coating and Luer lock adapter

6006	6 Fr (Ped.)	12
------	-------------	----



8887660903

Female catheter, 6" length

8887660903	14 Fr	50
------------	-------	----

Dover™ Urethral Catheters (continued)

Latex Urethral Catheters

Code	Size (Fr)	Ship Case
Red Rubber Coudé, 14" length, hydrophilic coating		
8403	14 Fr	12
8404	16 Fr	12



Red Rubber, 14" length, hydrophilic coating, smooth rounded tip		
8408	8 Fr	12
8410	10 Fr	12
8412	12 Fr	12
8414	14 Fr	12
8416	16 Fr	12
8418	18 Fr	12
8420	20 Fr	12
8422	22 Fr	12
8424	24 Fr	12



Red Rubber Robinson, 14" length, radiopaque, smooth rounded tip		
8887660085	8 Fr	100
8887660101	10 Fr	100
8887660127	12 Fr	100
8887660143	14 Fr	100
8887660168	16 Fr	100
8887660184	18 Fr	100
8887660218	20 Fr	100
8887660221	22 Fr	100
8887660242	24 Fr	100
8887660263	26 Fr	100
8887660284	28 Fr	100
8887660305	30 Fr	100



Dover™ Urethral Trays

Open Trays

All trays feature: CSR wrap, underpad, blue nitrile exam gloves, lubricating jelly, specimen container



3141

Code	Size (Fr)	Material	Antiseptic	Cotton Balls	Ship Case
------	-----------	----------	------------	--------------	-----------

Hydrogel coated catheters

3217	14 Fr	Red Rubber Latex	PVP Solution	•	20
------	-------	------------------	--------------	---	----

Urethral Catheter

3141	14 Fr	Vinyl	PVP Solution	•	20
3215	16 Fr	Red Rubber Latex	PVP Solution	•	20
406622	14 Fr	Vinyl	BZK Swabs		20
8887600057	14 Fr	Red Rubber Latex	PVP Solution	•	20
75000	14 Fr	Vinyl	BZK Swabs		20
75005	14 Fr	Red Rubber Latex	BZK Swabs		20
75030	14 Fr	Vinyl	PVP Swabs		20
75035	14 Fr	Red Rubber Latex	PVP Swabs		20

Add-A-Cath

Code	Antiseptic	Exam Gloves	Lubricating Jelly	Specimen Container	Underpad/ Drape	CSR Wrap	Ship Case
------	------------	-------------	-------------------	--------------------	-----------------	----------	-----------

Open Urethral Tray

3305	PVP Solution	•	•	•	•	•	20
3307	PVP Solution	•	•	•	•	•	20
75010	PVP Swab Sticks	•	•	•	•		20
75020	BZK Swab Sticks	•	•	•	•		20

Closed Trays

Trays include (unless indicated otherwise): CSR wrap, underpad, blue nitrile exam gloves, antiseptic, lubricating jelly, specimen container, 1500 mL collection bag



3411

Code	Size (Fr)	Material	Antiseptic	Ship Case
------	-----------	----------	------------	-----------

Hydrogel Coated Catheters

3410	14 Fr	Red Rubber Latex	PVP Solution	20
3408	16 Fr	Red Rubber Latex	PVP Solution	20

Clear Vinyl Catheters

3450	14 Fr	Vinyl	PVP Solution	20
2480	14 Fr	Vinyl	BZK Swabs	20

Female Catheterization Kit with 15 mL Specimen Tube and Flip Cap

3411	8 Fr	Vinyl	PVP Swabs	25
------	------	-------	-----------	----

Code	Size (Fr)	Material	Ship Case
------	-----------	----------	-----------

Urethral Catheter with 1000 mL bag only

3170	14 Fr	Vinyl	50
------	-------	-------	----

Dover™ Ancillary Products

Male External Catheters (Latex)



Code	Size	Attachment	Ship Case
------	------	------------	-----------

Uri-Drain Male External Catheter

Latex rolled sheath, molded one piece, individually packaged

8884732300	1.3"	Foam Strap	144
8884734600	1.1"	Foam Strap	144
8884732500	1.3"	Two Sided Adhesive	144
8884732000	1.1"	Two Sided Adhesive	144

Texas Catheter

Soft, lightweight latex with hard plastic tip

8884730300	1.3"	Foam Strap	144
8884731300	1.3"	Foam Strap, 12/box	144

Code	Size	Description	Ship Case
------	------	-------------	-----------

Adhesive Strap

8884730600	5.25" x 1"	Foam Strap	500
------------	------------	------------	-----

Leg Bags



Code	Small 9 oz.	Medium 17 oz. / 500 mL	Large 25 oz. / 730 mL	Latex	Extension Tubing	Reusable	Ship Case
3433		•			•		50
145514		•					50
145516			•				50
8887601121		•					20
8887601139			•				20
8884733100	•			•		•	50
8884733200		•		•		•	50
8884733900			•	•		•	50

Code	Description	Latex	Ship Case
------	-------------	-------	-----------

Leg Bag Accessory

8884731900	18" Urinary Extension Tube	•	24
------------	----------------------------	---	----

Dover™ Ancillary Products

Irrigation Trays



Code	1200 cc Tray	500 cc Container	Antiseptic	CSR Wrap	Sterile Saline	Lid	Ship Case
------	--------------	------------------	------------	----------	----------------	-----	-----------

60 cc Bulb Syringe Tray

67800	•	•				•	20
67803	•	•	•	•			20

60 cc Piston Syringe Tray

68800	•	•				•	20
68803	•	•	•	•			20
68820	•				•	•	20
3685	•	•	•	•			20



Syringes Only, 60 cc

Code	Description	Ship Case
67000	Sterile Bulb Syringe	50
68000	Sterile Piston Syringe	50

Midstream Catch Kits



Code	Capacity	Description	Patient Label	BZK Wipes	Castile Soap	Pneumatic Capable	Funnel	Screw Cap	Ship Case
22000	4.5 oz	Sterile	•	•			•		24
22015	4.5 oz	Sterile	•		•		•		24
24000	4.5 oz	Wide mouth container	•	•				•	50
24015	4.5 oz	Wide mouth container	•		•			•	50
24035	7 oz	Wide mouth container	•		•			•	50
25000	4.5 oz	Protective cover	•	•				•	100
25015	4.5 oz	Protective cover	•		•			•	100
5205	4.5 oz	Specimen container	•		•				48
5204	4.5 oz	Specimen container	•	•					48
2090SA	4 oz	Collection kit	•		•	•			100
4090SA	4 oz	Collection kit			•	•			100
2040SA	2 oz	Collection kit	•		•	•			100
2110SA		Calculi Strainer only							100

Dover™ Ancillary Products

Specimen Containers



2500SA

Code	Capacity	Lid	Ship Case
------	----------	-----	-----------

Measuring System

4014	26 oz		100
2400SA	700 cc	•	60
2500SA	800 cc		100



145500

Code	Capacity	Sterile	Ship Case
------	----------	---------	-----------

Pediatric Urine Collector

145500	100 mL		500
145501	100 mL	•	300

Patient Care Products (Latex Free)



145540

Code	Size (Fr)	Description	Tube length	Sterile	Pre-lubricated	Castile Soap	Underpad	Ship Case
145524	24 Fr	Flatus Bag	19"	•	•			50
145617		Vaginal Irrigation (Douche) Bag	60"	•				50
145540	24 Fr	Enema Bag	60"	•		•	•	50
155730	18 Fr	Rectal Tube	20"					50
155731	24 Fr	Rectal Tube	20"	•	•			50

Index

Code	HCPCS*	Latex	Page	Code	HCPCS	Latex	Page
1106	*	No	13	20512C	A4340	No	4
1416	A4338	Yes	8	20514C	A4340	No	4
1416C	A4340	Yes	9	20516C	A4340	No	4
1418C	A4340	Yes	9	20518C	A4340	No	4
1419	A4338	Yes	8	20520C	A4340	No	4
1420	A4338	Yes	8	20522C	A4340	No	4
1420C	A4340	Yes	9	20524C	A4340	No	4
1422	A4338	Yes	8	2090SA	A9270	No	20
1422C	A4340	Yes	9	2101	A4312	Yes	12
1424	A4338	Yes	8	2110SA	A9270	No	20
1424C	A4340	Yes	9	2116IC	A4312	No	3
1426	A4338	Yes	8	2118IC	A4312	No	3
1427	A4338	Yes	8	22000	A9270	No	20
1431	A4338	Yes	8	22015	A9270	No	20
145500	A9270	No	21	23018C	A4340	No	4
145501	A9270	No	21	23020C	A4340	No	4
145514	A4358	No	19	23022C	A4340	No	4
145516	A4358	No	19	23024C	A4340	No	4
145524	A9270	No	21	23026C	A4340	No	4
145540	A4458	No	21	24000	A9270	No	20
145617	A9270	Yes	21	2400SA	A9270	No	21
1512C	A4340	Yes	9	24015	A9270	No	20
1514C	A4340	Yes	9	24035	A9270	No	20
1516C	A4340	Yes	9	2480	A4353	No	18
1518C	A4340	Yes	9	25000	A9270	No	20
1520C	A4340	Yes	9	2500SA	A9270	No	21
1522C	A4340	Yes	9	25015	A9270	No	20
1524C	A4340	Yes	9	2718	A4346	Yes	8
155730	A9270	No	21	2720	A4346	Yes	8
155731	A9270	No	21	2722	A4346	Yes	8
1600	A9270	No	13	2723	A4346	Yes	8
1612	A4338	Yes	8	2819	A4346	Yes	8
1612C	A4340	Yes	9	2821	A4346	Yes	8
1614	A4338	Yes	8	2822	A4346	Yes	8
1614C	A4340	Yes	9	2824	A4346	Yes	8
1616	A4338	Yes	8	2826	A4346	Yes	8
1616C	A4340	Yes	9	3044	A4312	Yes	12
1618C	A4340	Yes	9	3141	A4353	No	18
1619	A4338	Yes	8	3170	*	No	18
1620C	A4340	Yes	9	3215	A4353	Yes	18
1621	A4338	Yes	8	3217	A4353	Yes	18
1622	A4338	Yes	8	3305	A4353	No	18
1622C	A4340	Yes	9	3307	A4353	No	18
1624	A4338	Yes	8	3408	A4353	Yes	18
1624C	A4340	Yes	9	3410	A4353	Yes	18
1626	A4338	Yes	8	3411	A9270	No	18
1628	A4338	Yes	8	3433	A4358	No	19
1631	A4338	Yes	8	3450	A4353	No	18
1814	A4338	Yes	8	3502	A4357	No	15
1817	A4338	Yes	8	3512	A4357	No	15
1818	A4338	Yes	8	3512N	A4357	No	15
2004LL	A4315	No	6	3512V	A4357	No	15
2040SA	A9270	No	20	3513CF	A4357	No	15

Index

Code	HCPCS*	Latex	Page	Code	HCPCS	Latex	Page
3515	A4354	Yes	15	40522L	A4344	No	4
3531	A4354	Yes	15	40524L	A4344	No	4
3532	A4354	No	15	406622	A4353	No	18
3558	A4338	Yes	9	40905A	A9270	No	20
3560	A4338	Yes	9	5027	A4310	No	13
3563	A4338	Yes	9	5029	A4310	No	13
3565	A4338	Yes	9	50514T	*	No	4
3567	A4338	Yes	9	50516ICT	*	No	1
3570	A4338	Yes	9	50516T	*	No	4
3573	A4338	Yes	9	50518ICT	*	No	1
3575	A4338	Yes	9	50518T	*	No	4
3577	A4338	Yes	9	5204	A9270	No	20
3601	A4338	Yes	9	5205	A9270	No	20
3607	A4338	Yes	9	54456	A4333	No	13
3611	A4338	Yes	9	6006	A4351	No	16
3614	A4338	Yes	9	6010	A4344	No	4
3618	A4338	Yes	9	605122IC	A4344	No	1
3619	A4338	Yes	9	605148IC	A4344	No	1
3623	A4338	Yes	9	605163IC	A4344	No	1
3631	A4338	Yes	9	605189IC	A4344	No	1
3685	A4320	No	20	605205IC	A4344	No	1
3716	A4314	Yes	12	605221IC	A4344	No	1
3718	A4314	Yes	12	605247IC	A4344	No	1
400608	A4351	No	16	605262IC	A4344	No	1
400610	A4351	No	16	6108	A4344	No	4
400612	A4351	No	16	6144ICLL	A4314	No	3
400614	A4351	No	16	6144LL	A4315	No	7
400616	A4351	No	16	6146ICLL	A4314	No	3
400618	A4351	No	16	6146LL	A4315	No	7
4014	A9270	No	21	6148ICLL	A4314	No	3
402712	A4338	Yes	8	6148LL	A4315	No	7
402714	A4338	Yes	8	6164LL	A4314	Yes	12
402716	A4338	Yes	8	6165LL	A4314	Yes	12
402718	A4338	Yes	8	6166LL	A4314	Yes	12
402720	A4338	Yes	8	6251	A4357	No	15
402722	A4338	Yes	8	6261	A4357	No	15
402724	A4338	Yes	8	630187IC	A4344	No	1
402726	A4338	Yes	8	630203IC	A4344	No	1
402728	A4338	Yes	8	630229IC	A4344	No	1
402730	A4338	Yes	8	630245IC	A4344	No	1
403712	A4346	Yes	8	630260IC	A4344	No	1
403714	A4346	Yes	8	6308LL	A4357	No	15
403716	A4346	Yes	8	665183IC	A4346	No	1
403718	A4346	Yes	8	665209IC	A4346	No	1
403720	A4346	Yes	8	665225IC	A4346	No	1
403722	A4346	Yes	8	665241IC	A4346	No	1
403724	A4346	Yes	8	67000	A4322	No	20
403726	A4346	Yes	8	67800	A4320	No	20
403728	A4346	Yes	8	67803	A4320	No	20
403730	A4346	Yes	8	68000	A4322	No	20
40516L	A4344	No	4	68800	A4320	No	20
40518L	A4344	No	4	68803	A4320	No	20
40520L	A4344	No	4	68820	A4320	No	20

Index

Code	HCPCS*	Latex	Page	Code	HCPCS	Latex	Page
7000LL	A4315	No	14	8887605163	A4344	No	4
7004ICLL	A4314	No	2	8887605189	A4344	No	4
7006ICLL	A4314	No	2	8887605205	A4344	No	4
7006LL	A4315	No	6	8887605221	A4344	No	4
7008ICLL	A4314	No	2	8887605247	A4344	No	4
7008LL	A4315	No	6	8887605262	A4344	No	4
7014LL	A4314	Yes	10	8887605282	*	No	4
7016LL	A4314	Yes	10	8887605303	*	No	4
7018LL	A4314	Yes	10	8887623125	A4338	Yes	9
7100	A4310	No	13	8887624149	A4338	Yes	9
75000	A4353	No	18	8887630142	*	No	4
75005	A4353	Yes	18	8887630163	*	No	4
75010	A4310	No	18	8887630187	A4344	No	4
75020	A4310	No	18	8887630203	A4344	No	4
75030	A4353	No	18	8887630229	A4344	No	4
75035	A4353	Yes	18	8887630245	A4344	No	4
76000	A4310	No	13	8887630260	A4344	No	4
76010	A4310	No	13	8887630284	*	No	4
76012	A4310	No	13	8887630305	*	No	4
76020	A4310	No	13	8887636002	*	No	13
76030	A4310	No	13	8887648221	A4338	Yes	9
8403	A4352	Yes	17	8887648247	A4338	Yes	9
8404	A4352	Yes	17	8887649229	A4346	Yes	9
8408	A4351	Yes	17	8887649245	A4346	Yes	9
8410	A4351	Yes	17	8887654088	A4338	Yes	8
8412	A4351	Yes	17	8887654104	A4338	Yes	8
8414	A4351	Yes	17	8887660085	A4351	Yes	17
8416	A4351	Yes	17	8887660101	A4351	Yes	17
8418	A4351	Yes	17	8887660127	A4351	Yes	17
8420	A4351	Yes	17	8887660143	A4351	Yes	17
8422	A4351	Yes	17	8887660168	A4351	Yes	17
8424	A4351	Yes	17	8887660184	A4351	Yes	17
8884730300	A4349	Yes	19	8887660218	A4351	Yes	17
8884730600	A4333	No	19	8887660221	A4351	Yes	17
8884731300	A4349	Yes	19	8887660242	A4351	Yes	17
8884731900	A4331	No	19	8887660263	A4351	Yes	17
8884732000	A4349	Yes	19	8887660284	A4351	Yes	17
8884732300	A4349	Yes	19	8887660305	A4351	Yes	17
8884732500	A4349	Yes	19	8887660903	A4351	No	16
8884733100	A4358	Yes	19	8887664161	A4346	No	4
8884733200	A4358	Yes	19	8887664182	A4346	No	4
8884733900	A4358	Yes	19	8887664203	A4346	No	4
8884734600	A4349	Yes	19	8887664224	A4346	No	4
8887600057	A4353	Yes	18	8887664245	A4346	No	4
8887600149	A4334	No	13	8887664266	A4346	No	4
8887600776	A4312	No	6	8887664287	A4346	No	4
8887600784	A4312	No	6	8887665161	*	No	4
8887601121	A4358	No	19	8887665183	A4346	No	4
8887601139	A4358	No	19	8887665209	A4346	No	4
8887603085	A4344	No	4	8887665225	A4346	No	4
8887603101	A4344	No	4	8887665241	A4346	No	4
8887605122	A4344	No	4	8887665262	*	No	4
8887605148	A4344	No	4	8887665283	*	No	4

Index

Code	HCPCS*	Latex	Page	Code	HCPCS	Latex	Page
8887688169	A4346	Yes	8	P4P16UL	A4314	Yes	10
8887688185	A4346	Yes	8	P4P16ULD	*	Yes	10
8887688201	A4346	Yes	8	P4P16XS	A4315	No	2
8887688227	A4346	Yes	8	P4P16XSD	*	No	2
8887688243	A4346	Yes	8	P4P16XTS	A4315	No	2
8887688268	A4346	Yes	8	P4P16XTSD	*	Yes	2
8887689167	A4346	Yes	8	P4P18S	A4315	No	5
8887689183	A4346	Yes	8	P4P18SD	*	No	5
8887689209	A4346	Yes	8	P4P18TS	A4315	No	5
8887689225	A4346	Yes	8	P4P18TSD	*	Yes	5
8887689241	A4346	Yes	8	P4P18UL	A4314	Yes	10
8887689266	A4346	Yes	8	P4P18ULD	*	Yes	10
8888492017	A4351	No	16	P4P18XS	A4315	No	2
8888492025	A4351	No	16	P4P18XSD	*	No	2
8888492033	A4351	No	16	P4P18XTS	A4315	No	2
8888492041	A4351	No	16	P4P18XTSD	*	Yes	2
8888492058	A4351	No	16	P4PA	A4354	No	14
8888492066	A4351	No	16	PA	A4357	No	15
8888492120	A4351	No	16	PA14S	A4315	No	7
8888492122	A4351	No	16	PA14UL	A4314	Yes	11
8888492124	A4351	No	16	PA16S	A4315	No	7
8888492126	A4351	No	16	PA16UL	A4314	Yes	11
8888492128	A4351	No	16	PA18S	A4315	No	7
8888492130	A4351	No	16	PA18UL	A4314	Yes	11
95516	A4335	No	13	PAA	A4354	No	15
P2A	A4357	No	14	PAB	A4357	No	15
P2A16S	A4315	No	5	PP	A4357	No	15
P2A16UL	A4314	Yes	10	PP14S	A4315	No	7
P2A18S	A4315	No	5	PP14SD	*	No	7
P2A18UL	A4314	Yes	10	PP14UL	A4314	Yes	11
P2AA	A4354	No	14	PP14ULD	*	Yes	11
P2K	A4357	No	14	PP14XS	A4315	No	3
P4A	A4357	No	14	PP14XSD	*	No	3
P4A14S	A4315	No	5	PP16S	A4315	No	7
P4A14UL	A4314	Yes	10	PP16SD	*	No	7
P4A16S	A4315	No	5	PP16UCL	*	Yes	11
P4A16UL	A4314	Yes	10	PP16UL	A4314	Yes	11
P4A18S	A4315	No	5	PP16ULD	*	Yes	11
P4AA	A4354	No	14	PP16XS	A4315	No	3
P4P	A4357	No	14	PP16XSD	*	No	3
P4P14S	A4315	No	5	PP18S	A4315	No	7
P4P14SD	*	No	5	PP18SD	*	No	7
P4P14TS	A4315	No	5	PP18UL	A4314	Yes	11
P4P14TSD	*	Yes	5	PP18ULD	*	Yes	11
P4P14UL	A4314	Yes	10	PP18XS	A4315	No	3
P4P14ULD	*	Yes	10	PP18XSD	*	No	3
P4P14XS	A4315	No	2	PPA	A4354	No	15
P4P14XSD	*	No	2				
P4P16S	A4315	No	5				
P4P16SD	*	No	5				
P4P16TS	A4315	No	5				
P4P16TSD	*	Yes	5				
P4P16UCL	*	Yes	10				

* HCPCS code pending verification by PDAC. Please call the Reimbursement Hot Line (800-346-7197 x6076) for update.

The Centers for Medicare and Medicaid Services was the primary source of any HCPCS information contained in this document. Covidien does not represent or guarantee that any HCPCS information contained in this document is complete, accurate, or applicable to any particular patient or third-party payer and disclaims all liability for any consequence resulting from any reliance on such information. The final billing-related decisions for any product must be made by the health care provider considering the medical necessity of the product furnished as well as the requirements of third-party payers and any local, state, or federal laws and regulations that apply to the product furnished. Covidien is providing this information in an educational capacity with the understanding that Covidien is not engaged in rendering legal, accounting, or other professional services. Covidien encourages all health care providers to consult with their own advisors regarding coding and payment.

*HCPCS code pending verification by PDAC. Please call the Reimbursement Hot Line (800-346-7197 x6076) for update.

The products shown in this catalog may not be licensed or registered in all parts of the world.
Please contact your sales representative for availability.

COVIDIEN, COVIDIEN with logo and Covidien logo are U.S. and internationally registered trademarks of Covidien AG.
Other brands are trademarks of a Covidien company. ©2010 Covidien AG.

H6200URO – 15M – 1110

