


ETHICON

Coated VICRYL*

(Polyglactin 910) Suture



Coated VICRYL*

(Polyglactin 910) Braided Suture

Braided

For easy handling and secure knot tying.

Coated

For smooth passage through tissue and easy knot tie down.

Synthetic

For minimal tissue reaction.

Absorbable

By hydrolysis. Absorption of VICRYL* Suture is essentially complete between 56 and 70 days.

Tensile Strength

Post implantation		Approximate % original strength remaining
14 days		75%
21 days	6-0 and larger	50%
21 days	7-0 and smaller	40%
28 days	6-0 and larger	25%

All of the original tensile strength is lost by five weeks post implantation.

Colour

Violet or undyed.

Range

Gauge sizes 8-0 to 2 (U.S.P.) Supplied as needled sutures and ligatures.

Sterilisation










Ethylene Oxide.

VICRYL* Monofilament











Available in gauge sizes 9-0 and 10-0 for Cataract Surgery.

These sutures comply with the requirements of U.S. Pharmacopoeia (U.S.P.) and The European Pharmacopoeia (PhEUR).




Coated VICRYL* (Polyglactin 910) Suture



Needle length and description	Length	Colour	Gauge				
⊙ 1/2 circle ETHIGUARD* Blunt Point			2-0 3	0 3.5	1 4	2 5	4 6
MHB-1 	31mm	75cm	violet	W9984			
CTB-1 	36mm	75cm	violet	W9985			
CTB-1 heavy 	36mm	75cm 90cm	violet violet	W9993	W9994	W9995	
CTB 	40mm	75cm	violet		W9988	W9989	
CTB heavy 	40mm	90cm	violet		W9996	W9997	
BT-3 	45mm	90cm	violet			W9998	W9999
CTXB 	48mm	90cm	violet		W9990	W9991	
○ 3/8 circle Blunt Point			2-0 3	0 3.5	1 4	2 5	4 6
BP-5 	65mm	100cm	violet		W9391		
⊙ 3/8 circle Taper Point			4-0 1.5	3-0 2	2-0 3	0 3.5	1 4
C-1 	13mm	60cm	violet	W9067			

Coated VICRYL* (Polyglactin 910) Suture





Needle length and description	Length	Colour	Gauge							
☉ 3/8 circle Taper Point			4-0 1.5	3-0 2	2-0 3	0 3.5	1 4	2 5	4 6	
BB 	17mm	45cm	violet	W9074						
☾ 1/2 circle Taper Point plus			6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4	
RB-2 	13mm	45cm	undyed	W9981	W9982					
RB-1 plus 	17mm	75cm	violet		W9105	W9106				
SH-2 plus 	20mm	75cm	violet			W9113	W9114			
SH-1 plus 	22mm	8x45cm CR	violet			V794E	V793E	V792E		
JB VISI-BLACK* 	26mm	70cm	violet				W9122	W9123		
SH plus 	26mm	8x45cm CR	violet				V774E	V775E		
SH plus 	26mm	75cm	violet				W9120	W9121		
MH-1 	31mm	75cm 100cm	violet violet					W9636 W9637		
MH-1 plus 	31mm	75cm	violet				W9130	W9136	W9138	W9213

Coated VICRYL* (Polyglactin 910) Suture








Needle length and description	Length	Colour	Gauge							
			6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4	
⊖ 1/2 circle Taper Point plus										
MH plus 	36mm	8x70cm CR	violet						V778G	
MH plus 	36mm	75cm 90cm	violet violet					W9140 W9440	W9141 W9441	
LH 	40mm	75cm	violet					W9150	W9154	

Needle length and description	Length	Colour	Gauge							
			2-0 3	0 3.5	1 4	2 5	4 6	5 7	7 9	
⊙ 1/2 circle Taper Point (heavier wire diameters)										
UCL 	13mm	70cm	violet	V114H						
MO-5 	31mm	75cm	violet		W9215	W9216				
CT-1 	36mm	8x45cm CR	violet		V740G	V741G				
CT 	40mm	75cm 90cm	violet violet		W9230 W9430	W9231 W9431				











Coated VICRYL* (Polyglactin 910) Suture


Needle length and description	Length	Colour	Gauge				
◎ 1/2 circle Taper Point (heavier wire diameters)			2-0 3	0 3.5	1 4	2 5	4 6
MO-45 	45mm	75cm	violet		W9244	W9245	W9246
CTX 	48mm	75cm 90cm	violet violet		W9250 W9450	W9251 W9451	W9252 W9452
TP-1 	65mm	75cm	violet	W9158			
MO-80 	80mm	100cm	violet		W9289		

Coated VICRYL* (Polyglactin 910) Suture





Needle length and description	Length	Colour	Gauge				
⊙ 5/8 circle Taper Point			2-0 3	0 3.5	1 4	2 5	4 6
 <div style="display: flex; justify-content: space-between; width: 100%;"> UR-5 36mm 70cm violet W9160 </div>							
 <div style="display: flex; justify-content: space-between; width: 100%;"> UR-4 40mm 75cm violet W9165 </div>							
⊙ J Needle Taper Point				3-0 2	2-0 3	0 3.5	1 4
 <div style="display: flex; justify-content: space-between; width: 100%;"> ASH-1 22mm 75cm violet W9223 </div>							
 <div style="display: flex; justify-content: space-between; width: 100%;"> ASH-35 36mm 75cm violet W9220 W9221 </div>							
⊙ 3/8 circle CC Needle			6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
 <div style="display: flex; justify-content: space-between; width: 100%;"> CC-1 13mm 75cm violet JV1012 JV1013 JV1014 </div>							
⊙ 1/2 circle CC Needle			6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
 <div style="display: flex; justify-content: space-between; width: 100%;"> CC-20 20mm 75cm violet W9118 </div>							
⊙ 3/8 circle TAPERCUT*			6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
 <div style="display: flex; justify-content: space-between; width: 100%;"> V-4 17mm 75cm undyed W9824 W9825 W9826 </div>							

Coated VICRYL* (Polyglactin 910) Suture







Needle length and description	Length	Colour	Gauge						
			4-0 1.5	3-0 2	2-0 3	0 3.5	1 4	2 5	
▼ 1/2 circle TAPERCUT*									
V-6 	22mm	70cm 75cm	undyed violet	V231H		W9828			
V-7 	26mm	75cm	violet		W9180	W9181			
V-7 heavy 	26mm	75cm	violet			W9350			
V-30 	31mm	75cm	violet			W9360	W9361	W9362	
V-34 	36mm	75cm 90cm	violet violet		V516H	W9363 W9463	W9364 W9464		
V-34 heavy 	36mm	75cm 90cm	violet violet					W9365 W9465	W9357
V-37 	40mm	75cm 75cm 90cm	violet undyed violet			W9900 W9466	W9901 W9467	W9368 W9468	
V-37 heavy 	40mm	75cm 75cm	violet undyed			W9359 W9904	W9372 W9905	W9373 W9906	
V-39 	45mm	75cm	violet			W9375	W9376	W9377	W9378
V-80 	80mm	75cm	violet				W9369	W9379	

Needle length and description	Length	Colour	Gauge						
			4-0 1.5	3-0 2	2-0 3	0 3.5	1 4	2 5	
▼ J Needle TAPERCUT*									
V-38 	31mm	75cm	violet					W9335	




Coated VICRYL* (Polyglactin 910) Suture





Needle length and description	Length	Colour	Gauge							
			4-0 1.5	3-0 2	2-0 3	0 3.5	1 4	2 5	4 6	
▼ 1/2 circle Reverse Cutting										
 <p style="text-align: center;">OS-4</p>	22mm	75cm	violet						W9287	
 <p style="text-align: center;">CP</p>	40mm	75cm 90cm	violet violet				W9320 W9420	W9321 W9421		
 <p style="text-align: center;">CPX</p>	48mm	75cm 90cm	violet violet				W9295	W9296 W9496	W9297 W9497	
▲ 1/2 circle TROCAR POINT										
 <p style="text-align: center;">TR-40</p>	40mm	75cm	violet					W9340		






Coated VICRYL* (Polyglactin 910) Suture

Needle length and description	Length	Colour	Gauge				
▼ 3/8 circle Reverse Cutting MultiPass* PRIME* Needle			6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
P-1 	11mm	45cm	undyed	MPV489H	MPV490H		
P-3 	13mm	45cm	undyed	MPV492H	MPV493H	MPV494H	
PS-3 	16mm	45cm	undyed		MPV500H		
PS-2 	19mm	45cm	undyed			MPV496H	MPV497H
PS-1 	24mm	45cm	undyed				MPV683H MPV684H
PS 	26mm	70cm	undyed				MPV9526H MPV9527H












Coated VICRYL* (Polyglactin 910) Suture

Needle length and description		Length	Colour	Gauge						
▼ 3/8 circle Reverse Cutting PRIME* Needle				6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4
 P-1	11mm	45cm	undyed	W9500T	W9501T					
 PS-2	19mm	45cm	undyed		W9514T	W9515T	W9516T			
 PS	26mm	45cm 75cm	undyed undyed				W9525T W9526T	W9527T		












▼ 3/8 circle Reverse Cutting				6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4
 FS-2	19mm	45cm 45cm	violet undyed			W9386 W9951				
 FS	26mm	45cm 75cm	violet undyed				W9388 W9890			
 FSLX	36mm	75cm	violet					W9390		
 LS-1	45mm	75cm	violet					W9282		

▲ 3/8 circle Conventional Cutting PRIME* Needle				6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4
 PC-1	13mm	45cm	undyed	W9831T	W9832T					
 PC-3	16mm	45cm	undyed		W9505T	W9506T	W9507T			
 PC-5	19mm	45cm 75cm	undyed undyed			W9510T W9570T	W9511T W9571T			
 PC-25	26mm	45cm 75cm	undyed undyed			W9520T W9580T	W9521T W9581T	W9522T W9582T		
 PCLX	36mm	75cm	undyed				W9531T	W9532T		


Coated VICRYL* (Polyglactin 910) Suture


Needle length and description		Length	Colour	Gauge				
▲ 3/8 circle Conventional Cutting				6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
FS-3 CONV 	16mm	75cm	violet		W9442	W9443	W9444	
▲ 1/2 circle Conventional Cutting				6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
X-1 CONV 	22mm	45cm 45cm	violet undyed			W9471	W9472 W9741	
▲ 5/8 circle Conventional Cutting				6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
UX-25 	26mm	75cm	undyed				W9730	
▼ Compound Curve Reverse Cutting - PRIME* Needle				6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
PS-4C 	16mm	45cm	undyed		W9436	W9437		
▼ Straight Cutting				6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
SC-1 	13mm	45cm	violet			W9410		
KS-45 	45mm	45cm	undyed			W9711		
KS-55 	55mm	75cm	undyed					W9713
KS 	60mm	75cm 75cm	violet undyed				W9719	W9418 W9718
KS fine 	60mm	75cm 75cm	violet undyed				W9415 W9717	
▼ 5/8 circle TAPERCUT* - P Needle				6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
UV-17 	17mm	75cm	undyed			W9509T		
▼ 5/8 circle TAPERCUT*				6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
UV-17 	17mm	75cm	undyed			W9820		

Coated VICRYL* (Polyglactin 910) Suture

Needle length and description	Length	Colour	Gauge						
▼ 3/8 circle MICRO-POINT spatula			8-0 0.4	7-0 0.5	6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
TG140-8 	6.5mm	15cm	violet	W9558					
TG140-8 	6.5mm	30cm 45cm	violet violet	W9560 W9559	W9561				
▼ 1/2 circle MICRO-POINT spatula			8-0 0.4	7-0 0.5	6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
TG175-8 	7mm	30cm 45cm	violet violet	W9564	W9565				
▼ 3/8 circle MICRO-POINT Reverse Cutting			8-0 0.4	7-0 0.5	6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
G-1 	11mm	45cm	violet		W9537				
▼ 1/2 circle MICRO-POINT Reverse Cutting			8-0 0.4	7-0 0.5	6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
G-7 	8mm	30cm	violet	W9545					
▼ 1/4 circle SPATULATED			8-0 0.4	7-0 0.5	6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
S-14 	8mm	20cm 45cm 45cm	violet violet undyed		W9562 W9552 W9752	W9553 W9753			
S-4 	8mm	45cm	violet				V572G		
▼ 3/8 circle SPATULATED			8-0 0.4	7-0 0.5	6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
SPAT-10 	11mm	45cm	undyed		W9760	W9761	W9762		
▼ 1/2 circle SPATULATED			8-0 0.4	7-0 0.5	6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
S-22 	8mm	20cm 45cm 45cm	violet violet undyed		W9563 W9566 W9756	W9567			
▼ 5/8 circle TAPERCUT			8-0 0.4	7-0 0.5	6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
UV-10 	10mm	30cm	violet		W9172				
▼ 3/8 circle Reverse Cutting			8-0 0.4	7-0 0.5	6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
G-3 	13mm	45cm	violet			JV557	JV559		

Coated VICRYL* (Polyglactin 910) Suture

Needle length and description	Length	Colour	Gauge							
⊙ 3/8 circle Taper Point			8-0 0.4	7-0 0.5	6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3	
BV-1 	9.3mm	30cm	violet	W9577		W9575				

Needle length and description	Length	Colour	Gauge						
⊙ Ski Needle Taper Point			8-0 0.4	7-0 0.5	6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
SKI-22 VISI-BLACK* 	22mm	20cm	violet					E9903S	E9902S

⊙ Straight Taper Point			8-0 0.4	7-0 0.5	6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
ST-4 <u>VISI-BLACK*</u>	19mm	20cm	violet				MIC103H		

Coated VICRYL* (Polyglactin 910) Suture

Ligatures and Non-needed Sutures

1.5m lengths

Gauge		violet 1.5m lengths 12 sachets per box
5-0	1	W9022
4-0	1.5	W9023
3-0	2	W9024
2-0	3	W9025
0	3.5	W9026
1	4	W9027
2	5	W9028

SUTUPAK* Pre-Cut Sutures

Gauge		violet 6 x 45cm lengths	violet 10 x 50cm lengths	violet 5 x 70cm lengths
4-0	1.5		V9048H	V1214E
3-0	2	V1225H	V9049H	V1215E
2-0	3	V1226H	V9050H	V1216E
0	3.5	V1227E	V9051H	V1217E
1	4	V1228E	V9052E	V1218H

LIGAPAK* Dispenser Reel (radio opaque)







Gauge		violet 2.5m per reel 12 reels per box
3-0	2	W9030
2-0	3	W9020
0	3.5	W9000
1	4	W9001

ENDOLOOP* Ligature

Pre-tied ligatures with Cannula for use during endoscopic surgery

Gauge		violet 45cm 12 sachets per box	Cannula length
0	3.5	EJ10C	29cm
0	3.5	MIC511G	38cm

Monofilament VICRYL* (Polyglactin 910) Suture

Needle length and description	Length	Colour	Gauge						
 3/8 circle CS ULTIMA* spatula			10-0 0.2	9-0 0.3	8-0 0.4	7-0 0.5	6-0 0.7	5-0 1	
CS140-6 	6.5mm	10cm	violet	V960G					
 3/8 circle Advanced MICRO-POINT spatula				10-0 0.2	9-0 0.3	8-0 0.4	7-0 0.5	6-0 0.7	5-0 1
TG140-6 	6.5mm	10cm	violet	W1702	W1703				
 3/8 circle Taper Point				10-0 0.2	9-0 0.3	8-0 0.4	7-0 0.5	6-0 0.7	5-0 1
BV-1 	9.3mm	15cm	violet	W9579					