


ETHICON

PDS* Plus

Antibacterial (Polydioxanone) Suture



PDS* Plus

Antibacterial (Polydioxanone) Suture

All the benefits of PDS*II (Polydioxanone) Monofilament Suture, with the additional benefit of an antibacterial.

Ideal for -

- Abdominal closure
- Fascia
- Pancreatic tissue
- Small bowel, colon
- Oesophagus
- Tendon
- Joint capsule
- Ligament

STRENGTH RETENTION












3-0 and larger	4-0 and smaller
80% at 2 weeks	60% at 2 weeks
70% at 4 weeks	40% at 4 weeks
60% at 6 weeks	35% at 6 weeks

Absorption complete at 182 - 232 days









Providing -






- **Protection** of suture from bacterial colonisation
- **Strength** when you need it most
- Wound support that you can **trust**
- **Security** of extended tissue support

PDS* Plus Antibacterial (Polydioxanone) Suture



Needle length and description	Length	Colour	Gauge				
☉ 1/2 circle ETHIGUARD* Blunt Point			3-0 2	2-0 3	0 3.5	1 4	2 5
MHB-1 	31mm	70cm	violet	PDP9959H			
BT-3 	45mm	90cm	violet			PDP9965H	
CTXB 	48mm	150cm loop	violet		PDP9966T	PDP9967T	
☉ 3/8 circle Taper Point			7-0 0.5	6-0 0.7	5-0 1	4-0 1.5	3-0 2
BV-175-6 	8mm	45cm	violet	PDP9091H			
BV-1 	9.3mm	70cm	violet	PDP1702H			
BV 	11mm	45cm 70cm	violet violet	PDP9093H PDP1002H			
C-1 	13mm	70cm	violet	PDP1721H			
C-1 	13mm	70cm	violet		PDP1013H		
C-1 	13mm	70cm	violet	PDP1721H			
BB 	17mm	70cm	violet		PDP9073H	PDP9077H	
☉ 1/2 circle Taper Point plus			6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
RB-2 	13mm	45cm	violet	PDP9100H	PDP9101H	PDP9102H	


PDS* Plus Antibacterial (Polydioxanone) Suture



Needle length and description	Length	Colour	Gauge							
			6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4	
⊖ 1/2 circle Taper Point plus										
 RB-2 13mm	70cm	violet		PDP9201H						
 RB-1 plus 17mm	90cm	violet		PDP9108H	PDP9109H					
 SH-2 plus 20mm	70cm	violet			PDP9115H	PDP9116H				
 SH 26mm	70cm	violet							PDP9205H	
 SH plus 26mm	70cm	violet					PDP9124H	PDP317H		
 MH-1 plus 31mm	70cm	violet				PDP9131H	PDP9132H	PDP9133H	PDP2017H	PDP9211H
 MH-1 plus 31mm	90cm	violet					PDP9134H			
 LH 40mm	70cm	violet					PDP9152H	PDP9151H		

Needle length and description	Length	Colour	Gauge							
			2-0 3	0 3.5	1 4	2 5	4 6	5 7	7 9	
⊙ 1/2 circle Taper Point (heavier wire diameters)										
 CT-2 26mm	70cm	violet	PDP333H	PDP334H	PDP335H					
 CT-1 36mm	70cm	violet	PDP339H	PDP340H	PDP341H					
 CT 40mm	70cm 90cm	violet violet		PDP9233H	PDP9234H					
 MO-2 40mm	150cm loop	violet		PDP9236T	PDP9237T					
 MO-45 45mm	90cm	violet				PDP9248H				







PDS* Plus Antibacterial (Polydioxanone) Suture


Needle length and description	Length	Colour	Gauge				
⊙ 1/2 circle Taper Point (heavier wire diameters)			2-0 3	0 3.5	1 4	2 5	4 6
 <p style="text-align: center;">CTX</p> <p style="text-align: right;">48mm</p>	90cm	violet		PDP9255H	PDP9254H		
 <p style="text-align: center;">CTX heavy</p> <p style="text-align: right;">48mm</p>	150cm loop	violet		PDP9261T	PDP9262T		




Needle length and description	Length	Colour	Gauge				
⊙ J Needle Taper Point			3-0 2	2-0 3	0 3.5	1 4	2 5
 <p style="text-align: center;">ASH-35</p> <p style="text-align: right;">36mm</p>	70cm	violet		PDP9219H			

Needle length and description	Length	Colour	Gauge						
⊙ 3/8 circle CC Needle			7-0 0.5	6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5
 <p style="text-align: center;">CC</p> <p style="text-align: right;">9.3mm</p>	70cm	violet		PDP1712H					
 <p style="text-align: center;">CC-1</p> <p style="text-align: right;">13mm</p>	70cm	violet		PDP1032H					








PDS* Plus Antibacterial (Polydioxanone) Suture



Needle length and description	Length	Colour	Gauge						
▼ 1/2 circle TAPERCUT*			3-0 2	2-0 3	0 3.5	1 4	2 5	4 6	5 7
 V-7 26mm	70cm	violet	PDP9179H	PDP9184H					
 V-34 36mm	70cm 90cm	violet violet		PDP9380H		PDP9381H	PDP9352H		
 V-34 heavy 36mm	90cm	violet				PDP9355H	PDP9382H		
 V-37 40mm	90cm	violet				PDP9366H	PDP9370H		
 V-39 45mm	90cm	violet					PDP9385H		
 V-40 48mm	90cm	violet				PDP9395H	PDP9396H		


▼ J Needle TAPERCUT*			3-0 2	2-0 3	0 3.5	1 4	2 5	4 6	5 7
 V-38 31mm	70cm	violet				PDP9333H	PDP9334H		



Needle length and description	Length	Colour	Gauge				
▼ 1/2 circle Reverse Cutting			3-0 2	2-0 3	0 3.5	1 4	2 5
 OS-4 22mm	70cm	violet				PDP9286H	
 CP 40mm	90cm	violet				PDP9324H	
 CPX 48mm	90cm	violet				PDP9299H	




PDS* Plus Antibacterial (Polydioxanone) Suture

Needle length and description		Length	Colour	Gauge				
▼ 3/8 circle Reverse Cutting MultiPass* PRIME* Needle				6.0 0.7	5-0 1.0	4-0 1.5	3-0 2.0	2.0 3.0
 P-6	8mm	45cm	undyed	PDP9860H				
 P-1	11mm	45cm 45cm	undyed violet	PDP489H	PDP490H PDP487H			
 P-3	11mm	45cm	undyed		PDP493H	PDP494H		
 PS-3	16mm	45cm	undyed		PDP500H	PDP501H		
 PS-2	19mm	45cm	undyed			PDP496H	PDP497H	
 PS-1	24mm	45cm	undyed			PDPZ682H	PDP683H	
 PS	26mm	70cm	undyed			PDP9624H		PDP9626H


▲ 3/8 circle Conventional Cutting MultiPass* PRIME* Needle				6.0 0.7	5-0 1.0	4-0 1.5	3-0 2.0	2.0 3.0
 PC-3	16mm	45cm	undyed		PDP9865H	PDP9867H		
 PC-5	19mm	45cm	undyed			PDP9615H		


▼ 1/2 circle Reverse Cutting MultiPass* PRIME* Needle				6.0 0.7	5-0 1.0	4-0 1.5	3-0 2.0	2.0 3.0
 PS-4	16mm	45cm	undyed			PDP507H		



▼ 3/8 circle Reverse Cutting				6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4
 FS-2	19mm	45cm	undyed			PDP9950H				
 FS	26mm	45cm	undyed			PDP9957H				

▲ 3/8 circle Conventional Cutting PRIME* Needle				6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4
 PC-1	13mm	45cm	undyed	PDP9861H	PDP9863H					
 PC-25	26mm	45cm	undyed				PDP9625H			
 PCLX	36mm	70cm	undyed					PDP9631H		

PDS* Plus Antibacterial (Polydioxanone) Suture

Needle length and description	Length	Colour	Gauge						
▲ 3/8 circle Conventional Cutting			6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4
 FS-3 CONV 16mm	70cm	undyed		PDP9733H	PDP9734H				

Needle length and description	Length	Colour	Gauge				
▲ 1/2 circle Conventional Cutting			6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
 CX-1 22mm	45cm	undyed				PDP9740H	

Needle length and description	Length	Colour	Gauge				
▼ Straight Cutting			6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
 KS 60mm	70cm	undyed					PDP9715H
 KS fine 60mm	70cm	undyed			PDP622H	PDP9714H	