


ETHICON

PROLENE*

Monofilament Polypropylene Suture



PROLENE*

Monofilament Polypropylene Suture

Non Absorbable

Provides permanent tensile strength retention in tissue, even in the presence of infection.

Exceptionally smooth

For easy passage through tissue.

Monofilament

Smooth surface will not support bacterial growth.

Synthetic

Elicits only minimal tissue reaction.

Pliable

For easy knot tying.

Colour

Blue or clear.

Range

Gauge sizes 10-0 to 1 (U.S.P.). Supplied as needled sutures.

Sterilisation

Ethylene Oxide.



As with other synthetic sutures, knot security requires the standard technique of square ties with additional throws as indicated by surgical circumstances and the experience of the operator.














Caution:

Care should be taken to avoid damaging the surface of the material with surgical instruments, as this could lead to fracture of the material.














These sutures comply with the requirements of U.S. Pharmacopoeia (U.S.P.) and The European Pharmacopoeia (PhEUR).

PROLENE* Monofilament Polypropylene Suture












Needle length and description	Length	Colour	Gauge							
⊙ 1/2 circle ETHIGUARD* Blunt Point			2-0 3	0 3.5	1 4	2 5	4 6	5 7	7 9	
 CTB <div style="float: right; margin-right: 50px;">40mm</div>	100cm	blue								W8998
 BT-3 <div style="float: right; margin-right: 50px;">45mm</div>	100cm	blue								W8999

⊙ 3/8 circle Taper Point			8-0 0.4	7-0 0.5	6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
 BV130-5 <div style="float: right; margin-right: 50px;">6.5mm</div>	13cm	blue	W2775						
 BV130-5 MultiPass* Needles <div style="float: right; margin-right: 50px;">6.5mm</div>	45cm	blue	W2777						
 BV130-5 <div style="float: right; margin-right: 50px;">6.5mm</div>	60cm	blue	F1890						
 BV175-6 VISI-BLACK* <div style="float: right; margin-right: 50px;">8mm</div>	60cm	blue	8841H	8335H					
 BV175-6 MultiPass* Needles <div style="float: right; margin-right: 50px;">8mm</div>	60cm	blue	8741H	8735H					
 BV175-7 VISI-BLACK* <div style="float: right; margin-right: 50px;">8mm</div>	60cm	blue	8967H	8966H					
 BV175-7 MultiPass* Needles <div style="float: right; margin-right: 50px;">8mm</div>	60cm	blue	8767H	8766H					
 BV175-8 MultiPass* Needles <div style="float: right; margin-right: 50px;">9.3mm</div>	60cm	blue		W8702					
 BV175-8 STRAIGHT PACK <div style="float: right; margin-right: 50px;">9.3mm</div>	4x60cm	blue		W8702S					
 BV-1 VISI-BLACK* <div style="float: right; margin-right: 50px;">9.3mm</div>	60cm	blue		W8304	W8305				
 BV-1 <div style="float: right; margin-right: 50px;">9.3mm</div>	45cm 60cm	blue blue			W8711 W8712				
 BV <div style="float: right; margin-right: 50px;">11mm</div>	60cm	blue			W8597				
 C-1 VISI-BLACK* <div style="float: right; margin-right: 50px;">13mm</div>	60cm 75cm	blue blue			W8307 W8316	W8321			

PROLENE* Monofilament Polypropylene Suture




Needle length and description	Length	Colour	Gauge						
☉ 3/8 circle Taper Point			8-0 0.4	7-0 0.5	6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
 C-1 13mm	60cm 75cm	blue blue		W8725	W8718 W8706				
 C-1 13mm	75cm	blue				W8890			
☉ 1/2 circle Taper Point			8-0 0.4	7-0 0.5	6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
 TF-6 6.5mm	75cm	blue		EH7813E					
 RB-2 VISI-BLACK* 13mm	75cm	blue				W8310			
 RB-2 13mm	75cm	blue				W8710			
 RB-1 VISI-BLACK* 17mm	90cm	blue				W8330	W8329		
 RB-1 17mm	90cm	blue				W8556	W8557	W8558	
 SH-2 VISI-BLACK* 20mm	90cm	blue					W8340		
 SH-2 20mm	90cm	blue					W8761		
 SH-1 22mm	75cm	blue						W8770	
 SH VISI-BLACK* 26mm	90cm	blue					W8355	W8354	
 SH 26mm	75cm 90cm 120cm	blue blue blue				F1821	W8521	W8522 W8534	
 JMH-1 VISI-BLACK* 31mm	120cm	blue						X1295H	X1296H











PROLENE* Monofilament Polypropylene Suture



Needle length and description		Length	Colour	Gauge					
⊙ 1/2 circle Taper Point				8-0 0.4	7-0 0.5	6-0 0.7	5-0 1	3-0 2	2-0 3
MH-1 	31mm	75cm	blue						W295
MH-1 	31mm	90cm	blue					W8525	W8526
⊙ 1/2 circle Taper Point (heavier wire diameters)				2-0 3	0 3.5	1 4	2 5	5 7	7 9
MO-5 	31mm	100cm	blue		W8430				
CT 	40mm	100cm	blue			W742			
MO-2 	40mm	100cm	blue		W8434				
TP-1 	65mm	100cm	blue			W8401			
⊙ Straight Taper Point				2-0 3	0 3.5	1 4	2 5	5 7	7 9
ST-70 	70mm	75cm	blue	W8400					
⊙ 3/8 circle CC Needle				8-0 0.4	7-0 0.5	6-0 0.7	5-0 1	3-0 2	2-0 3
CC175-6 VISI-BLACK* 	8mm	60cm	blue	W8101					
CC MultiPass* Needles 	9.3mm	60cm 75cm	blue blue	W8703	W8704	W8807 F1832			
CC HEMO-SEAL* MultiPass* Needles 	9.3mm	60cm	blue				W8662		
CC-11 MultiPass* Needles 	11mm	60cm	blue		W8801	W8802			

PROLENE*

Monofilament Polypropylene Suture









Needle length and description		Length	Colour	Gauge						
	3/8 circle CC Needle			8-0 0.4	7-0 0.5	6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
	13mm	60cm 75cm	blue blue		W8813	W8707 W8814				
	13mm	75cm	blue				W8721			

	1/2 circle CC Needle			8-0 0.4	7-0 0.5	6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
	11mm	60cm	blue				W8803			
	13mm	60cm 75cm	blue blue			W8815	W8816			
	13mm	60cm 75cm	blue blue				W8664	W8665		
	16mm	90cm	blue				W8830	W8831		
	16mm	90cm	blue						W8667	
	20mm	90cm	blue					W8840		
	26mm	90cm 120cm	blue blue					W8845	W8844 W8895	W8843 W8896
	31mm	90cm	blue						W8849	W8850
	40mm	90cm	blue						W8851	W8852

	MULTI-CURVE CC Needle			8-0 0.4	7-0 0.5	6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
	9.3mm	75cm	blue		EH8021E	EH8025H				



PROLENE*


Monofilament Polypropylene Suture

Needle length and description	Length	Colour	Gauge						
Ⓜ MULTI-CURVE CC Needle			8-0 0.4	7-0 0.5	6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
ACC-1 VISI-BLACK* 	13mm	75cm	blue			EH8030H	EH8036H		
Ⓜ 3/8 circle TAPERCUT*			8-0 0.4	7-0 0.5	6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
V-4 	17mm	75cm	blue			F2829			
V-26 	26mm	75cm 120cm	blue blue				F2845 F2850		
Ⓜ 1/2 circle TAPERCUT*			4-0 1.5	3-0 2	2-0 3	0 3.5	1 4	2 5	4 6
V-5 	17mm	90cm	blue	W8935	W8936	W8937			
V-7 	26mm	90cm	blue		W8976	W8977			
V-34 	36mm	100cm	blue				W486		
V-39 	45mm	100cm	blue					W8450	
Ⓜ J Needle TAPERCUT*			4-0 1.5	3-0 2	2-0 3	0 3.5	1 4	2 5	4 6
V-38 	31mm	100cm	blue				W982		

PROLENE*






Monofilament Polypropylene Suture

Needle length and description	Length	Colour	Gauge						
▼ 1/2 circle Reverse Cutting			4-0 1.5	3-0 2	2-0 3	0 3.5	1 4	2 5	4 6
CP 	40mm	100cm	blue				W8470		
CPX 	48mm	100cm	blue				W8475		








Needle length and description	Length	Colour	Gauge						
▲ 1/2 circle TROCAR POINT			4-0 1.5	3-0 2	2-0 3	0 3.5	1 4	2 5	4 6
TR-30 	31mm	100cm	blue				W8440		



PROLENE*




Monofilament Polypropylene Suture

Needle length and description	Length	Colour	Gauge						
▼ 3/8 circle Reverse Cutting MultiPass* PRIME* Needle			7-0 0.5	6-0 0.7	5-0 1.0	4-0 1.5	3-0 2.0	2-0 3.0	
P-1 	11mm	45cm 75cm	blue blue	MPP2854H	MPP8697H				
P-3 	13mm	45cm 45cm 75cm 75cm	blue clear blue clear		MPP8695H	MPP8698H MPP8605H MPP2831H MPP2858H	MPP8699H MPP8604H MPP2847H MPP2860H		
PS-3 	16mm	45cm 75cm 75cm	blue blue clear		MPP2832H	MPP8681H MPP2833H MPP2859H	MPP7764H		
PS-2 	19mm	45cm 75cm 75cm	blue blue clear			MPP8686H	MPP8682H MPP2834H MPP2861H	MPP8687H MPP2835H MPP2862H	
PS-1 	24mm	45cm	blue				MPP7758H	MPP7757H	MPP7756H





PROLENE* Monofilament Polypropylene Suture



Needle length and description		Length	Colour	Gauge						
▼ 3/8 circle Reverse Cutting PRIME* Needle				6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4
 P-6	8mm	45cm	blue	W8870T						
 P-1	11mm	45cm	blue	W8868T						
 P-3	13mm	45cm 45cm 75cm	blue clear clear	W8871T	W8872T F2858	8604G F2860				
 PS-3	16mm	45cm 75cm	blue clear		W8882T F2859	W8884T				
 PS-2	19mm	45cm 75cm	blue clear		W8014T	W8015T F2861	F2862			
 PS	26mm	45cm	blue				W8025T	W8026T		
 2 locking buttons	26mm	45cm	blue				W8625			


▼ 3/8 circle Reverse Cutting				6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4
 FS-2	19mm	45cm	blue			W8683				
 FS	26mm	45cm	blue				W8684			

▲ 3/8 circle Conventional Cutting PRIME* Needle				6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4
 PC-3	16mm	45cm	blue	W8005T	W8006T	W8007T				
 PC-5	19mm	45cm	blue		W8010T	W8011T	W8878T			
 PC-25	26mm	45cm 45cm	blue clear	W8886T	W8881T	W8020T W534	W8021T			










PROLENE* Monofilament Polypropylene Suture


Needle length and description		Length	Colour	Gauge					
▲ 3/8 circle Conventional Cutting				6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3	1 4
FS-3 CONV 	16mm	75cm	blue				W8549		
PCLX 	36mm	75cm	blue				W626		
PCLX 	36mm	75cm	blue				W8626		
CEC 135-40 	40mm	100cm	blue					W8689	



▲ 3/8 circle Conventional Cutting (Hand Held)				6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3	1 4
KP-3 CONV 	90mm	2x50cm	blue						W2989
KP-3 CONV 	90mm	90cm	blue						W987


▲ 1/2 circle Conventional Cutting				6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3	1 4
X-1 CONV 	22mm	45cm	blue				W8731		



PROLENE* Monofilament Polypropylene Suture


Needle length and description	Length	Colour	Gauge							
			6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5		
▼ Straight Cutting										
 KS-55	55mm	90cm	blue						W8571	
 KS	60mm	75cm	blue						W631	
2 beads & 2 collars  KS	60mm	75cm	blue						W8631	
2 locking buttons  KS	60mm	75cm	blue						W8623	
 KS fine	60mm	75cm	blue					W630		
2 beads & 2 collars  KS fine	60mm	75cm	blue					W8630		
2 locking buttons  KS fine	60mm	75cm	blue					W8621		
 KS fine	60mm	75cm	blue					W621		
 KS-75	75mm	100cm	blue					W751	W752	W753

Needle length and description	Length	Colour	Gauge							
			10-0 0.2	9-0 0.3	8-0 0.4	7-0 0.5	6-0 0.7	5-0 1	4-0 1.5	
▼ 3/8 circle Advanced MICRO-POINT spatula										
 TG140-6	6.5mm	30cm	blue	W1777						

Needle length and description	Length	Colour	Gauge							
			10-0 0.2	9-0 0.3	8-0 0.4	7-0 0.5	6-0 0.7	5-0 1	4-0 1.5	
▼ 3/8 circle MICRO-POINT spatula										
 TG140-8	6.5mm	23cm	blue	W1711	W1709					
 TG140-8	6.5mm	30cm	blue	W1710						




Needle length and description	Length	Colour	Gauge							
			10-0 0.2	9-0 0.3	8-0 0.4	7-0 0.5	6-0 0.7	5-0 1	4-0 1.5	
▼ 1/4 circle TAPERCUT*										
 CIF-4	13mm	20cm	blue	788G						

Needle length and description	Length	Colour	Gauge							
			10-0 0.2	9-0 0.3	8-0 0.4	7-0 0.5	6-0 0.7	5-0 1	4-0 1.5	
▼ Straight MICRO-POINT										
 STC-6	16mm	15cm loop	blue	EH7900G						
 STC-6	16mm	23cm	blue	W1713						

Needle length and description	Length	Colour	Gauge							
			10-0 0.2	9-0 0.3	8-0 0.4	7-0 0.5	6-0 0.7	5-0 1	4-0 1.5	
▼ 1/4 circle Spatula										
 CTC-6L	13mm	10cm loop	blue	EH7896G						

PROLENE*

Monofilament Polypropylene Suture

Needle length and description	Length	Colour	Gauge							
⊙ 3/8 circle Taper Point			10-0 0.2	9-0 0.3	8-0 0.4	7-0 0.5	6-0 0.7	5-0 1	4-0 1.5	
BV75-3 	3.8mm	13cm	blue	W2790						
BV100-4 	5mm	13cm	blue	W2794	W2783					
BV130-4 	5mm	13cm	blue		W2780					

PROLENE*

Monofilament Polypropylene Suture

ENDOLOOP* Ligature

Pre-tied ligatures with Cannula for use during endoscopic surgery

Gauge		blue 45cm
2-0	3	EH464G